


**COSMOSIL**



High Performance  
Liquid Chromatography

**Genore**  
chromatografia

ul. Dynarska 1/23, 01-493 Warszawa  
tel. 22 40 107 34 do 35, fax: -36  
e-mail: [info@genore.pl](mailto:info@genore.pl)  
[www.genore.pl](http://www.genore.pl)

## CORPORATE PROFILE

Nacalai Tesque dates back to 1846 when the company's founder Mansuke Nakarai opened Nakarai Mansuke Shoten, Ltd., an apothecary selling traditional Japanese and Chinese medicines. In 1958, this company's reagent department became an independent company, Nakarai Chemicals, Ltd.

The company has since dedicated itself to expanding its corporate base and has strived to be an enterprise that our customers always rely on, while taking pride in its contribution to scientific and industrial development.


Making the most of this 30-year history and as a step toward the future, Nakarai Chemicals changed its corporate name to Nacalai Tesque, Inc. in 1988. At Nacalai Tesque, we have fostered a corporate commitment to the pursuit of reliable quality and the creation of products of real value, while serving as a vital link between humanity and science. Centering around research chemicals, the fields of our activities include fine chemicals, diagnostics and related laboratory equipment and supplies.

The pace of scientific and technological progress in every industrial field is rapidly accelerating, and all business partners and affiliates are required to provide even more diversified and advanced expertise.

It is our corporate policy to strive for our lofty ideals for excellence while respecting our long history and tradition. We consider it our mission to maintain close contact with our customers by offering reliable quality in all our products, information and services, and by making full use of the knowledge and experience of our staff.

### The 30<sup>th</sup> Anniversary of COSMOSIL HPLC Column

In 2010, Nacalai Tesque reached another important milestone in the company's history with the 30<sup>th</sup> anniversary of COSMOSIL brand HPLC column. For the past 30 years, Nacalai Tesque has provided innovative and reliable chromatography products to the analytical industry. Our goal has not changed since the first COSMOSIL HPLC column was developed in 1980; it is to focus on customer needs and make customer's separations as successful as possible by offering the highest level of products and technical support. Keeping this goal in mind, we look forward to serving the industry for the next 30 years and beyond.



# CONTENTS

CORPORATE PROFILE .....	1
CONTENTS .....	2-3
GENERAL INFORMATION .....	4

## I. COSMOSIL HPLC COLUMN

1. COSMOSIL/COSMOGEL .....	6
2. Column selection guide .....	7
3. COSMOSIL silica packing materials .....	8-9
4. Performance guarantee .....	10-11
1) Quality guarantee of packing materials .....	10-11
2) Quality guarantee of COSMOSIL packed columns .....	11
5. Reversed phase chromatography .....	12-35
1) Octadecyl types .....	12-19
C <sub>18</sub> -MS-II .....	14-15
C <sub>18</sub> -AR-II .....	16-17
C <sub>18</sub> -PAQ .....	18-19
2) Special columns .....	20-29
Cholester .....	21-23
πNAP .....	24-25
PYE .....	26-27
NPE .....	28
PBB-R .....	29
3) Other columns .....	30-32
CN-MS .....	31
C <sub>22</sub> -AR-II, C <sub>8</sub> -MS, C <sub>4</sub> -MS, TMS-MS, PE-MS .....	32
4) Silica based preparative columns .....	33
15C <sub>18</sub> -MS-II, 15C <sub>18</sub> -AR-II, 15C <sub>18</sub> -PAQ .....	33
5) Ultra-Fast Liquid Chromatography (UFLC) .....	34-35
2.5C <sub>18</sub> -MS-II, 2.5Cholester .....	34-35
6. Normal phase chromatography .....	36-37
SL-II .....	36-37
7. Hydrophilic interaction chromatography .....	38-39
HILIC .....	38-39
8. Saccharide analysis .....	40-43
Sugar-D .....	41-43
NH <sub>2</sub> -MS .....	43
9. Protein separation wide pore columns .....	44-53
1) Reversed phase chromatography .....	44-47
Protein-R .....	44-45
C <sub>18</sub> -AR-300, C <sub>8</sub> -AR-300, C <sub>4</sub> -AR-300, Ph-AR-300 .....	46-47
2) Gel filtration chromatography .....	48-49
Diol-120-II, Diol-300-II .....	48-49
3) Ion-exchange chromatography .....	50-52
DEAE, QA, CM, SP .....	50-52
4) Hydrophobic interaction chromatography .....	53
HIC .....	53
10. Special columns for fullerenes .....	54-59
Buckyprep .....	56
Buckyprep-M .....	57
PBB .....	58
PYE .....	59
NPE .....	59
11. Special columns for carbon nanotubes .....	60
CNT-300, CNT-1000, CNT-2000 .....	60
12. Conventional columns versus high performance columns .....	61

## II. LIQUID CHROMATOGRAPHY RELATED PRODUCT

1. Packing materials for column chromatography .....	64-71
C <sub>18</sub> -OPN .....	65-67
C <sub>18</sub> -PREP .....	68
SL-II-PREP .....	69
Silica Gel (spherical • neutral) .....	70
Silica Gel (for column chromatograph) .....	71

2. Liquid chromatography related products (for mobile phase) .....	72
Ion-pair reagents .....	72-73
3. Liquid chromatography related products (for pretreatment) .....	74-75
Cosmonice filter .....	74-75
Cosmospin filter .....	74-75
COSMOSIL HPLC accessories .....	76

## III. APPLICATION DATA

1. COSMOSIL Chromatogram Index .....	78
2. Application data of substances in Japanese Pharmacopoeia, 15 <sup>th</sup> version .....	79
3. COSMOSIL Application Data .....	80-149
1) Drugs .....	80-112
2) Crude Drugs .....	113-117
3) Natural Compounds .....	118-120
4) Pesticides .....	121-122
5) Food Additives .....	123-124
6) Vitamins .....	125-127
7) Metabolites .....	128-129
8) Carbohydrates .....	130-132
9) Lipids .....	133-134
10) Nucleic Acid Related Substances .....	135
11) Amino Acids, Peptides and Proteins .....	136-142
12) The others .....	143-149
4. Reference list .....	150-151

## IV. TECHNICAL NOTE

1. Selectivity of packing materials in reversed phase liquid chromatography .....	154-160
2. Preparation of mobile phase for HPLC .....	161-163
3. Sample pretreatment for HPLC .....	164-167
4. Baseline noise in gradient elution .....	168
5. Troubleshooting for increased pressure .....	169-171
6. Effect of guard column .....	172
7. Troubleshooting for normal phase chromatography .....	173-174
8. Inner diameter of column (scale down and scale up) .....	175

## V. INDEX

1. Sample name .....	178-187
2. Packing material's name .....	188-202



# COSMOSIL HPLC COLUMN

## General information

### General ordering information

When placing an order with us or making an enquiry, please contact our International Business Development Group or your local distributor. Please clearly identify the product in question when submitting your enquiry.

The speed of innovation is accelerating. We always have brand new or improved columns not listed here. There are also many other products Nacalai Tesque can supply. Therefore we urge you to make enquiries.

### Product description and availability

Please note that the product specifications are subjected to changes and the manufacturing of some product may be stopped. Please consult the table on page 11 for cross-reference information on old products and their newer and better equivalents.

### Column identification

At the end of each section, the COSMOSIL and COSMOGEL packed columns are listed in a way that the particle size, stationary phase, column size of the packing material can be easily determined.

#### 38019-81 COSMOSIL 5 C<sub>18</sub>-MS-II 4.6 mm I.D. x 150 mm

(1)        (2)        (3)        (4)        (5)

When placing an order, please clearly indicate the product number (1), product name (2), particle size (3), type of stationary phase (4) and column size (5).

### Warranty claims

The manufacturer will replace defective columns if notified within 2 weeks of receipt of the product by the customer under the following conditions:

- 1) Column abnormalities are due to accidents in shipping or rough handling.
- 2) The number of effective plates of the column is considerably lower than the minimum guaranteed theoretical plate number documented in the inspection report that accompanies each column.

Please contact the International Business Development Section of Nacalai Tesque or your local distributor for additional information.

### Terms and conditions of sale

Terms are subject to conditions set forth by the authorized Nacalai Tesque distributors in each country.

### Not for clinical use

Nacalai Tesque products are not intended for clinical use. While clinical applications may be shown, these products are not validated for clinical use.

1. COSMOSIL/COSMOGEL .....	6
2. Column selection guide .....	7
3. COSMOSIL silica packing materials .....	8
4. Performance guarantee .....	10
1) Quality guarantee of packing materials .....	10
2) Quality guarantee of COSMOSIL packed columns .....	11
5. Reversed phase chromatography .....	12
1) Octadecyl types .....	12
C <sub>18</sub> -MS-II .....	14
C <sub>18</sub> -AR-II .....	16
C <sub>18</sub> -PAQ .....	18
2) Special columns .....	20
Cholester .....	21
πNAP .....	24
PYE .....	26
NPE .....	28
PBB-R .....	29
3) Other columns .....	30
CN-MS .....	31
C <sub>22</sub> -AR-II, C <sub>8</sub> -MS, C <sub>4</sub> -MS, TMS-MS, PE-MS .....	32
4) Silica based preparative columns .....	33
15C <sub>18</sub> -MS-II, 15C <sub>18</sub> -AR-II, 15C <sub>18</sub> -PAQ .....	33
5) Ultra-Fast Liquid Chromatography (UFLC) .....	34
2.5C <sub>18</sub> -MS-II, 2.5Cholester .....	34
6. Normal phase chromatography .....	36
SL-II .....	36
7. Hydrophilic interaction chromatography .....	38
HILIC .....	38
8. Saccharide analysis .....	40
Sugar-D .....	41
NH <sub>2</sub> -MS .....	43
9. Protein separation wide pore columns .....	44
1) Reversed phase chromatography .....	44
Protein-R .....	44
C <sub>18</sub> -AR-300, C <sub>8</sub> -AR-300, C <sub>4</sub> -AR-300, Ph-AR-300 .....	46
2) Gel filtration chromatography .....	48
Diol-120-II, Diol-300-II .....	48
3) Ion-exchange chromatography .....	50
DEAE, QA, CM, SP .....	50
4) Hydrophobic interaction chromatography .....	53
HIC .....	53
10. Special columns for fullerenes .....	54
Buckyprep .....	56
Buckyprep-M .....	57
PBB .....	58
PYE .....	59
NPE .....	59
11. Special columns for carbon nanotubes .....	60
CNT-300, CNT-1000, CNT-2000 .....	60
12. Conventional columns versus high performance columns .....	61



# 1. COSMOSIL/COSMOGEL


## General description of the COSMOSIL/COSMOGEL packing materials

Sample	Separation mode	Packing material	Stationary phase	Special features and applications	USP category	Page	
Organic compounds (low M.W.)	Reversed phase	C <sub>18</sub> -MS-II	Octadecyl group	Multi-purpose C <sub>18</sub> column. Monofunctional silylation on ultra-pure silica gel for separation of the widest range of compounds	L1	14	
		C <sub>18</sub> -AR-II		Multi-purpose C <sub>18</sub> column using ultrapure silica gel. Features strong acid resistance and suitable for a wide range of separation	L1	16	
		C <sub>18</sub> -PAQ		Reversed phase column, compatible with 100% water based mobile phases	L1	18	
		Cholester	Cholesteryl group	Usable under the same condition as C <sub>18</sub> . Unique rigid cholesteryl structure improves separation		21	
		πNAP	Naphthylethyl group	Stronger π-π interaction than phenyl column		24	
		PYE	Pyrenylethyl group	The most powerful π-π interaction		26	
		NPE	Nitrophenylethyl group	Separation utilizing π-π interaction and Dipole-dipole interaction		28	
		PBB-R	Pentabromobenzyl group	Separation utilizing dispersion force		29	
		CN-MS	Cyanopropyl group	Great reproducibility using isocratic elution mode	L10	31	
		C <sub>22</sub> -AR-II	Docosyl group	Alkyl chain columns except C <sub>18</sub> column			32
		C <sub>8</sub> -MS	Octyl group		L7		
		C <sub>4</sub> -MS	Butyl group		L26		
		TMS-MS	Trimethyl group		L13		
		PE-MS	Phenylethyl group	π-π interaction	L11		
			Normal phase	SL-II	**	Normal Phase chromatography with nonpolar organic solvents	L3
	Hydrophilic interaction	HILIC	Triazole	Retains highly polar compounds that would not be retained in C <sub>18</sub> column		38	
Monosaccharides Oligosaccharides	Hydrophilic interaction	Sugar-D	Secondary/Tertiary amine	A novel stationary phase for saccharide separation. Extended column life and increased stability. Alternative to aminopropyl type		41	
		NH <sub>2</sub> -MS	Aminopropyl group	Primary amine bonded column		43	
Proteins	Reversed phase	Protein-R	Octadecyl group	The most suitable reverse phased column for proteins		44	
		C <sub>18</sub> -AR-300	Octadecyl group	Wide pore type reversed phase columns with high acid resistance recommended for the separation of proteins, polypeptides, nucleic acids and other large molecules	L1	46	
		C <sub>8</sub> -AR-300	Octyl group		L7		
		C <sub>4</sub> -AR-300	Butyl group		L26		
		Ph-AR-300	Phenyl group		L11		
	Gel permeation	Diol-120-II	Diol group	Silica-based gel filtration column for high speed separation of proteins and water soluble polymer	L20	48	
		Diol-300-II					
	Ion-exchange	DEAE	Diethylaminoethyl type	Weak anion-exchange		50	
		QA	Quarternary ammonium type	Strong anion-exchange	L23		
		CM	Carboxymethyl type	Weak cation-exchange			
SP		Sulfopropyl type	Strong cation-exchange				
	Hydrophobic interaction	HIC	**	Hydrophobic interaction chromatography column for protein separation		53	
Fullerenes	**	Buckyprep	Pyrenylpropyl group	Standard column for fullerenes separation		56	
		Buckyprep-M	Phenothiazinyl group	Designed to separate metallofullerenes		57	
		PBB	Pentabromobenzyl group	Designed for the preparative separation of fullerenes using carbon disulphide, o-dichlorobenzene and toluene		58	
		NPE	Nitrophenylethyl group	Separation of derivatized fullerenes		59	
		PYE	Pyrenylethyl group	Separation of fullerenes and structural isomer			
Carbon nanotubes	Gel permeation	CNT-300 <b>NEW</b>	Hydrophilic group (neutral)	Separation of soluble carbon nanotubes		60	
		CNT-1000 <b>NEW</b>					
		CNT-2000 <b>NEW</b>					

For old type columns, please refer to page 61.

# 2. Column selection guide

## COSMOSIL columns selection guide



## COSMOSIL columns selection guide

### Organic compounds (low M.W.)

Octadecyl group bonded column (C<sub>18</sub>, ODS) are recommended as first-choice columns for separations of organic compounds (low M.W.). If there is not enough separation or no retention using COSMOSIL C<sub>18</sub> columns, COSMOSIL series offer many kinds of specialty columns.

- Medicines
- Crude drugs
- Natural compounds
- Pesticides
- Food additives
- Vitamins
- Lipids etc.

### Saccharides

- COSMOSIL Sugar-D is recommended for the separation of monosaccharides and oligosaccharides as a first-choice column.
- For the separation of sugar derivatives, COSMOSIL C<sub>18</sub>-PAQ is suitable as well.

### Proteins

- Please select based on the separation mode. Please refer to page 44.

### Fullerenes

- COSMOSIL Buckyprep is most suitable for the separation of fullerenes.

# 3. COSMOSIL silica packing material

## Introduction

Superior HPLC columns can be produced only with excellent packing materials and superb packing technique. COSMOSIL columns are well known for their high efficiency and high-resolution separations. Based on spherical, totally porous silica, COSMOSIL columns provide enhanced chemical and mechanical stability as well as very high surface coverage.

The selection of the C<sub>18</sub> chemistries available enables the chromatographer to tailor separation to special applications. The ultra pure silica based MS-II series with widely extended pH range are developed for improved separation of basic compounds. The C<sub>18</sub> AR-II phase provides increased acid resistance. Four unique bonded chemistries are available for COSMOSIL specialty columns: Cholester, πNAP, PYE and HILIC. These specialty columns may improve the separation compared with conventional columns. Five highly effective phases for fullerene separation are also available: Buckyprep, Buckyprep-M, PBB, PYE and NPE. COSMOGEL packing materials are non-silica based and provide superior in ion exchange columns.






Figure. Microscopic photograph of the silica gel

## Packing material and view showing a frame format of column



## Raw material silica gel

COSMOSIL is based on ultra pure porous spherical silica gel (purity: 99.99% or higher). Low-purity silica gel contains metal impurity which may cause interference in the separation, especially for metal coordination compounds.



### Metal coordination compounds

The compounds, which have 2 or more hydroxyl groups, carboxyl groups, carbonyl groups and/or amino groups in the adjacent positions, can complex with metal impurity, which results in peak tailing.




Table. Metal content percentage (ppm)

Packing material	Metal content (ppm)						Total
	Al	Ca	Fe	Mg	Na	Ti	
Ultra pure silica gel	4	10	11	3	12	<1	41
Normal silica gel	32	450	26	88	56	134	786

## Stationary phase construction

While C<sub>18</sub> columns are most widely used in reversed phase HPLC, it is important to distinguish between two very different bonded phase formats. Monomeric type C<sub>18</sub> format incorporates the bonding of the C<sub>18</sub> alkyl chain to a single silica atom on the silica gel backbone. Monomeric type columns such as the COSMOSIL C<sub>18</sub>-MS-II and the MS series have excellent synthesis reproducibility, very good lot-to-lot reproducibility and short mobile phase equilibration times. On the other hand, the polymeric C<sub>18</sub> format incorporates a tri-functional silylation procedure whereby the octadecyl group is bonded to 2 or 3 silica atoms on the silica gel backbone. This increases silylation results in far greater column stability particularly in acidic mobile phase conditions. Stereo recognition capability is also greater than that of the monofunctional silylation type C<sub>18</sub> columns. The polymeric format is offered in the AR-II and the entire AR-300 series of COSMOSIL columns. Please refer to product descriptions and application chromatograms for selection guidance.




Figure. Diagrams of different stationary phase constructions (before end-capping treatment)


## End-capping treatment

The silanols (Si-OH groups) on the silica surface provided bonding site for stationary phases. However, part of the silanol groups remain uncapped even after the end-capping treatment, they cause peak tailing for basic compounds. COSMOSIL packing materials for reversed phase chromatography are of near-perfectly end-capped residual silanol groups.

- R-NH<sup>+</sup> (Basic compounds)
- Si-O<sup>-</sup> (Free silanol group)

Basic compound can form ionic bonds with residual free silanols. The ionic bonding causes peak tailing of basic compounds if a silica-based column is not perfectly end-capped.

### Effect of residual silanol



## Synthesis reproducibility

By using strictly selected silica gel and constant synthesis conditions, the chemically bonded type column retains a variance of the capacity factor (k') between synthetic lots of within ±10% and a variance of the separation factor (α) of within ± 5%. The figures below show in graphic form the lot inspection results of synthesized packing material (COSMOSIL 5C<sub>18</sub>-MS-II). Figure 1 shows the variance of stationary phase (octadecyl group) introduced volume which is the basic indicator of the quality of the packing material. Figure 2 shows the end-capping efficiency of the packing material. The variance among the lots is reduced to the minimum in the COSMOSIL packed columns.





Figure 1. Variance of the combining volume between silica gel and C<sub>18</sub>

Figure 2. Variance of end capping efficiency of the packing material

# 4. Performance guarantee

## 1) Quality guarantee of packing materials

The strict quality control system of Nacalai Tesque supports the customers with an individual "Inspection Report" which accompanies each and every COSMOSIL and COSMOGEL Packed Column (except guard columns) and an additional "Certificate of Analysis" for the COSMOSIL 5C<sub>18</sub>-MS-II and 5C<sub>18</sub>-AR-II (4.6 mm I.D. x 150 mm and 4.6 mm I.D. x 250 mm).


### Validated columns

Product name	Product number	Column size
COSMOSIL 5C <sub>18</sub> -MS-II	38019-81	4.6 mm I.D. x 150 mm
	38020-41	4.6 mm I.D. x 250 mm
COSMOSIL 5C <sub>18</sub> -AR-II	38144-31	4.6 mm I.D. x 150 mm
	38145-21	4.6 mm I.D. x 250 mm
COSMOSIL Cholester	05976-61	4.6 mm I.D. x 150 mm
	05977-51	4.6 mm I.D. x 250 mm
COSMOSIL HILIC	07056-51	4.6 mm I.D. x 150 mm
	07057-41	4.6 mm I.D. x 250 mm

### COSMOSIL certificate of analysis


Validate terms of the physical properties of the silica gel, the carbon content, polar selectivity, hydrophobicity, silanol capacity, steric selectivity, inactive degree to basic and chelating compounds.

Certificate of Analysis		MS- II
COSMOSIL 5C <sub>18</sub> -MS- II		GEL Lot No. 29
Base silicagel material	Specification	Results
Median Particle Size		
50% cum vol [μm]	4.3-4.6	4.5
Surface Area [m <sup>2</sup> /g]	320-350	321
Pore Volume [ml/g]	0.9-1.1	0.9
Median Pore Diameter [nm]	11.0-13.0	11.8
Carbon content [%]	15.5-17.5	15.6
Atomic Emission [ppm]		
Al	≤ 5	1.5
Fe	≤ 20	8.3
Ti	≤ 0.5	0.1
Na	≤ 20	1.7
Chromatographic Results		
TEST A		
α [k'(Methyl Benzoate)/k'(Benzene)]	≤ 0.80	0.73
α [k'(Toluene)/k'(Benzene)]	≥ 1.60	1.76
TEST B		
α [k'(Caffeine)/k'(Phenol)]	≤ 0.44	0.40
TEST C		
α [k'(N-Acetyl Procainamide) /k'(Benzyl Alcohol)]	≤ 0.60	0.6
α [N(N'-Acetyl Procainamide) /N(Benzyl Alcohol)]	≥ 0.30	0.38
TEST D		
α [k'(8-Quinolol) /k'(Propylbenzene)]	≤ 0.14	0.09
α [N(8-Quinolol) /N(Propylbenzene)]	≥ 0.15	0.17
Nacalai Tesque Inc. Kyoto, Japan		
Approved : Quality Control Dept.		
Name: <i>K. Iwaguchi</i>	Date: 2009.8.21	00-3




- 3 lot numbers
- 3 lot numbers are available.

Please contact us. (info.intl@nacalai.com)



## 2) Quality guarantee of COSMOSIL packed columns

Inspection report contains data of number of theoretical plates (N), peak asymmetry (s), capacity factor (k') and separation factor (α) and pressure.



Inspection report





# 5. Reversed phase chromatography

## 1) Octadecyl types

### Introduction

The reversed phase HPLC column is most commonly used because of the high theoretical plate number, excellent separation characteristics, reproducibility, affordable cost and ease of use. Columns packed with the octadecyl group bonded type silica gel (C<sub>18</sub>, ODS) are the most widely employed. We offer three types of octadecyl group bonded columns: COSMOSIL C<sub>18</sub>-MS-II, C<sub>18</sub>-AR-II and C<sub>18</sub>-PAQ, each of which has a different separation property.


### Material characteristics

Packing material	C <sub>18</sub> -MS-II	C <sub>18</sub> -AR-II	C <sub>18</sub> -PAQ
Silica gel	High purity porous spherical silica		
Average particle size	3, 5, 15 μm*	3, 5, 15 μm	5, 15 μm
Average pore size	approx. 120 Å		
Specific surface area	approx. 300 m <sup>2</sup> /g		
Stationary phase			
	Octadecyl group		
Bonding type	Monomeric	Polymeric	
Main interaction	Hydrophobic interaction		
End-capping treatment	Near-perfect treatment		
Carbon content	approx. 16%	approx. 17%	approx. 11%
pH range	2 ~ 10**	1.5 ~ 7.5**	2 ~ 7.5
Feature	This phase is recommended for most applications but particularly effective for basic compounds.	This phase is recommended for separations requiring acidic mobile phase conditions. It also shows superior molecular shape selectivity to monomeric type C <sub>18</sub> columns.	This phase is designed to offer superior retention of polar compounds and excellent reproducibility in highly aqueous mobile phases, even in 100% aqueous.

\*For 2.5C<sub>18</sub>-MS-II, please refer to page 34.

\*\*Optimum pH range of columns based on silica gel is between 2 and 7.5.

### Difference of separation property



### Column selection based on application data

We prepare following application data to help you select separation conditions.

#### • COSMOSIL Chromatogram Index

More than 6,100 single compound elution profiles with full chromatographic condition description are available. They are not only an incredible help for chromatographers, but also can be used as references in choosing conditions for similar compounds.

These data are available at our web site: <http://www.nacalai.com>

#### • COSMOSIL Application Data

COSMOSIL Application Data is now available on our website. The online version includes more than 1,000 application data using COSMOSIL columns. The online data are searchable by name of sample and column. If you have any questions regarding the application data or separations of compounds not listed here, please feel free to e-mail us at [info.intl@nacalai.com](mailto:info.intl@nacalai.com).


#### • Application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version (246 data)

We prepare data of drugs using three kinds of C<sub>18</sub> columns that are specified in HPLC analysis in application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version. The data are available at our web site.

<http://www.nacalai.co.jp/en/cosmosil/TheJP15.htm>, or type "COSMOSIL Japanese Pharmacopoeia" at a search site.

### Column selection by mobile phase

- If a mobile phase is determined, use the following chart to select an appropriate COSMOSIL column.
- Refer to application data above for choosing a mobile phase of new analysis.
- Adjustment of pH is required for dissociative compounds.
- Generally acidic mobile phase is suitable for acidic compounds, and neutral mobile phase is suitable for basic compounds.
- If you are not sure about the mobile phase, try C<sub>18</sub>-MS-II first.




# C<sub>18</sub>-MS-II

The COSMOSIL 5C<sub>18</sub>-MS-II is a well-balanced column with better basic performance such as sharper peaks of basic compounds and chelating compounds, and high theoretical plate number. It is the most popular HPLC column because we produce very consistent products and minimize variation from lot-to-lot. Furthermore we provide abundant application data with the column, so it will help you to choose an analysis condition of your sample.

## Separation property

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Acetophenone (0.05µg)  
 2; Methyl Benzoate (0.5µg)  
 3; Benzene (2.0µg)  
 4; Toluene (2.0µg)



Column	Hydrophobic interaction α(Toluene/Benzene)	Pressure (MPa)	Theoretical plates number(Toluene)
COSMOSIL 5C <sub>18</sub> -MS-II	1.96	8.3	14300
A company C <sub>18</sub>	1.99	13.0	16800
B company C <sub>18</sub>	1.94	8.0	14000
C company C <sub>18</sub>	1.69	11.2	5600
D company C <sub>18</sub>	1.84	10.5	14200


## Analysis of basic compounds and metal coordination compounds

The COSMOSIL 5C<sub>18</sub>-MS-II column, taking advantage of a new end-capping treatment, can replace the original COSMOSIL C<sub>18</sub> (ODS) column. A new end-capping treatment with polar groups for "shield effect" has significantly improved peak shape for basic compounds. Ultra pure silica gel with low trace-metal content is used for COSMOSIL columns; thus the columns provide excellent peak shapes for chelating compounds.

### Basic compounds

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample:  
 1; Procainamide (0.38µg)  
 2; N-Acetylprocainamide (0.25µg)  
 3; Benzylalcohol (5.63µg)



### Metal coordination compounds

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 20mmol/l Phosphoric Acid = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV240nm

Sample: Oxine-copper



## Validation

The strict quality control system of Nacalai Tesque supports the customers with an individual "Inspection Report" which accompanies each and every COSMOSIL and COSMOGEL Packed Column (except guard columns) and an additional "Certificate of Analysis" for the COSMOSIL 5C<sub>18</sub>-MS-II (4.6 mm I.D. x 150 mm and 4.6 mm I.D. x 250 mm). For more informations, please refer to page 10.


## Application data

### Parasympatholytic agents

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 0.2AUFS

Sample:  
 1; Atropine Sulfate  
 2; Scopolamine Hydrobromide




### Analgesic antipyretic drugs

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample:  
 1; Acetoaminophen  
 2; Antipyrine  
 3; 4-Aminoantipyrine  
 4; Phenacetin



## Ordering information

### Analytical / Preparative column (Particle size: 5 µm)

#### COSMOSIL 5C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
1.0× 50	02824-31	4.6×100	38018-91
1.0×150	02896-01	4.6×150	38019-81
2.0× 30	05876-71	4.6×150 3 lots set*	09397-73
2.0× 50	04355-21	4.6×250	38020-41
2.0×100	05597-31	6.0×150	38021-31
2.0×150	38025-91	6.0×250	38022-21
2.0×250	05761-61	10× 50	05789-21
3.0×100	05458-51	10×150	34355-91
3.0×150	34245-31	10×250	38023-11
3.0×250	34254-11	20×150	05091-41
4.6× 30	34341-61	20×250	38024-01
4.6× 50	38017-01	28×250	05760-71

#### COSMOSIL 5C<sub>18</sub>-MS-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38014-31
4.6×10 Cartridge**	38015-89
10×20	38016-11
20×20	05790-81
20×50	34371-71
28×50	34347-01

\* For 4.6×150 3 lots set, please refer to page 11.  
 \*\* 3 cartridges included, needs a holder.  
 Please refer to page 76.

### Preparative column (Particle size: 15 µm)

#### COSMOSIL 15C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	34525-61
50×250	05886-41
50×500	34531-71

#### COSMOSIL 15C<sub>18</sub>-MS-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	05885-51
50×50	34527-41

### Fast LC column (Particle size: 3 µm)

#### COSMOSIL 3C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number
2.0× 50	05514-01
4.6× 10	38065-71
4.6× 50	38066-61
4.6×100	38067-51

For more informations, please refer to page 33 for 15C<sub>18</sub>-MS-II and page 34 for 2.5C<sub>18</sub>-MS-II.  
 For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.

## C<sub>18</sub>-AR-II

The COSMOSIL 5C<sub>18</sub>-AR-II packed column features a polymeric type of C<sub>18</sub> reversed phase material. The column employs an ultra-pure silica gel low in metal impurities and has near-perfect end-capping. In addition, it has stronger acid resistance than the COSMOSIL 5C<sub>18</sub>-AR. The COSMOSIL 5C<sub>18</sub>-AR-II column is especially effective for the separation of chelating compounds as well as both acidic and basic compounds.

### Acid resistance

The acidic resistance of COSMOSIL 5C<sub>18</sub>-AR-II is much improved compared with commercially available monomeric type octadecyl stationary phases. It retains high performance even in case of acidic mobile phases commonly used to separate acidic compounds and peptides.




Figure.  
Decomposition test in 0.1% Trifluoroacetic acid solution at 60°C.  
Capacity factor(k') = Naphthalene,  
Mobile phase: 70% Methanol


### Application data

#### Peptides

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
B; 0.05%TFA-Acetonitrile  
B conc. 10→40% 20min Linear gradient  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV220nm

Sample: 1; Oxytocin (0.8μg)  
2; Angiotensin II(Human) (0.8μg)  
3; Angiotensin I(Human) (0.8μg)  
4; Substance P (0.8μg)




NACALAI TESQUE, INC  
AP-0351

#### Organic acids

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 20/80  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.5AUFS

Sample: 1; Gallic Acid (0.63μg)  
2; Protocatechuic Acid (0.63μg)  
3; Gentisic Acid (0.63μg)  
4; Phthalic Acid (0.63μg)




NACALAI TESQUE, INC  
AP-0159

#### Salicylic acid esters

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.2AUFS

Sample: 1; Methyl Salicylate (2.3μg)  
2; Ethyl Salicylate (2.6μg)  
3; Propyl Salicylate (2.3μg)  
4; Butyl Salicylate (2.6μg)




NACALAI TESQUE, INC  
AP-0165

#### Parabens

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
Column size: 4.6mm I.D.-150mm  
Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
Flow rate: 1.0 ml/min  
Temperature: 30°C  
Detection: UV254nm, 0.12AUFS


Sample: 1; Methyl *p*-Hydroxybenzoate (0.125μg)  
2; Ethyl *p*-Hydroxybenzoate (0.125μg)  
3; Propyl *p*-Hydroxybenzoate (0.125μg)  
4; Butyl *p*-Hydroxybenzoate (0.125μg)



NACALAI TESQUE, INC  
AP-0099

## LC/MS Application data

### Identification of herbal medicine constituents by LC/MS



### Ordering information

#### Analytical / Preparative column (Particle size: 5 μm)

##### COSMOSIL 5C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
1.0× 50	02955-21	4.6×100	38143-41
1.0×150	02951-61	4.6×150	38144-31
2.0× 30	05098-71	4.6×150 3 lots set*	09396-83
2.0× 50	34400-81	4.6×250	38145-21
2.0×100	34469-11	6.0×150	38146-11
2.0×150	37992-51	6.0×250	38147-01
2.0×250	05272-71	10× 50	05369-21
3.0×100	05791-71	10×150	34350-41
3.0×150	38028-61	10×250	38149-81
3.0×250	38029-51	20×150	34316-01
4.6× 30	05877-61	20×250	38150-41
4.6× 50	38142-51	28×250	34362-91

##### COSMOSIL 5C<sub>18</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38141-61
4.6×10 Cartridge**	38008-89
10×20	38148-91
20×20	34458-51
20×50	34479-81
28×50	34363-81

\* For 4.6×150 3 lots set, please refer to page 11.  
\*\* 3 cartridges included, needs a holder.  
Please refer to page 76.

#### Preparative column (Particle size: 15 μm)

##### COSMOSIL 15C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	37978-51
50×250	38058-71
50×500	05884-61

##### COSMOSIL 15C<sub>18</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	38030-11
50×50	38057-81

#### Fast LC column (Particle size: 3 μm)

##### COSMOSIL 3C<sub>18</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number
2.0× 50	05478-91
4.6× 10	38068-41
4.6× 50	38069-31
4.6×100	38070-91


For more information on 15C<sub>18</sub>-AR-II, please refer to page 33.

For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.

# C<sub>18</sub>-PAQ

The COSMOSIL 5C<sub>18</sub>-PAQ maintains stable retention time even in 100% aqueous mobile phases. The new polymeric linking style gives this column a strong acidic resistance so that it is compatible with mobile phases of acidic pH that can permanently damage conventional octadecyl stationary phases. This type is especially good for separation of hydrophilic compounds.


## Stable performance under 100% aqueous conditions




The figure shows the change of retention time for thymine with 100% aqueous mobile phase (20 mmol/l phosphate buffer pH7). The sample was analyzed 4 times (1 hour). Flow of mobile phase was then stopped for 1 hour. The sample was analyzed under the same condition again after 1 hour. The conventional C<sub>18</sub> column showed change of retention time, but COSMOSIL 5C<sub>18</sub>-PAQ maintained stable retention time.

## Application data


### • Organic acids




### • Nucleobases and Nucleosides




### • dNTPs




### • Catecholamines



### • Ascorbic acid



### • 2-Phenylethyl glycosides



## Ordering information

### • Analytical / Preparative column (Particle size: 5 µm)

#### COSMOSIL 5C<sub>18</sub>-PAQ Packed Column

Column size I.D. x length (mm)	Product number
1.0× 50	05792-61
1.0×150	05793-51
2.0× 30	05878-51
2.0× 50	05794-41
2.0×100	05470-71
2.0×150	34449-71
2.0×250	05795-31
3.0×100	05796-21
3.0×150	05797-11
3.0×250	05798-01
4.6× 30	05879-41
4.6× 50	34451-21

#### COSMOSIL 5C<sub>18</sub>-PAQ Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	02484-91
10×20	34457-61
20×20	05803-11
20×50	05804-01
28×50	34455-81

### • Preparative column (Particle size: 15 µm)

#### COSMOSIL 15C<sub>18</sub>-PAQ Packed Column

Column size I.D. x length (mm)	Product number
28×250	05888-21
50×250	05890-71
50×500	05891-61

#### COSMOSIL 15C<sub>18</sub>-PAQ Guard Column

Column size I.D. x length (mm)	Product number
28×50	05887-31
50×50	05889-11

For more information on 15C<sub>18</sub>-PAQ, please refer to page 33.  
For flow rate and device of semi-micro columns, or preparative columns, please refer to page 175.








## 2) Special columns

### Introduction


Reversed phase HPLC columns have been widely used because of their superior resolution, high theoretical plate number and ease of use. Since hydrophobic interaction is the dominant separation mechanism in reversed phase chromatography, conventional stationary phases such as C<sub>18</sub> and C<sub>8</sub> do not offer optimum selectivity for compounds with similar hydrophobicity. COSMOSIL offers a broad selection of columns with unique stationary phases for separation of these difficult analytes. These columns offer improved separation of structurally similar compounds that are difficult to analyze with a C<sub>18</sub> type column.

### Material characteristics

Packing material	Cholester	πNAP	PYE	NPE	PBB-R
Silica gel	High purity porous spherical silica				
Average particle size	5 μm*	5 μm			
Average pore size	approx. 120 Å				
Specific surface area	approx. 300 m <sup>2</sup> /g				
Stationary phase					
	Cholesteryl group	Naphtylethyl group	Pyrenylethyl group	Nitrophenylethyl group	Pentabromobenzyl group
Bonding type	Monomeric type				
Main interaction	-Hydrophobic interaction -Molecular shape selectivity	-Hydrophobic interaction -π-π interaction	-Hydrophobic interaction -π-π interaction -Dispersion force -Charge-transfer interaction	-Hydrophobic interaction -π-π interaction -Dipole-dipole interaction	-Hydrophobic interaction -Dispersion force
End-capping treatment	Near-perfect treatment				
Carbon content	approx. 20%	approx. 11%	approx. 18%	approx. 9%	approx. 8%
Feature	-Usable under condition the same as C <sub>18</sub> -High molecular sharp selectivity	-Stronger π-π interaction than phenyl column	-Strongest π-π interaction	-Dipole-dipole interaction	-Dispersion force interaction

\* For 2.5 μm particle size, please refer to page 34.

### Column selection guide




## Cholester

COSMOSIL Cholester is a reversed phase HPLC column with cholesteryl groups bonded silica packing material, which provides equivalent hydrophobicity like traditional alkyl group bonded silica packing materials (C<sub>18</sub>, C<sub>30</sub>). However, Cholester offers strong molecular shape selectivity for hydrophobic compounds to yield unique and reproducible separation patterns following the same analytical conditions used with other ODS columns.

### Hydrophobic interaction

Figure shows the comparison of hydrophobic interactions with competitive C<sub>18</sub> columns. Cholester provides the same hydrophobicity as alkyl group bonded types (C<sub>18</sub>, C<sub>30</sub>). It is not necessary to change the analytical conditions when replacing C<sub>18</sub> or C<sub>30</sub> columns with Cholester.



### Molecular shape selectivity

The stationary phase of Cholester has very rigid structures and can distinguish different molecular shapes. Cholester offers improved separation for structurally similar compounds that are difficult to analyze with alkyl group bonded materials (C<sub>18</sub> and C<sub>30</sub>). As in the following example Cholester retains planar Triphenylene longer than stereoscopic o-Terphenyl.






Fig. Comparison of molecular sharp selectivity




### Improvement in separation

COSMOSIL Cholester provides enhanced selectivity over traditional C<sub>18</sub> columns and offers greater performance in separating isomers or other closely related compounds. COSMOSIL Cholester is ideal for method development and serves as an excellent alternative to traditional C<sub>18</sub> columns. The figure below shows analytical data of chloroacetophenone isomers. These isomers are difficult to separate with C<sub>18</sub> and C<sub>30</sub>, but they are well resolved by COSMOSIL Cholester.



### Efficiency of preparative separation

The figure below shows the comparison of efficiency of preparative separation with a C<sub>18</sub> column. Both columns show good separation. However, sample loading capacity for preparative separations can be affected by a slight difference in separation ability. COSMOSIL Cholester can load 4 times of sample volume compared with C<sub>18</sub> columns.




### Application data

#### • Catechins

##### COSMOSIL Application Data


Column: 4.6mm I.D. x 150mm  
 Mobile phase: A: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 10/90  
 B: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 30/70  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1: Galocatechin(GC) (0.80mg/ml)  
 2: Caffeine (0.08mg/ml)  
 3: Epigallocatechin(EGC) (0.80mg/ml)  
 4: Catechin (C) (0.40mg/ml)  
 5: Epicatechin(EC) (0.40mg/ml)  
 6: Epigallocatechin gallate(EGCG) (0.20mg/ml)  
 7: Galocatechin gallate(GCG) (0.40mg/ml)  
 8: Epicatechin gallate(EGC) (0.20mg/ml)  
 9: Catechin gallate(CG) (0.20mg/ml)  
 Injection Vol. 1.0µl



#### • Flavones

##### COSMOSIL Application Data

Column: 4.6mm I.D. x 150mm  
 Mobile phase: A: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 20/80  
 B: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 70/30  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1: Fisetin (0.25µg)  
 2: Myricetin (0.20µg)  
 3: 7,8-Dihydroxyflavone (0.05µg)  
 4: Luteolin (0.20µg)  
 5: Quercetin (0.20µg)  
 6: 7-Hydroxyflavone (0.10µg)  
 7: Baicalein (0.05µg)  
 8: 6-Hydroxyflavone (0.05µg)  
 9: Flavone (0.05µg)  
 10: Chrysin (0.05µg)  
 11: 6-Methoxyflavone (0.05µg)  
 12: 3-Hydroxyflavone (0.25µg)  
 13: 5-Hydroxyflavone (0.05µg)



### Ordering information

#### • Analytical / Preparative column (Particle size: 5 µm)

##### COSMOSIL Cholester Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	05968-71
1.0×250	05969-61
2.0× 30	08565-51
2.0× 50	06352-91
2.0×100	06948-01
2.0×150	05971-11
2.0×250	05972-01
3.0×150	05973-91
3.0×250	05974-81

##### COSMOSIL Cholester Guard Column


Column size I.D. x length (mm)	Product number
4.6×10	05976-61
4.6×150 3 lots set*	07970-03
10×20	05977-51
20×20	05978-41
20×50	05979-31
20×100	08011-91
20×150	05980-91
20×250	05981-81
28×50	05982-71
28×250	05983-61
28×250	05985-41

\* For 4.6×150 3 lots set, please refer to page 11.  
 \* For 2.5 µm particle size, please refer to page 34

# πNAP

COSMOSIL πNAP is a reversed phase HPLC column with naphthylethyl group bonded silica packing material. The naphthylethyl group is composed of two fused aromatic rings and forms strong π-π interactions with unsaturated compounds. This column offers improved separation of compounds such as positional isomers that are difficult to analyze with alkyl group bonded materials.

## Comparison of π-π interactions




COSMOSIL πNAP shows stronger π-π interactions than phenyl columns.


Figure. Comparison of π-π interaction

## Application data


### • Tolunitriles




### • Adrenal Cortical Hormones



### • Berberines



### • Isoflavones



## Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

### COSMOSIL πNAP Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	08076-61
1.0×250	08077-51
2.0× 30	08566-41
2.0× 50	08567-31
2.0×100	08299-51
2.0×150	08078-41
2.0×250	08079-31
3.0×150	08080-91


### COSMOSIL πNAP Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	08082-71
10×20	08087-21
20×20	08090-61
20×50	08091-51
28×50	08094-21

# PYE


COSMOSIL PYE column is a reversed phase column with 2-(1-pyrenyl) ethyl groups bonded silica packing material. This column utilizes  $\pi$ - $\pi$  interactions originating from the planar pyrene ring structure to separate structural isomers.

## Comparison $\pi$ - $\pi$ interaction




## Application data


### • Sterols




### • Phthalonitriles



### • Tocopherols



### • Capsaicins



## Attention

1. Methanol is recommended as a mobile phase for COSMOSIL PYE column. Acetonitrile is not recommended because it has many  $\pi$  electrons and interferes  $\pi$ - $\pi$  interactions between a sample and the stationary phase.
2. The stationary phase of COSMOSIL PYE, nitrophenyl group, has a large UV absorption. When the stationary phase detaches from silica gel and elutes, even a slight quantity can be detected and causes baseline noise. In such a case, wash the column with tetrahydrofuran. Detachment of a small amount of the stationary phase does not deteriorate a column's separation ability.
3. COSMOSIL PYE column is not suitable for gradient analysis.

## Ordering information

• Analytical / Preparative column (Particle size: 5 µm)

COSMOSIL 5PYE Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	02851-71
2.0×150	38042-61
2.0×250	34450-31

COSMOSIL 5PYE Guard Column


Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×150	37837-91	4.6×10	37903-11
4.6×250	37989-11	10×20	38041-71
10×250	37996-11	20×20	05867-91
20×250	38044-41	20×50	34475-21

## NPE

COSMOSIL NPE column is a reversed phase column with nitrophenylethyl groups bonded silica packing material. This column provides unique retention characteristics, slightly different from the COSMOSIL PYE column, utilizing both dipole-dipole and  $\pi$ - $\pi$  interactions.

### Selectivity for dipole-dipole interactions

COSMOSIL NPE strongly retains 1,8-dinitronaphthalene because of the strong dipole formed by the two nitro groups positioned on the same side of naphthalene.



### Attention

1. Methanol is recommended as a mobile phase for COSMOSIL NPE column. Acetonitrile is not recommended because it has many  $\pi$  electrons and interferes  $\pi$ - $\pi$  interactions between a sample and the stationary phase.
2. The stationary phase of COSMOSIL NPE, nitrophenyl group, has a large UV absorption. When the stationary phase detaches from silica gel and elutes, even a slight quantity can be detected and causes baseline noise. In such a case, wash the column with tetrahydrofuran. Detachment of a small amount of the stationary phase does not deteriorate a column's separation ability.
3. COSMOSIL NPE column is not suitable for gradient analysis.

### Ordering information

- Analytical / Preparative column (Particle size: 5  $\mu$ m)

#### COSMOSIL 5NPE Packed Column

Column size I.D. x length (mm)	Product number
1.0x150	05897-01
2.0x150	34328-51
2.0x250	34379-91

#### COSMOSIL 5NPE Guard Column

Column size I.D. x length (mm)	Product number
4.6x10	37902-21
4.6x250	37990-71
10x20	38045-31
20x20	05868-81
20x50	05869-71


## PBB-R

COSMOSIL PBB-R is a reversed phase column with pentabromobenzyl groups bonded silica packing material. This column provides unique selectivity for structurally similar compounds utilizing the dispersion force interaction. The dispersion force interaction of COSMOSIL PBB-R makes it useful for separation of structural isomers differing only by a double bond.

### Application data


- Halogen exchange reaction products

COSMOSIL PBB-R strongly retains iodine atom which has a large dispersion force, than bromine atom. So it can separate halogen exchange reaction products that are difficult to analyze with C<sub>18</sub> column.



- Surfactant agents

C<sub>18</sub> column can not separate Triton X-100 mixture, because (-OCH<sub>2</sub>CH<sub>2</sub>-) group has little hydrophobicity. However, COSMOSIL PBB-R can separate them because it distinguishes difference in the dispersion force, which depends on its molecular weight.



### Ordering information

- Analytical / Preparative column (Particle size: 5  $\mu$ m)

#### COSMOSIL 5PBB-R Packed Column

Column size I.D. x length (mm)	Product number
1.0x150	05899-81
2.0x150	05900-31
2.0x250	05904-91

#### COSMOSIL 5PBB-R Guard Column

Column size I.D. x length (mm)	Product number
4.6x10	05697-21
4.6x250	05698-11
10x20	05699-01
20x20	05700-51
20x50	05700-51









### 3) Other columns

#### Introduction


COSMOSIL alkyl type columns also include stationary phases with packing materials of C<sub>22</sub> (Dococyl group), C<sub>8</sub> (Octyl group), C<sub>4</sub> (Butyl group) and TMS (Trimethyl group). The order of retention force by hydrophobicity of each packing material is C<sub>22</sub>>C<sub>8</sub>>C<sub>4</sub>>TMS. The columns having lower hydrophobicity than C<sub>18</sub> or C<sub>22</sub> are effective for separation of high hydrophobic compounds and compounds with big difference in hydrophobic. Hydrophobicity of C<sub>22</sub> is about the same as C<sub>18</sub>. However, stereoselectivity of C<sub>22</sub> is higher than C<sub>18</sub>, and so C<sub>22</sub> may provide better separation.

The COSMOSIL PE-MS columns (phenylethyl group) and the COSMOSIL CN-MS columns (cyanopropyl group) provide a secondary separation mode ( $\pi$ - $\pi$  interaction). These columns are recommended when the other alkyl chain stationary phases do not offer optimum selectivity for structurally similar compounds.

#### Material characteristics

Packing material	C <sub>22</sub> -AR-II	C <sub>8</sub> -MS	C <sub>4</sub> -MS	TMS-MS	PE-MS	CN-MS
Silica gel	High purity porous spherical silica					
Average particle size	5 $\mu$ m					
Average pore size	approx. 120 Å					
Specific surface area	approx. 300 m <sup>2</sup> /g					
Stationary phase						
	Dococyl group	Octyl group	Butyl group	Trimethyl group	Phenylethyl group	Cyanopropyl group
Bonding type	Polymeric	Monomeric				
Main interaction	Hydrophobic interaction				Hydrophobic interaction $\pi$ - $\pi$ interaction	
End-capping treatment	Near-perfect treatment					
Carbon content	approx. 19%	approx. 10%	approx. 7%	approx. 5%	approx. 10%	approx. 7%


#### Different of separation characteristic



## CN-MS


Gradient elution is commonly used for the samples containing both polar and non-polar compounds. However, gradient elution may cause reproducibility problem depending on the gradient mixer and pump, and need an equilibration time for each analysis. COSMOSIL 5CN-MS offers rapid analysis and great reproducibility using isocratic elution mode.

#### Rapid analysis of samples containing both polar and non-polar compounds



#### Application data

##### • Acetoaminophen



#### Ordering information

##### • Analytical / Preparative column (Particle size: 5 $\mu$ m)

###### COSMOSIL 5CN-MS Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	38233-61
4.6×100	38234-51
4.6×150	38235-41
4.6×250	38236-31

###### COSMOSIL 5CN-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38231-81
10×20	38232-71

# C<sub>22</sub>-AR-II, C<sub>8</sub>-MS, C<sub>4</sub>-MS, TMS-MS, PE-MS

## Ordering information

- Analytical / Preparative column (Particle size: 5 µm)

### COSMOSIL 5C<sub>22</sub>-AR-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	05848-41	6.0×150	05850-91
4.6×100	05849-31	6.0×250	05851-81
4.6×150	04598-51	10×250	04969-91
4.6×250	04599-41	20×250	05183-41

### COSMOSIL 5C<sub>22</sub>-AR-II Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	04881-21
10×20	05554-81

### COSMOSIL 5C<sub>8</sub>-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38153-11	6.0×150	38157-71
4.6×100	38154-01	6.0×250	38158-61
4.6×150	38155-91	10×250	38159-51
4.6×250	38156-81	20×250	38160-11

### COSMOSIL 5C<sub>8</sub>-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38151-31
10×20	38152-21

### COSMOSIL 5C<sub>4</sub>-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38163-81	6.0×150	38167-41
4.6×100	38164-71	6.0×250	38168-31
4.6×150	38165-61	10×250	38169-21
4.6×250	38166-51	20×250	38170-81

### COSMOSIL 5C<sub>4</sub>-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38161-01
10×20	38162-91

### COSMOSIL 5TMS-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38173-51	6.0×150	38177-11
4.6×100	38174-41	6.0×250	38178-01
4.6×150	38175-31	10×250	38179-91
4.6×250	38176-21	20×250	38180-51

### COSMOSIL 5TMS-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38171-71
10×20	38172-61

### COSMOSIL 5PE-MS Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6× 50	38183-21	6.0×150	38187-81
4.6×100	38184-11	6.0×250	38188-71
4.6×150	38185-01	10×250	38189-61
4.6×250	38186-91	20×250	38190-21

### COSMOSIL 5PE-MS Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	38181-41
10×20	38182-31

## 4) Silica based preparative columns

# 15C<sub>18</sub>-MS-II, 15C<sub>18</sub>-AR-II, 15C<sub>18</sub>-PAQ

COSMOSIL series is available in 10 mm I.D. and 20 mm I.D. for semi-preparative applications and 28 mm I.D. and 50 mm I.D. for preparative scales. For column sizes and packing materials not listed below, contact either your local distributor or the manufacturer directly.


## Material characteristics

Packing material	15C <sub>18</sub> -MS-II	15C <sub>18</sub> -AR-II	15C <sub>18</sub> -PAQ
Silica gel	High purity porous spherical silica		
Average particle size	15 µm		
Average pore size	approx. 120 Å		
Specific surface area	approx. 300 m <sup>2</sup> /g		
Stationary phase	Octadecyl group (please refer to page 12.)		
Bonding type	Monomeric	Polymeric	
Main interaction	Hydrophobic interaction		
End-capping treatment	Near-perfect treatment		
pH range	2 ~ 10	1.5 ~ 7.5	2 ~ 7.5
Carbon content	approx. 16%	approx. 17%	approx. 11%
Feature	This phase is recommended for most of applications but particularly effective for basic organic compounds.	This phase is recommended for the separations requiring acidic mobile phase conditions. It also shows superior molecular shape selectivity to monomeric type C <sub>18</sub> columns.	This phase is designed to offer superior retention of polar compounds and excellent reproducibility in highly aqueous mobile phases, even in 100% aqueous.

## Application data

- Preparative separation using 50 mm I.D. column

Carbazole is extracted from anthracene oil (coal tar) and required high purity because it is often used for analytical applications. Following is the preparative separation of carbazole using a 50 mm I.D. COSMOSIL 15C<sub>18</sub>-AR-II.



Please refer to TECHNICAL NOTE 8, Inner diameter of column (scale down and scale up) at page 175.

## Ordering information

Please refer to page 15 for 15C<sub>18</sub>-MS-II, page 17 for 15C<sub>18</sub>-AR-II and page 19 for 15C<sub>18</sub>-PAQ.



## 5) Ultra-Fast Liquid Chromatography (UFLC)

2.5C<sub>18</sub>-MS-II, 2.5Cholester

Ultra-Fast Liquid Chromatography (UFLC) columns filled with sub-2  $\mu\text{m}$  particles have become widely used. Smaller particle columns generate higher pressures, while provide higher theoretical plate number. Therefore the column length should be shorter, and lead to lower resolution. COSMOSIL 2.5Cholester columns provide enhanced selectivity over traditional C<sub>18</sub> materials, and greater performance in separating isomers or other closely related compounds.


## Material characteristics

Packing material	C <sub>18</sub> -MS-II	Cholester
Silica gel	High purity porous spherical silica	
Average particle size	2.5 $\mu\text{m}$	
Average pore size	approx. 130 Å	
Specific surface area	approx. 330 m <sup>2</sup> /g	
Stationary phase	Octadecyl group	Cholesteryl group
Bonding type	Monomeric	
Main interaction	Hydrophobic interaction	Hydrophobic interaction Molecular shape selectivity
End-capping treatment	Near-perfect treatment	
Carbon content	approx. 18%	approx. 21%
Features	•First choice of reversed phase column	•The same mobile phase as C <sub>18</sub> •High molecular shape selectivity

## Ultra-Fast Liquid Chromatography


Very fast and efficient separation can be achieved using 2.5  $\mu\text{m}$  particles.

Note: Ultra high pressure liquid chromatography system or some modification of HPLC system is required for UFLC analysis.



## Comparison of analytical pressure


COSMOSIL 2.5Cholester and 2.5MS-II are reversed phase columns with 2.5  $\mu\text{m}$  particle size of silica packing material. It can be used under lower pressure than competitor's 2  $\mu\text{m}$  columns.




## Application data

## • Saikosaponins

COSMOSIL 2.5Cholester offers improved resolution for compounds difficult to analyze with C<sub>18</sub> without changing analytical condition.



## • Catechins



## Ordering information

• Analytical column (Particle size: 2.5  $\mu\text{m}$ )COSMOSIL 2.5C<sub>18</sub>-MS-II Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
2.0× 50	08994-31	3.0× 50	08997-01
2.0× 75	08995-21	3.0× 75	08998-91
2.0×100	08996-11	3.0×100	08999-81

## COSMOSIL 2.5Cholester Packed Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
2.0× 50	09000-01	3.0× 50	09049-91
2.0× 75	09047-11	3.0× 75	09050-51
2.0×100	09048-01	3.0×100	09051-41

# 6. Normal phase chromatography


## SL-II

Ultra-pure silica gel of more than 99.99% purity is used for the COSMOSIL SL-II packed column series. This column provides improved separation and reproducibility for compounds with carbonyl or phenol hydroxyl groups, which are often problematic to separate using conventional silica gel columns.

### Material characteristics


Packing material	SL-II
Silica gel	High purity porous spherical silica
Average particle size	3,5,15 μm
Average pore size	approx. 120 Å
Specific surface area	approx. 300 m <sup>2</sup> /g
Feature	•High purity silica gel (>99.99%) with specially treatment •Suitable for preparative separation

### Comparison with old type




### Comparison with medium-pressure column

COSMOSIL SL-II offers sharper peak compared with packing materials for medium-pressure liquid chromatography and open chromatography.



### Scaling up from analytical to preparative separation



Please refer to TECHNICAL NOTE 8, Inner diameter of column (scale down and scale up) at page 175

### Ordering information

- Analytical / Preparative column (Particle size: 5 μm)

#### COSMOSIL 5SL-II Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37999-81
4.6×100	38000-01
4.6×150	38001-91
4.6×250	38002-81

#### COSMOSIL 5SL-II Guard Column

Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×10	37997-01	6.0×150	38003-71
10×20	37998-91	6.0×250	38004-61
20×20	05874-91	10×250	38005-51
20×50	05875-81	20×250	38006-41
28×50	34359-51	28×250	34358-61

- Preparative column (Particle size : 15 μm)

#### COSMOSIL 15SL-II Packed Column

Column size I.D. x length (mm)	Product number
28×250	05893-41
50×250	05895-21
50×500	05896-11

#### COSMOSIL 15SL-II Guard Column

Column size I.D. x length (mm)	Product number
28×50	05892-51
50×50	05894-31

- Columns of high number of theoretical plates (Particle size : 3 μm)

#### COSMOSIL 3SL-II Packed Column


Column size I.D. x length (mm)	Product number
4.6× 10	38059-61
4.6× 50	38060-21
4.6×100	38061-11

# 7. Hydrophilic interaction chromatography

## HILIC


COSMOSIL HILIC is a new column for hydrophilic interaction chromatography with Triazole bonded silica packing material. The hydrophilic interaction chromatography is a variation of normal phase chromatography where a polar stationary phase is used with a mobile phase which contains a high concentration of water miscible organic solvent and a low concentration of aqueous eluent. The main retention mechanism is the partitioning of the polar analytes between the polar stationary and the non-polar mobile phase. As it is also called "aqueous normal phase", the elution order is similar to that of normal phase and the sample elution is in the order of increasing hydrophilicity. Without using ion-pair reagent COSMOSIL HILIC retains highly polar analytes that would not be retained in reversed phase chromatography. It also shows a weak anion-exchange mechanism with the positively charged stationary phase, thus acidic compound is strongly retained.

### Material characteristics

Packing material	HILIC
Silica gel	High purity porous spherical silica
Average particle size	5 µm
Average pore size	approx. 120 Å
Specific surface area	approx. 300 m <sup>2</sup> /g
Stationary phase	 Triazole
Interactions	Hydrophilic interaction, anion exchange
Object substance	Hydrophilic compounds, acidic compounds
Feature	Suitable for non-retaining by C <sub>18</sub>


### Comparison with C<sub>18</sub>

COSMOSIL HILIC can separate glycine and glycyglycine without ion-pair reagent. Although C<sub>18</sub> column can separate them with ion-pair reagents, there are some disadvantages such as column equilibration, preparation of mobile phase and column deterioration.




### Application data

#### Melamine related compounds





#### BCAA (amino acid branched-chain)



### Optimizing analytical conditions

COSMOSIL HILIC Chromatogram Index, which includes 154 chromatograms using COSMOSIL HILIC, is now available online at <http://www.nacalai.co.jp/en/cosmosil/>. This index is useful for optimizing analytical conditions for hydrophilic interaction chromatography.



### Ordering information

#### Analytical / Preparative column (Particle size: 5 µm)

##### COSMOSIL HILIC Packed Column

Column size I.D. x length (mm)	Product number
1.0×150	07869-11
1.0×250	07870-71
2.0× 30	08568-21
2.0× 50	07052-91
2.0×100	08569-11
2.0×150	07054-71
2.0×250	07489-91
3.0×150	07871-61
3.0×250	07872-51

Column size I.D. x length (mm)	Product number
4.6×150	07056-51
4.6×150 3 lots set*	09385-23
4.6×250	07057-41
10×250	07059-21
20×250	07060-81
28×250	07875-21

\* for 4.6×150 3 lots set, please refer to page 11.

##### COSMOSIL HILIC Guard Column


Column size I.D. x length (mm)	Product number
4.6×10	07055-61
10×20	07058-31
20×20	07854-91
20×50	07873-41
28×50	07874-31

# 8. Saccharide analysis

## Introduction


Saccharides are not retained on standard C<sub>18</sub> columns because of the low hydrophobicity of compounds. COSMOSIL Sugar-D and NH<sub>2</sub>-MS are specifically designed for separation of saccharides. COSMOSIL C<sub>18</sub>-PAQ is recommended for hydrophobic glycosides or saccharide derivatives.

## Material characteristics

Packing material	Sugar-D	NH <sub>2</sub> -MS
Silica gel	High purity spherical silica	
Average particle size	5 μm	
Average pore size	-	approx. 120 Å
Specific surface area	-	approx. 300 m <sup>2</sup> /g
Stationary phase	-	
Bonding type	Secondary/Tertiary amine	Aminopropyl
Object substance	Monosaccharides, oligosaccharides	
End capping treatment	-	Near-perfect treatment
Carbon content	-	approx. 4 %
Feature	<ul style="list-style-type: none"> <li>First choice of saccharide analysis</li> <li>High durability</li> <li>Good quantitative analysis</li> </ul>	<ul style="list-style-type: none"> <li>Different selectivity from Sugar-D</li> </ul>

## Comparison of retention


The conventional aminopropyl column is slightly more retentive than Sugar-D. The retention time of Sugar-D can be adjusted by increasing the concentration of acetonitrile in the mobile phase by 5%-10%.



# Sugar-D

Conventionally aminopropyl bonded stationary phases are used for liquid chromatographic analysis of mono- and oligosaccharides. General shortcomings of the conventional aminopropyl bonded phases are tailing and adsorption of certain saccharides and the general low durability (short active life) of these columns. These problems are addressed and solved by the novel COSMOSIL Sugar-D, resulting in better (sharper) separation and much improved durability.

## Durability




Decomposition condition  
 Solution: Water  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature

Column: 4.6 mm I.D. x 250 mm  
 Mobile phase: Acetonitrile : Water = 70 : 30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI  
 Sample: Maltose

## Adsorption characteristics

Certain types of saccharides such as arabinose or galactose are partially or temporarily adsorbed on conventional aminopropyl stationary phases causing tailing or no elution at all. COSMOSIL Sugar-D provides superior separation and high recovery for these saccharides.




## Application data

### • Mono- and Oligosaccharides

#### COSMOSIL Application Data

Column: Sugar-D  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1: Rhamnose (10μg)  
 2: Fructose (10μg)  
 3: Glucose (10μg)  
 4: Sucrose (10μg)  
 5: Maltose (10μg)  
 6: Raffinose (10μg)




### • Polyols

#### COSMOSIL Application Data

Column: Sugar-D  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1: Glycerol (10μg)  
 2: meso-Erythritol (10μg)  
 3: Xylitol (10μg)  
 4: Glucitol (10μg)  
 5: Maltitol (10μg)  
 6: Inositol (10μg)







# 9. Protein separation wide pore columns

## 1) Reversed phase chromatography

### Protein-R


COSMOSIL Protein-R is a reversed phased HPLC column designed specifically for peptide and protein separations.

#### Material characteristics


Packing material	Protein-R
Silica gel	High purity spherical porous silica
Average particle size	5 µm
Average pore size	approx. 300 Å
Specific surface area	approx. 150 m <sup>2</sup> /g
Stationary phase	Octadecyl group
Bonding type	Polymeric
Main Interactions	Hydrophobic interaction
End capping treatment	Near-perfect treatment
Feature	•High recovery rate •Acid-resistant

#### Comparison of separation

Protein-R shows sharper peaks for proteins than conventional C<sub>18</sub> wide pore columns.




#### Recovery rate



The figure left shows recovery rates for proteins using different columns. Protein-R shows a higher recovery rate than C<sub>4</sub>-300 and a much higher recovery rate than C<sub>18</sub>-300.

#### Durability


The figure below shows durability against acidic mobile phase of various columns. Protein-R shows a higher acid durability than C<sub>4</sub>-300.



Degradation test with 0.1%-Trifluoroacetic Acid at 60°C  
 (K): Napthalene in the mobile phase (Methanol : Water = 50 : 50 )

#### Application of peptide separation

5C<sub>18</sub>-MS-II (pore size 120 Å) shows better separation of low molecular weight proteins, but Protein-R shows better separation of high-molecular weight proteins.



#### Ordering information

• Analytical / Preparative column (Particle size: 5 µm)

COSMOSIL Protein-R Packed Column

Column size I.D. x length (mm)	Product number
2.0×150	06514-71
4.6× 50	06525-31
4.6×150	06526-21
4.6×250	06527-11





COSMOSIL Protein-R Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	06518-31
10×20	06528-01
20×20	08692-81


# C<sub>18</sub>-AR-300, C<sub>8</sub>-AR-300, C<sub>4</sub>-AR-300, Ph-AR-300

COSMOSIL offers a variety of stationary phases with wide-pore silica gel material for separations of polypeptides and proteins.

## Material characteristics

Packing material	5C <sub>18</sub> -AR-300	5C <sub>8</sub> -AR-300	5C <sub>4</sub> -AR-300	5Ph-AR-300
Silica gel	High purity spherical porous silica			
Average particle size	5 μm			
Average pore size	approx. 300 Å			
Specific surface area	approx. 150 m <sup>2</sup> /g			
Stationary phase				
	Octadecyl group	Octyl group	Butyl group	Phenyl group
Bonding type	Polymeric			
Main Interaction	Hydrophobic interaction			Hydrophobic interaction π-π interaction
End capping treatment	Near-perfect treatment			
Carbon content	approx. 12%	approx. 7%	approx. 6%	approx. 7%

## Comparison of separation



## Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

### COSMOSIL 5C<sub>18</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37911-01
4.6×150	37913-81
4.6×250	37914-71

Column size I.D. x length (mm)	Product number
10×150	37917-41
10×250	37918-31
20×150	37919-21
20×250	37920-81

### COSMOSIL 5C<sub>18</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37910-11
10×20	37965-11

### COSMOSIL 5C<sub>8</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37951-81
4.6×150	37953-61
4.6×250	37954-51

Column size I.D. x length (mm)	Product number
10×150	34345-21
10×250	34247-11
20×150	05861-51
20×250	34364-71

### COSMOSIL 5C<sub>8</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37950-91
10×20	34464-61

### COSMOSIL 5C<sub>4</sub>-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37956-31
4.6×150	37958-11
4.6×250	37959-01

Column size I.D. x length (mm)	Product number
10×150	34249-91
10×250	38047-11
20×150	34477-01
20×250	38048-01

### COSMOSIL 5C<sub>4</sub>-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37955-41
10×20	05862-41

### COSMOSIL 5Ph-AR-300 Packed Column

Column size I.D. x length (mm)	Product number
4.6× 50	37961-51
4.6×150	37963-31
4.6×250	37964-21

Column size I.D. x length (mm)	Product number
10×150	05865-11
10×250	34267-51
20×150	05866-01
20×250	34468-21

### COSMOSIL 5Ph-AR-300 Guard Column

Column size I.D. x length (mm)	Product number
4.6×10	37960-61
10×20	34268-41

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index



## 2) Gel filtration chromatography


# Diol-120-II, Diol-300-II

COSMOSIL Diol-120-II and Diol-300-II gel filtration columns are ideal for the size-based separation of proteins and other water soluble polymers. The separation MW range is 5,000 - 700,000 daltons for proteins and 300 - 300,000 daltons for water soluble polymers when Diol-120-II and Diol-300-II are used in series.

### Material characteristics

Packing material	5Diol-120-II	5Diol-300-II
Silica gel	High purity spherical porous silica	
Average particle size	5 µm	
Average pore size	approx. 120 Å	approx. 300 Å
Specific surface area	Diol group	
Object substance	Proteins, water soluble polymers	
Flow rate	0.5-1.0 (ml/min)	
Selection of pore size (Protein)	5,000-100,000 Da	10,000-700,000 Da
Selection of pore size (Water soluble polyme)	300-30,000 Da	500-300,000 Da


### Calibration curve of proteins



Column: COSMOSIL 5Diol-II 7.5mm I.D. x 600mm  
 Mobile phase: 20mmol/l Phosphate Buffer (pH 7.0) + 100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0ml/min  
 Temperature: 30°C


Sample	M.W.	Sample	M.W.
Thyroglobulin	660,000	Peroxidase	40,000
Catalase	250,000	Carbonic Anhydrase	30,000
Glucose Oxidase	160,000	α-Chymotrypsinogen A	25,700
Uricase	128,000	α-Chymotrypsin	25,200
Choline Oxidase	95,000	Trypsinogen	24,000
Transferrin	85,000	Trypsin (bovine)	23,300
Conalbumin	77,500	Myoglobin	17,000
Malate Dehydrogenase	70,000	Lysozyme	14,300
α-Glucosidase	68,500	Ribonuclease A	13,700
Albumin (BSA)	66,000	Cytochrome C	12,400
α-Amylase	52,500	Aprotinin	6,500
Fetuin	48,000	Gly-Gly	132
Albumin (Ovalbumin)	45,000		

### Comparison of separation property




### Application data


#### E.Coli Proteins




#### Water Soluble Polystyrene



#### Polyethylene Glycol



#### Pullulan



### Ordering information

#### Analytical column (Particle size: 5 µm)

COSMOSIL 5Diol-120-II Packed Column		COSMOSIL 5Diol-120-II Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5x300	38050-51	7.5x50	38049-91
7.5x600	38051-41		

COSMOSIL 5Diol-300-II Packed Column		COSMOSIL 5Diol-300-II Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5x300	38053-21	7.5x50	38052-31
7.5x600	38054-11		

### 3) Ion-exchange chromatography

## DEAE, QA, CM, SP

The packing materials for COSMOGEL ion-exchange glass columns are based on hydrophilic polymethacrylate 10 µm particles with a 1000 Å pore size. COSMOGEL packed columns are available with DEAE, a weak anion exchanger; QA, a strong anion exchanger; CM, a weak cation exchanger; and SP, a strong cation exchanger. The availability of four different ion exchangers provides chromatographers with the flexibility of column selection based on charge differences of samples.

#### Material characteristics

Packing material	DEAE	QA	CM	SP
Type	Diethylaminoethyl type Weak anion exchange	Quaternary ammonium type Strong anion exchange	Carboxymethyl type Weak cation exchange	Sulfopropyl type Strong cation exchange
Gel	Totally porous spherical hydrophilic polymer			
Average particle size	10 µm			
Average pore size	approx. 1,000 Å			
Functional group	N <sup>+</sup> H(C <sub>2</sub> H <sub>5</sub> ) <sub>2</sub>	N <sup>+</sup> H(CH <sub>3</sub> ) <sub>3</sub>	COO <sup>-</sup>	SO <sub>3</sub> <sup>-</sup>
Counter ion	Cl <sup>-</sup>	Cl <sup>-</sup>	Na <sup>+</sup>	Na <sup>+</sup>
Capacity (meq/g)	0.6	0.4	0.3	0.4
50% Ionization pH	10.8	11.0	5.7	2.6
pH range	<11	Whole area	>4	Whole area
Flow rate (appropriate)	7.5 mm I.D. / 8.0 mm I.D. ; 0.5 ~ 1.0 ml/min		20 mm I.D. ; 4.0 ~ 6.0 ml/min	
Flow rate (maximum)	7.5 mm I.D. / 8.0 mm I.D. ; 1.5 ml/min		20 mm I.D. ; 8.0 ml/min	
Pressure (maximum)	1.5 MPa			
Temperature	10 ~ 50°C			

#### Collection rate of proteins

##### Anion exchange columns

Protein	Collection rate(%)	
	DEAE	QA
Catalase	93	92
BSA	98	97
Ovalbumin	100	103
β-Lactoglobulin	95	102
Myoglobin	97	101
Transferrin	105	102

Sample: 0.02 M Tris-HCl buffer (pH 8.2)  
Elution: Tris-HCl buffer (pH 8.2) + 0.5 M NaCl

##### Cation exchange columns



Pretein	Collection rate(%)	
	CM	SP
γ-Globulin	97	96
α-Chymotrypsinogen	96	103
Myoglobin	99	95
Cytochrome C	98	94

Sample: 0.02 M Tris-HCl buffer (pH 6.5)  
Elution: Tris-HCl buffer (pH 6.5) + 0.5 M NaCl

### Difference between anion exchange columns and cation exchange columns

- Separation of proteins on anion exchange columns


The higher the negative charge, the longer the sample is retained on an anion exchange column. As shown below, with a weak alkaline mobile phase, the lower the isoelectric point, the longer the sample is retained.



### Sample loading capacity and resolution

Up to 1 mg of sample can be well separated on an 8 mm I.D. column.

If the resolution is high enough, 10 mg of sample can be separated.



## Selection of the mobile phase

Generally, anion-exchange columns are operated with the mobile phase pH at least one point higher than the isoelectric point (pI) of samples, while cation-exchange columns are operated with the mobile phase pH at least one point lower than the pI. The elution force of bivalent ions such as Bis-tris HCl is stronger than univalent ions such as Tris HCl.

Table. Buffer type and pH

Anion exchange (DEAE, QA)	pH	Cation exchange (CM, SP)
	4.0	Formic acid buffer
Piperazine buffer	5.0	Acetic acid buffer
Bis-Tris buffer	6.5	Phosphoric acid buffer
Tris buffer	8.0	HEPES buffer
Monoethanolamine buffer	9.5	

The initial mobile Phase (A) is 20-50 mmol/l of one of the above mentioned buffer solutions and the final mobile phase (B) is the mobile phase (A) with an addition of 20-600 mmol/l of salt.

## Selection of salts

High concentration of salts, generally NaCl, is used in elution buffers. When stronger elution buffer is needed, CaCl<sub>2</sub> or MgCl<sub>2</sub> can be used for DEAE and QA columns, and Na<sub>2</sub>SO<sub>4</sub> can be used for CM and SP columns.

## Selection of organic solvents

Water miscible solvents can be used when the elution is not strong enough. COSMOGEL columns can be used with up to 20% of water miscible solvents such as acetonitrile and 2-propanol.

## Ordering information

- Analytical column (Particle size: 10 µm)

### COSMOGEL Stainless Packed Column

Product name	Column size I.D. x length (mm)	Product number
COSMOGEL DEAE Stainless Packed Column	7.5×75	43371-91
COSMOGEL QA Stainless Packed Column		43373-71
COSMOGEL CM Stainless Packed Column		43375-51
COSMOGEL SP Stainless Packed Column		43377-31

### COSMOGEL Glass Packed Column

Product name	Column size I.D. x length (mm)	Product number
COSMOGEL DEAE Glass Packed Column	8.0×75	37845-81
COSMOGEL QA Glass Packed Column		37846-71
COSMOGEL CM Glass Packed Column		37844-91
COSMOGEL SP Glass Packed Column		37847-61

## 4) Hydrophobic interaction chromatography

### HIC

COSMOSIL 5HIC is designed for one step desalting and separation of proteins. Hydrophobic interaction chromatography (HIC) is an effective method for purification and separation of proteins (especially enzymes) based on differences in their surface hydrophobicity. Since this method does not use organic solvents like reversed phase chromatography, there is only a little loss in enzyme activity and the tertiary structure of proteins.


## Material characteristics

Packing material	HIC
Silica gel	High purity spherical porous silica
Average particle size	5 µm
Average pore size	approx. 300Å
Specific surface area	approx. 150m <sup>2</sup> /g
Main interaction	Hydrophobic interaction


## Application data

A buffer with high salt concentration, usually 1-2 mol/l of (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, is used as an initial mobile phase for adsorption of samples to a weakly hydrophobic stationary phase. The elution is done with a decreasing salt gradient.

- Separation of protein standards



- Separation of β-Glucosidase



## Ordering information

- Analytical column (Particle size: 5 µm)

### COSMOSIL 5HIC Packed Column






Column size I.D. x length (mm)	Product number
4.6×50	04263-21

# 10. Special column for fullerenes


## Introduction

Separation of fullerenes, especially preparative-scale separation, on conventional HPLC columns is always problematic due to the low solubility and small recovery of fullerenes. COSMOSIL offers a variety of columns designed for preparative-scale separation of fullerenes including higher fullerenes, metallofullerenes and fullerene derivatives.

## Material characteristics

Packing material	Buckyprep	Buckyprep-M	PBB	PYE	NPE
Silica gel	High purity porous spherical silica				
Average particle size	5 μm				
Average pore size	approx. 120Å				
Specific surface area	approx. 300 m <sup>2</sup> /g				
Stationary phase					
	Pyrenylpropyl group	Phenothiazinyl group	Pentabromobenzyl group	Pyrenylethyl group	Nitrophenylethyl group
Bonding type	Monomeric				
End capping treatment	Near-perfect	None	Near-perfect treatment		
Carbon content	approx. 17%	approx. 13%	approx. 8%	approx. 18%	approx. 9%
Feature	Standard column for fullerene separation.	Designed to separate metallofullerenes.	Designed for preparative separation of C <sub>60</sub> , C <sub>70</sub> .	Separation of fullerene and structural isomers.	Separation of fullerene derivatives

## Comparison of retention in toluene




## Solubility and boiling point of each solvent for C<sub>60</sub>

Solubility and boiling point of each solvent for C<sub>60</sub>

Solvent	mg/ml	b.p.(°C)
Methanol	0.001	64.5
Acetonitrile	0.018	81.8
n-Hexane	0.046	68.7
Toluene	3.2	111
Chlorobenzene*	7.0	132
Carbon disulfide	12	46.3
o-Dichlorobenzene*	27	180
1,2,4-Trichlorobenzene	21.3	213

\*: R.S.Ruoff, et al., J.Phys.Chem.,97,3379(1993)

## Comparison of separation property



## Suggested solvents

Solvent	Feature
Chlorobenzene	Stronger eluent than toluene. Recommended for higher fullerenes.
o-Dichlorobenzene	Stronger eluent than chlorobenzene.
1,2,4-Trichlorobenzene	Strongest eluent. It can be used as a washing solvent for higher fullerenes. To wash a column, inject 3 ml of 1,2,4-trichlorobenzene to a 4.6 mm I.D. × 250 mm column and 50 ml to a 20 mm I.D. × 250 mm column after every operation.
n-Hexane	Weak eluent. Recommended for weakly retained fullerenes.
Acetonitrile	Weak eluent. Recommended for weakly retained fullerenes.


Note : Use them after filtration or distillation, if they are not for HPLC.

# Buckyprep


COSMOSIL Buckyprep is a pyrenylpropyl group bonded silica based column specifically designed for fullerene separation. The unparalleled separation capabilities have enabled COSMOSIL Buckyprep to become the world benchmark of HPLC column for fullerene separation. COSMOSIL Buckyprep retains fullerenes very strongly with a mobile phase of 100% toluene and exceeds the injection volume of a standard C<sub>18</sub> column by a factor of 35. Therefore, preparative-scale separation can be obtained with a 250 mm × 4.6 mm I.D. analytical column.

## Application data

### Higher fullerenes




### Detivatized fullerene



## Fullerene Chromatogram Index

A comprehensive index with more than 100 chromatograms for fullerene separation is available from Nacalai Tesque, Inc. Please feel free to contact us for a copy.



## Ordering information

### Analytical / Preparative column (Particle size: 5 µm)



COSMOSIL Buckyprep Packed Column		COSMOSIL Buckyprep Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37977-61	4.6×10	37983-71
10×250	37981-91	10×20	37984-61
20×250	37982-81	20×50	34374-41
28×250	34346-11	28×50	05871-21

# Buckyprep-M

COSMOSIL Buckyprep-M is a phenothiazinyl group bonded silica based column specifically designed for metallofullerene separation. Metallofullerenes are retained more strongly than other fullerenes on this column. COSMOSIL Buckyprep-M is also effective for the separation of higher fullerenes and fullerene derivatives.

## Application data

### Metallo fullerenes



## Ordering information

### Analytical / Preparative column (Particle size: 5 µm)


COSMOSIL Buckyprep-M Packed Column		COSMOSIL Buckyprep-M Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	04138-71	4.6×10	04139-61
10×250	04141-11	10×20	04140-21
20×250	04142-01	20×50	34474-31
28×250	05873-01	28×50	05872-11



# PBB


COSMOSIL PBB is a pentabromobenzyl group bonded silica based column useful for preparative-scale separation of fullerenes. It can be used with *o*-Dichlorobenzene, which has greater solubility for fullerenes than toluene. The loading capacity of COSMOSIL PBB for C<sub>60</sub> and C<sub>70</sub> can be three times greater than COSMOSIL Buckyrep.

## Separation of fullerenes with different mobile phases



## Preparative-scale separation

The loading capacity of COSMOSIL PBB for C<sub>60</sub> and C<sub>70</sub> can be three times greater than COSMOSIL Buckyrep.



## Ordering information


• Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL 5PBB Packed Column		COSMOSIL 5PBB Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37980-01	4.6×10	37987-31
10×250	37985-51	10×20	37988-21
20×250	37986-41	20×50	34375-31

# PYE • NPE

## Application data

• PCBM, PCBB



## Ordering information

• Analytical / Preparative column (Particle size: 5 μm)

COSMOSIL 5PYE Packed Column		COSMOSIL 5PYE Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×250	37989-11	4.6×10	37903-11
10×250	37996-11	10×20	38041-71
20×250	38044-41	20×50	34475-21
28×250	34300-91		

COSMOSIL 5NPE Packed Column		COSMOSIL 5NPE Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
4.6×150	37902-21	4.6×10	37904-01
4.6×250	37990-71	10×20	38045-31
10×250	05469-11	20×50	05869-71
20×250	38046-21		

# 11. Special columns for carbon nanotubes

## CNT-300, CNT-1000, CNT-2000

COSMOSIL CNT series are ideal for separation of soluble carbon nanotubes based on sizes. COSMOSIL CNT series are packed with hydrophilic group-bonded silica packing material. The columns are specially designed to avoid adsorption of carbon nanotubes to silica support and thus ensure high resolution and maximum recovery of carbon nanotubes. COSMOSIL CNT series are available in three different pore sizes, 300 Å, 1000 Å and 2000 Å, respectively.


### Material characteristics

Packing material	CNT-300	CNT-1000	CNT-2000
Silica gel	High purity porous spherical silica		
Average particle size	5 µm		
Average pore size	approx. 300 Å	approx. 1000 Å	approx. 2000 Å
Stationary phase	Hydrophilic group (neutral)		
pH range	2.0-7.5		
Pressure	15 MPa and below		

### Application data

#### • Carbon nanotubes

COSMOSIL CNT columns offered improved separation for DNA wrapped carbon nanotubes by connecting three columns with different pore sizes.



### Ordering information

#### • Analytical column (Particle size: 5 µm)

COSMOSIL CNT-300 Packed Column		COSMOSIL CNT-300 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09195-71	7.5×50	09194-81

COSMOSIL CNT-1000 Packed Column		COSMOSIL CNT-1000 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09197-51	7.5×50	09196-61

COSMOSIL CNT-2000 Packed Column		COSMOSIL CNT-2000 Guard Column	
Column size I.D. x length (mm)	Product number	Column size I.D. x length (mm)	Product number
7.5×300	09199-31	7.5×50	09198-41

# 12. Conventional columns versus high performance columns

### Introduction

A period of more than 30 years has passed since the first COSMOSIL 5C<sub>18</sub> columns were developed and offered for sale. Continuous technical improvement has made many of these columns obsolete and not of the highest quality and performance available any more. However, many long-term users continue to employ these older conventional columns for routine analysis and quality control. Nevertheless, the manufacture of these older columns will eventually cease and we strongly urge customers to replace the conventional columns with their higher performance equivalents outlined in the table below. For additional information, contact the manufacturer or your local distributor directly.

Conventional columns (old)		High performance columns (new)
5C <sub>18</sub> -AR	→	5C <sub>18</sub> -AR-II
5C <sub>18</sub>	→	5C <sub>18</sub> -MS-II
5C <sub>18</sub> -MS	→	5C <sub>18</sub> -MS-II
5C <sub>18</sub> -P	→	5C <sub>18</sub> -PAQ
5C <sub>18</sub> -P-MS	→	5C <sub>18</sub> -PAQ
5C <sub>8</sub>	→	5C <sub>8</sub> -MS
5TMS	→	5TMS-MS
5PE	→	5PE-MS
5CN-R	→	5CN-MS
5NH <sub>2</sub>	→	5NH <sub>2</sub> -MS
5C <sub>18</sub> -300	→	5C <sub>18</sub> -AR-300
5C <sub>8</sub> -300	→	5C <sub>8</sub> -AR-300
5C <sub>4</sub> -300	→	5C <sub>4</sub> -AR-300
5SL	→	5SL-II

### Ordering information

Product name	Column size I.D. x length (mm)	Product number
COSMOSIL 5C <sub>18</sub> Packed Column	4.6×150	39047-81
	4.6×250	39265-21
COSMOSIL 5C <sub>18</sub> -MS Packed Column	4.6×150	37971-21
	4.6×250	37972-11
COSMOSIL 5C <sub>18</sub> -AR Packed Column	4.6×150	37861-61
	4.6×250	37862-51
COSMOSIL 5C <sub>18</sub> -P Packed Column	4.6×150	39103-31
	4.6×250	39280-11
COSMOSIL 5C <sub>18</sub> -P-MS Packed Column	4.6×150	37995-21
	4.6×250	37994-31
COSMOSIL 5C <sub>8</sub> Packed Column	4.6×150	39042-31
	4.6×250	39260-71
COSMOSIL 5TMS Packed Column	4.6×150	39057-51
	4.6×250	39275-91
COSMOSIL 5CN-R Packed Column	4.6×150	39114-91
	4.6×250	39285-61
COSMOSIL 5NH <sub>2</sub> Packed Column	4.6×150	39150-11
	4.6×250	39290-81
COSMOSIL 5C <sub>18</sub> -300 Packed Column	4.6×150	39607-41
	4.6×150	39037-11
COSMOSIL 5SL Packed Column	4.6×150	39037-11
	4.6×250	39255-51

For more information on other columns, please feel free to contact us.





# LIQUID CHROMATOGRAPHY RELATED PRODUCT

1. Packing materials for column chromatography .....	64
C <sub>18</sub> -OPN .....	65
C <sub>18</sub> -PREP .....	68
SL-II-PREP .....	69
Silica gel (spherical • neutral) .....	70
Silica gel (for column chromatograph) .....	71
2. Liquid chromatography related products (for mobile phase) .....	72
Ion-pair reagents .....	72
3. Liquid chromatography related products (for pretreatment) .....	74
Cosmonice filter .....	74
Cosmospin filter .....	74
COSMOSIL HPLC accessories .....	76



# 1. Packing materials for column chromatography

## Introduction


Open column chromatography is an excellent and easy technique for large-scale preparation and purification at low cost. COSMOSIL offers both normal and reversed phase packing materials based on totally porous spherical silica, which provides higher separation, less pressure and higher reproducibility than irregular silica.

## Material characteristics


Packing material	C <sub>18</sub> -OPN	C <sub>18</sub> -PREP	SL-II-PREP	Silica gel 60(neutral)
Silica gel	High purity porous spherical silica			
Average particle size	75 · 140 μm	40 · 75 · 140 μm	75 · 140 μm	
Average pore size	approx. 120 Å			approx. 60 Å
Specific surface area	approx. 300 m <sup>2</sup> /g			approx. 500 m <sup>2</sup> /g
Stationary phase	Octadecyl group		None	
Carbon content	—	approx. 19%	0%	
End-capping treatment	Treated	None	—	
Useful range	Open column chromatography / Flash column chromatography			
	Reversed phase chromatography		Normal phase chromatography	

For more informations on other silica gel, please refer to page 71.

## Selection guide (Reversed phase)



## Selection guide (Normal phase)



## C<sub>18</sub>-OPN

Conventional reversed phase C<sub>18</sub> packing materials are restricted to about 30-50% water in the mobile phase. The COSMOSIL C<sub>18</sub>-OPN is a new "Water-Wet" C<sub>18</sub> packing material developed for reversed phase open column chromatography. The C<sub>18</sub>-OPN material can be used in 100% aqueous effluents.

## Characteristic

The external surface of the C<sub>18</sub>-OPN gel is coated with hydrophilic group to increase wettability of the gel, and octadecyl group is bonded in the pore of the gel. This physical characteristic of the gel makes the reversed phase open column chromatography possible with 100 % water.




Figure 1. Structure of C<sub>18</sub>-OPN



Figure 2. Packing material in water

Left : C<sub>18</sub>OPN provides good resolution  
Right : C<sub>18</sub>PREP float up

## Application data

### • Separation of Theobromine and Theophylline




Figure shows the sample are clearly separated by reversed open column chromatography with 70% of water.

## Influence of particle size





Table. Comparison between 75 µm and 140 µm particle size silica.

## Flow rate




Since reversed phase chromatography generally employs high viscosity solvents such as water and methanol, the flow rate is lower than that of normal phase chromatography. The flow rate of reversed phase depends on the mobile phase composition. Figure indicates that the flow rate of the COSMOSIL 140C<sub>18</sub>-OPN (140 µm in particle size) is about 2.5 times higher than that of the COSMOSIL 75C<sub>18</sub>-OPN.

Figure. Concentration of methanol against flow rate  
 Column size: 10 mm I.D. x 180 mm bed height (gravitational liquid flow)

## Comparison of normal phase

- Separation of *p*-Cresol and *p*-Ethylphenol by normal and reversed phase mode.

Since the structural difference between *p*-Cresol and *p*-Ethylphenol is only one methylene group, it is difficult to separate such samples under normal phase condition. On the other hand, the samples are clearly separated under reversed phase condition with COSMOSIL C<sub>18</sub>-OPN packing material.



## Column size and required amount of packing material

Table. Column size and required amount of C<sub>18</sub>-OPN packing material

Column I.D. (mm)	Bed height (mm)	Amount of C <sub>18</sub> -OPN(g)
10	150	4
	250	7
20	150	17
	250	28
30	150	38
	250	63

## Reproducibility and washing methods

Wash the COSMOSIL C<sub>18</sub>-OPN packing material with tetrahydrofuran, chloroform or other solvents to remove the impurities. This packing material has excellent reproducibility and can be used repeatedly.

### "CAUTION"

Do not wash with basic solvents of pH 7 or more which will dissolve the silica gel or pH 2 or less which will cleave the C<sub>18</sub> stationary phase. Dry the packing material at 50°C or less. See end of this chapter for packing method.

## Ordering information

- COSMOSIL C<sub>18</sub>-OPN

Product name	Average particle size	Product number	PKG size
COSMOSIL 75C <sub>18</sub> -OPN	75 µm	37842-66	100 g
		37842-95	500 g
		37842-11	1 kg
COSMOSIL 140C <sub>18</sub> -OPN	140 µm	37878-16	100 g
		37878-45	500 g
		37878-61	1 kg

## C<sub>18</sub>-PREP

The large particle size C<sub>18</sub> bulk materials are widely used for lab to process scale purifications. COSMOSIL offers three different particle sizes of C<sub>18</sub> packing materials.

### Particle size, flow rate and theoretical plate number

Because reversed phase chromatography employs effluents of high viscosity such as methanol and water, the flow rate is lower than that of normal phase chromatography, which uses effluents of low viscosity such as hexane and ethyl acetate.




Figure 1. Flow rate against theoretical plate number  
Column size: 20 mm I.D. x 300 mm




Figure 2. Concentration of methanol against flow rate  
Column size: 10 mm I.D. x 180 mm bed height (gravitational liquid flow)


### Application data

#### • Vitamin E

##### COSMOSIL Application Data

Column: 40C<sub>18</sub>-PREP  
Column size: 20mm I.D.-300mm (Closed column)  
Mobile phase: Methanol  
Flow rate: 9.9 ml/min  
Temperature: Room temperature  
Detection: UV280nm

Sample: 1; DL- $\alpha$ -Tocopherol (5mg)  
2; DL- $\alpha$ -Tocopherol Acetate (5mg)




#### • Natural compounds

##### COSMOSIL Application Data

Column: 40C<sub>18</sub>-PREP  
Column size: 20mm I.D.-300mm (Closed column)  
Mobile phase: Methanol/ 0.05%TFA-H<sub>2</sub>O = 70/30  
Flow rate: 9.9 ml/min  
Temperature: Room temperature  
Detection: UV254nm

Sample: 1; Baicalin (40µg)  
2; Baicalein (120µg)  
3; Wogonin (40µg)



### Ordering information

#### • COSMOSIL C<sub>18</sub>-PREP

Product name	Average particle size	Product number	PKG size
COSMOSIL 40C <sub>18</sub> -PREP	40 µm	37932-86	100 g
		37932-15	500 g
		37932-31	1 kg
COSMOSIL 75C <sub>18</sub> -PREP	75 µm	37933-76	100 g
		37933-05	500 g
		37933-21	1 kg
COSMOSIL 140C <sub>18</sub> -PREP	140 µm	37934-66	100 g
		37934-95	500 g
		37934-11	1 kg

## SL-II-PREP

COSMOSIL SL-II-PREP is ultra pure silica gel packing material more than 99.99% purity. COSMOSIL SL-II-PREP provides improved separation and reproducibility for compounds with carbonyl or phenol hydroxyl groups, which are often problematic on conventional silica gel materials.

\*All chromatograms shown below are obtained with silica gel packed into stainless steel columns.

### Performance for chelating compounds


Highly purified silica gel of COSMOSIL SL-II-PREP enables separation of metal coordination compounds without adsorption.

#### • Metal Coordination Compounds

##### COSMOSIL Application Data

Column: 10mm I.D.-250mm  
Mobile phase: Hexane/Ethanol = 95/5  
Flow rate: 5.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Quinizarin  
2; p-Nitrobenzyl Alcohol




#### • Organic Acid and Amide

##### COSMOSIL Application Data

Column: 10mm I.D.-250mm  
Mobile phase: Hexane/Ethanol = 90/10  
Flow rate: 5.0 ml/min  
Temperature: 30°C  
Detection: UV254nm

Sample: 1; Salicylic Acid  
2; Salicylamide



### Ordering information

#### • COSMOSIL SL-II

Product name	Average particle size	Product number	PKG size
COSMOSIL 75SL-II-PREP	75 µm	38012-64	100 g
		38012-35	500 g
		38012-51	1 kg
COSMOSIL 140SL-II-PREP	140 µm	38013-54	100 g
		38013-25	500 g
		38013-41	1 kg



## Silica gel (spherical • neutral)

Since conventional silica gels are weakly acidic, some pH sensitive compounds may be decomposed during the purification by column chromatography with the acidic silica gels. The pH of Silica gel 60 (spherical • neutral) is adjusted to nearly neutral for the separation of not only pH sensitive compounds but also new compounds that the physical properties are still unknown.

\*All chromatograms shown below are obtained with silica gel packed into stainless and steel columns.

### Comparison with conventional silica gel

#### • Purification of Acetal



### Ordering information

#### • Silica gel (spherical • neutral)

Product name	Average particle size	Product number	PKG size
Silica gel 60 (spherical • neutral) for column chromatograph	75 µm	30511-64	100 g
		30511-35	500 g
		30511-51	1 kg
		30511-06	5 kg
		30511-22	25 kg
	140 µm	30518-94	100 g
		30518-65	500 g
		30518-81	1 kg
		30518-52	25 kg

## Silica gel (for column chromatograph)

### Ordering information

#### • Silica gel (spherical)

Product name	Particle size	Pore size	Grade	Product number	PKG size
Silica Gel 60, spherical	approx. 70 ~ 230 mesh	60 Å	SP	30731-71	1 kg
				30731-42	25 kg
	approx. 150 ~ 325 mesh		SP	30733-51	1 kg
				30733-22	25 kg
Silica Gel 120, spherical	approx. 70 ~ 230 mesh	120 Å	SP	30734-41	1 kg
	approx. 150 ~ 325 mesh			30735-31	1 kg

#### • Silica gel (irregular)

Product name	Particle size	Pore size	Grade	Product number	PKG size
Silica Gel 60	approx. 70 ~ 230 mesh	60 Å	SP	30724-55	500 g
				30724-71	1 kg
				30724-84	5 kg
				30724-42	25 kg
	approx. 230 ~ 400 mesh		SP	30721-85	500 g
				30721-01	1 kg
				30721-14	5 kg
				30721-72	25 kg



## 2. Liquid chromatography related products (for mobile phase)


### Ion-pair reagents

The use of ion pair reagents as mobile phase additives extends the applicability of reversed phase HPLC. Ionic or highly polar compounds are difficult to analyze by reversed phase using only organic solvent and buffer solution because of the short retention time. Ion pair reagents are strong hydrophobic ions which form neutral ion pairs with oppositely charged samples molecules, making the efficient ODS columns amenable to separate ionic or highly polar samples. Nacalai Tesque offers a broad range of ion pair reagents for pharmaceutical compounds and other highly polar materials.

#### General use of ion-pair reagents in the mobile phase

When using ion pair reagents, ample time should be allowed for establishing equilibrium and for cleaning the column.


When using ion pair reagents with an alkyl chain of C<sub>10</sub> or shorter, it typically takes 20 minutes for establishing equilibrium and 30 minutes for cleaning. It may take more than 1 hour to clean the column when using ion pair reagents with an alkyl chain longer than C<sub>10</sub>. Therefore, it is highly recommended to prepare a column for exclusive use with ion pair reagents.



\*For HILIC column, please refer to page 38.

#### Separation adjustment

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm  
 Sample: 1; L-Noradrenaline  
 2; L-Adrenaline  
 3; L-DOPA  
 4; Dopamine Hydrochloride




#### Application data

- Low-molecular-weight unsaturated carboxylic acids
- Amino acids

##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra-n-butylammonium Phosphate, 20mmol/l Phosphate buffer(pH7) = 10 /90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS


Sample: 1; Acrylic Acid (2.0µg)  
 2; Crotonic Acid (2.0µg)  
 3; Tiglic Acid (2.0µg)  
 4; Sorbic Acid (2.0µg)



##### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Sodium 1-Hexanesulfonate, 20mmol/l Phosphate buffer(pH3) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV215nm, 0.5AUFS

Sample: 1; Histidine (1.0µg)  
 2; Tyrosine (3.0µg)  
 3; Phenylalanine (1.0µg)  
 4; Tryptophan (0.5µg)



#### Ordering information

- for Basic samples

Product name	R:	Grade	Product number	PKG size
Sodium 1 - Butanesulfonate	C <sub>4</sub> H <sub>9</sub> -	SP	31331-94	5 g
Sodium 1 - Pentanesulfonate	C <sub>5</sub> H <sub>11</sub> -	SP	31730-64	5 g
			31730-22	25 g
Sodium 1 - Hexanesulfonate	C <sub>6</sub> H <sub>13</sub> -	SP	31529-24	5 g
			31529-82	25 g
Sodium 1 - Heptanesulfonate	C <sub>7</sub> H <sub>15</sub> -	SP	31528-34	5 g
			31528-92	25 g
Sodium 1 - Octanesulfonate	C <sub>8</sub> H <sub>17</sub> -	SP	31729-04	5 g
			31729-62	25 g
Sodium 1 - Nonanesulfonate	C <sub>9</sub> H <sub>19</sub> -	SP	31626-44	5 g
Sodium 1 - Decanesulfonate	C <sub>10</sub> H <sub>21</sub> -	SP	31429-34	5 g
Sodium 1 - Undecanesulfonate	C <sub>11</sub> H <sub>23</sub> -	SP	32030-04	5 g
Sodium 1 - Dodecane sulfonate	C <sub>12</sub> H <sub>25</sub> -	SP	31426-64	5 g
Sodium Lauryl Sulfate	**	SP	31623-32	25 g

#### 0.5M Solution

Sodium 1 - Butanesulfonate	C <sub>4</sub> H <sub>9</sub> -	SP	31332-84	5×10 ml
Sodium 1 - Hexanesulfonate	C <sub>6</sub> H <sub>13</sub> -	SP	31532-64	10 ml
			31532-06	5×10 ml
Sodium 1 - Octanesulfonate	C <sub>8</sub> H <sub>17</sub> -	SP	31733-34	10 ml
			31733-76	5×10 ml

- for Acid samples

Product name	X-:	Grade	Product number	PKG size
Tetra - n - butylammonium Bromide	-Br	SP	32824-72	25 g
Tetra - n - butylammonium Chloride	-Cl	EP	32935-51	1 g
			32935-64	5 g
			32935-22	25 g
Tetra - n - butylammonium Hydrogensulfate	-HSO <sub>4</sub>	GR	32924-62	25 g
Tetra - n - butylammonium Iodide	-I	SP	32905-54	5 g
			32905-12	25 g
Tetra - n - butylammonium Perchlorate	-ClO <sub>4</sub>	SP	32906-44	5 g
			32906-02	25 g
Tetra - n - butylammonium Phosphate	-H <sub>2</sub> PO <sub>4</sub>	SP	32929-54	5 g


#### 0.5M Solution

Tetra - n - butylammonium Phosphate	-H <sub>2</sub> PO <sub>4</sub>	SP	32926-26	10 ml
			32926-84	5×10 ml

### 3. Liquid chromatography related products (for pretreatment)

## Cosmonice filter

Injection of samples containing particulates (microparticles, precipitates, colloid substances) will clog HPLC columns, shorten injector life, and result in extensive maintenance on pumps. Cosmonice filters are used to remove particulates from samples and prolong the life of HPLC system components. There are two types of Cosmonice filters as stated below.



#### W series (Aqueous solution)

W series are installed with low adsorption hydrophilic durapore filter (polyvinylidenedifluoride, PVDF). W series can be used with both aqueous and organic solvents. They are best suited for prefiltration of protein and other biological samples.

#### S series (Organic solvents)

S series are installed Teflon filter (polytetrafluoroethylene, PTFE) with strong resistance to organic solvents, acids and alkalis. They are best suited for prefiltration of samples with aggressive organic solvents such as chloroform and tetrahydrofuran.

Please refer to TECHNICAL NOTE 3, Sample pretreatment for HPLC at page 164.

#### Ordering information

##### • Cosmonice filter

Product name	Diameter (mm)	Pore size (µm)	Process volume	Hold-up volume	Product number	PKG size
Cosmonice Filter W (Aqueous)	4	0.45	less than 1 ml	< 10 µl	06543-04	100 pkg
	13	0.45	0.5 - 10 ml	< 30 µl	06544-94	100 pkg
	25	0.45	3 - 50 ml	< 100 µl	06545-84	50 pkg
Cosmonice Filter S (Solvent)	4	0.45	less than 1 ml	< 10 µl	06541-24	100 pkg
	13	0.45	0.5 - 10 ml	< 30 µl	06542-14	100 pkg

[Connection] Inlet: luer-lock, Outlet: luer-slip, Connectable needles

## Cosmospin fileter

Cosmospin filters are used to remove fine particles and precipitates from samples by centrifugation. They utilize omnipore hydrophilic PTFE membrane filter, which has a wide range of chemical resistance. Cosmospin filters are the best choice for HPLC sample filtration. Two pore sizes, G (0.2 µm) and H (0.45 µm), are available.

Please refer to TECHNICAL NOTE 3, Sample pretreatment for HPLC at page 164.



#### Ordering information

##### • Cosmospin filter

Product name	Pore size (µm)	Maximum sample volume	Hold-up volume	Maximum centrifugal force	Rotor size (fixed-angle)	Filtration area	Color	Product number	PKG size
Cosmospin Filter G	0.2	0.4 ml	5 µl	5,000× g	1.5 ml	0.2 cm <sup>2</sup>	brown	06549-44	100 pkg
Cosmospin Filter H	0.45	0.4 ml	5 µl	5,000× g	1.5 ml	0.2 cm <sup>2</sup>	white	06540-34	100 pkg

Dimension: Diameter 10.6 mm x Length 45 mm Membrane: Omnipore Hydrophilic PTFE Sample reservoir and collection tube: Polypropylene

#### Chemical compatibility

Solvent	Cosmonice W series	Cosmonice S series	Cosmospin	Solvent	Cosmonice W series	Cosmonice S series	Cosmospin
Acetic acid, 98%	+	+	+	Hydrogen gas	+	+	+
Acetone	-	+	+	Hydrogen peroxide (3%)	+	+	
Acetonitrile	+	+	+	Hydraulic oil (5606)	+	+	+
Ammonia solution (6N)	+	+	+	Hypo (photo)	+	+	+
Ammonium hydroxide (conc.)	+	+	-	Isopropyl acetate	+	+	+
Amyl alcohol	+	+	+	Isopropyl alcohol	+	+	+
Benzene	+	+	-	Kerosene	+	+	+
Benzyl alcohol	+	+	-	Methanol	+	+	+
Boric acid	+		+	Methyl ethyl ketone	-	+	+
Butyl acetate		+		Methyl isobutyl ketone	+	+	-
Carbon tetrachloride	+	+	+	2-Methyl-1-propanol	+	+	+
Chloroform	+	+	+	Nitric acid (6N)	+	+	
Cyclohexanone	-	+	-	Nitrobenzene	+	+	-
Dichloromethane	+	+	-	Ozone gas	-	+	-
Dimethylacetamide	-	+	+	Paraldehyde		+	
Dimethylformamide	+	+	+	Pentane	+	+	-
Dimethylsulfoxide	-	+	-	Petroleum ether	+	+	
Dioxane	+	+	+	Phenol (water saturation)	+	+	-
DMSO	-	+	-	Phosphate buffer solution	+		+
Ethers	+	+	+	2-Propanol	+	+	+
Ethyl acetate	+	+	+	Pyridine	-	+	+
Ethyl alcohol	+	+	+	Seawater	+	+	+
Ethyl cello solve	+	+	+	Silicone oils	+	+	+
Ethylene glycol	+	+	+	Sodium hydroxide (conc.)	+	+	+
Formamide	+	+	+	Sulfuric acid (6N)		+	
Freon, TF or PCA solvent	+	+	+	Toluene	+	+	-
Gasoline	+	+	+	THF	-	+	-
Glycerine (Glycerol)	+	+	+	Trichloroacetic acid	+	+	+
Helium gas		+	+	Trichloroethane	+	+	-
Hexane	+	+	-	Trichloroethylene	+	+	-
Hydrochloride (6N)	+	+	+	TFA	+	+	-
Hydrofluoric acid	-	+	-	Xylene	+	+	+

+ : Recommended, - : Not recommended, (blank) : Not data available

# COSMOSIL HPLC accessories



## APPLICATION DATA

I. COSMOSIL HPLC column  
 II. Liquid chromatography related product  
 III. Application data  
 IV. Technical note  
 V. Index

### Ordering information


#### COSMOSIL Guard Cartridge Holder

Product name	Product number	PKG size
COSMOSIL Guard Cartridge Holder	38009-79	1 PKG



#### COSMOSIL Column Prefilter

Product name	Product number	PKG size
COSMOSIL Column Prefilter	39361-19	1 PKG



#### COSMOSIL Column Spare Filter for Prefilter

Product name	Product number	PKG size
COSMOSIL Column Spare Filter for Prefilter	39539-09	2 PKG




#### COSMOSIL Column Connecting Tube

Product name	Product number	PKG size
COSMOSIL Column Connecting Tube	37843-69	1 PKG



1. COSMOSIL Chromatogram Index .....	78
2. Application data of substances in Japanese Pharmacopoeia, 15 <sup>th</sup> version .....	79
3. COSMOSIL Application Data .....	80
1) Drugs .....	80
2) Crude Drugs .....	113
3) Natural Compounds .....	118
4) Pesticides .....	121
5) Food Additives .....	123
6) Vitamins .....	125
7) Metabolites .....	128
8) Carbohydrates .....	130
9) Lipids .....	133
10) Nucleic Acid Related Substances.....	135
11) Amino Acids, Peptides and Proteins.....	136
12) The others .....	143
4. Reference list .....	150



# 1. COSMOSIL Chromatogram Index

More than 6,100 single compound elution profiles with full chromatographic condition description are available. They are not only an incredible help for chromatographers, but also can be used as references in choosing conditions for similar compounds.

These data are available at our web site: <http://www.nacalai.com>

Category:

Column name:


Sample Name:  begins with

Molecular Formula:  contains (Keyword search)

CAS number:  (ex: 498-02-2)

Result/Page:

## Data example



# 2. Application data of substances in Japanese Pharmacopoeia, 15<sup>th</sup> version

We prepare data of drugs using three kinds of C<sub>18</sub> columns that are specified in HPLC analysis in Application Data of Substances in Japanese Pharmacopoeia, 15<sup>th</sup> version. The data are available at our web site.

<http://www.nacalai.co.jp/en/cosmosil/TheJP15.htm>, or type "Cosmosil Japanese Pharmacopoeia" at a search site.

## Interpretation of application data

THE JAPANESE PHARMACOPOEIA (Fifteenth Edition)

**1 Acetaminophen (Purity)**

Acetaminophen [103-90-2]  
 C<sub>9</sub>H<sub>9</sub>NO<sub>2</sub> 151.16  
*N*-(4-Hydroxyphenyl) acetamide

	5C <sub>18</sub> -MS-II	5C <sub>18</sub> -AR-II	5C <sub>18</sub> -PAQ
3 <sup>*</sup> Suitability	10.7	10.5	11.2
5 <sup>*</sup> Recommendation			recommendation

4<sup>\*</sup> HPLC chromatograms for 5C<sub>18</sub>-MS-II, 5C<sub>18</sub>-AR-II, and 5C<sub>18</sub>-PAQ columns.

6<sup>\*</sup> Condition and JP description

Column size: 4.6 mmID-150 mm  
 Mobile phase: Methanol : 0.05 mol/L potassium dihydrogenphosphate (pH 4.7)  
 = 20 : 80  
 Temperature: 40°C  
 Detection: UV 225 nm  
 Flow rate: 1.0 ml/min  
 Sample: 1: *p*-Aminophenol Hydrochloride (0.02 mg/ml) [Internal standard]  
 2: Acetaminophen (0.02 mg/ml) [Purity sample]  
 3: 4'-Acetoxycetamide (0.02 mg/ml) [Impurity]  
 Test solution: Mobile phase  
 Injection volume: 10 µL

Items described in The Japanese Pharmacopoeia  
 (Flow rate) Adjust the flow rate so that the retention time of Sample 2 is about 5 minutes.  
 (Selection of column) (4<sup>\*</sup>)  
 Use a column eluting off Sample 1 before Sample 2, and resolution (R<sub>s</sub>) between those 2 peaks should not be less than 7.  
 (Detection sensitivity)  
 Adjust the detection sensitivity so that the peak height of Sample 2 obtained from 10 µL of the standard solution is about 15% of the full scale.  
 (Time span of measurement)  
 About 5 times as long as the retention time of Sample 2 after the solvent peak.  
 (Point to notice)  
 \*2 0.05 mol/L potassium dihydrogenphosphate (pH 4.7)  
 Dissolve 0.50 g of potassium dihydrogenphosphate in 900 mL of water, adjust pH to exactly 4.7 with dilute sodium hydrochloride test solution, and add water to make 1000 mL solution.

NACALAI TESQUE, INC

- ① Substance name
- ② Substance information
- ③ Suitability (○ : suitable, × : unsuitable, \*\* : depend on condition)
- ④ HPLC chromatogram
- ⑤ Recommending column
- ⑥ Condition and JP description



### 3. COSMOSIL Application Data

1) Drugs	P80	(Application data of substances in Japanese Pharmacopoeia ...P80-109, P113-117)
2) Crude Drugs	P113	
3) Natural Compounds	P118	
4) Pesticides	P121	
5) Food Additives	P123	
6) Vitamins	P125	
7) Metabolites	P128	
8) Carbohydrates	P130	
9) Lipids	P133	
10) Nucleic Acid Related Substances	P135	
11) Amino Acids, Peptides and Proteins	P136	
12) The others	P143	


#### 1) Drugs

##### • Acrinol

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 4.6mmol/l Sodium *I*-Octanesulfonate, 65mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV268nm

Sample: 1; Acrinol (0.04µg)




AP-0457

##### • Azathioprine

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 25mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV296nm

Sample: 1; 6-Mercaptopurine Hydrate (0.16µg)  
 2; Azathioprine (0.16µg)  
 3; Benzoic Acid (9.6µg)




AP-0471

##### • Aztreonam

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Tetra-*n*-butylammonium Hydrogensulfate (pH3.0 with 0.5mol/l Na<sub>2</sub>HPO<sub>4</sub>) = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; 4-Aminobenzoic Acid (0.5µg)  
 2; Aztreonam (5.0µg)




AP-0474

##### • Acetaminophen

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH4.7 with NaOH) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV225nm

Sample: 1; *p*-Aminophenol Hydrochloride (0.2µg)  
 2; Acetaminophen (0.2µg)  
 3; 4'-Acetoxyacetanilide (0.2µg)



AP-0452


#### 1) Drugs

##### • Atenolol

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/Tetrahydrofuran/ 4.6mmol/l Sodium *I*-Octanesulfonate, 1.2mmol/l Tetra-*n*-butylammonium Hydrogensulfate, 25mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 9/1/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV226nm

Sample: 1; Atenolol (0.1µg)




AP-0468

##### • Amikacin Sulfate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH6.5 with KOH) = 72/28  
 Flow rate: 1.5 ml/min  
 Temperature: 35°C  
 Detection: UV340nm

Sample: 1; Amikacin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative (0.6µg)  
 2; Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative (0.6µg)




AP-0460

##### • Meglumine Sodium Amidotrizoate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 6.25mmol/l Tetra-*n*-butylammonium Phosphate, 50mmol/l K<sub>2</sub>HPO<sub>4</sub> (pH7.0 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Amidotrizoic Acid (0.25µg)  
 2; Acetizoic Acid (0.30µg)




AP-0739

##### • Amlodipine Besilate

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 30mmol/l KH<sub>2</sub>PO<sub>4</sub> = 65/35  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV237nm

Sample: 1; Amlodipine Besilate (0.56µg)  
 2; Isobutyl *p*-Hydroxybenzoate (0.60µg)




AP-1084

##### • Amoxicillin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 10mmol/l CH<sub>3</sub>COONa (pH4.5 with CH<sub>3</sub>COOH) = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Amoxicillin (3.0µg)




AP-0463

##### • Ampicillin

#### COSMOSIL Application Data

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 10% Acetonitrile/ 90% 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub> (pH5.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Ampicillin Sodium Salt (10µg)  
 2; Guaiacol Glycerol Ether (5µg)



AP-0464



## 1) Drugs

## ● Isoxsuprine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 32.6mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>,  
 18.4mmol/l Sodium *L*-Pentanesulfonate  
 (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV269nm

Sample: 1; Isoxsuprine Hydrochloride (2.0µg)




## ● Isoniazid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 10mmol/l-Sodium *L*-Tridecaneulfonate-  
 Methanol/ 50mmol/l Phosphate buffer  
 (pH2.5) =60/40  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV265nm

Sample: 1; Isonicotinic Acid (0.4µg)  
 2; Isoniazid (0.5µg)




## ● Idoxuridine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 13/87  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; 2'-Deoxyuridine (0.3µg)  
 2; 5-Iodouracil (1.2µg)  
 3; Idoxuridine (3.0µg)  
 4; Sulfathiazole (0.5µg)




## ● Idoxuridine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 4/96  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; 2'-Deoxyuridine (0.04µg)  
 2; 5-Iodouracil (0.12µg)  
 3; Idoxuridine (0.40µg)




## ● Ipratropium Bromides

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.5% H<sub>3</sub>PO<sub>4</sub>/  
 Methanesulfonic Acid = 120/1000/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Ipratropium Bromide Derivative  
 2; Ipratropium Bromide




## ● Indomethacin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.1% H<sub>3</sub>PO<sub>4</sub> = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; 4-Chlorobenzoic Acid (1.0µg)  
 2; Butyl *p*-Hydroxybenzoate (0.6µg)  
 3; Indometacin(Indomethacin) (1.0µg)




## 1) Drugs

## ● Estradiol Benzoate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV230nm

Sample: 1; Progesterone (0.15µg)  
 2; Estradiol Benzoate (0.65µg)




## ● Estriol

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 50/50  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV280nm

Sample: 1; Estriol (0.5µg)  
 2; Methyl Benzoate (0.5µg)




## ● Epirubicin Hydrochloride

**COSMOSIL Application Data**

Column: 5TMS-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 6.9mmol/l Sodium Lauryl Sulfate-  
 H<sub>2</sub>O/Acetonitrile/Methanol/H<sub>3</sub>PO<sub>4</sub>  
 = 540/290/170/1  
 Flow rate: 2.0 ml/min  
 Temperature: 35°C  
 Detection: UV254nm

Sample: 1; 2-Naphthalenesulfonic Acid Sodium(5.0µg)  
 2; Doxorubicin Hydrochloride (5.0µg)  
 3; Epirubicin Hydrochloride (10µg)  
 4; Daunorubicin Hydrochloride (5.0µg)




## ● Oxytocin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 100mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 B; Acetonitrile/H<sub>2</sub>O = 50/50  
 B conc. 30→60%(0→30min)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; (Arg<sup>8</sup>)-Vasopressin (5.0µg)  
 2; Oxytocin Acetate Salt (5.0µg)




## ● Omeprazole

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.9mmol/l Na<sub>2</sub>HPO<sub>4</sub>,  
 1.35mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH7.6 with 1% H<sub>3</sub>PO<sub>4</sub>)  
 = 11/29  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Omeprazole (0.10µg)  
 2; 1,2-Dinitrobenzene (0.25µg)




## ● Captopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.1% Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Captopril (2.6µg)




## 1) Drugs

## ● Gabexate Mesilate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ (0.1% Sodium Lauryl Sulfate/  
 0.5% Sodium *I*-Heptanesulfonate/  
 Acetic Acid = 200/20/1) = 71/29  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV245nm

Sample: 1; Ethyl *p*-Hydroxybenzoate (0.39µg)  
 2; Butyl *p*-Hydroxybenzoate (0.39µg)  
 3; Gabexate Mesilate (0.75µg)




AP-0673

## ● Camostat Mesilate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ (0.2% Sodium *I*-Heptane Sulfonate/  
 0.1% Sodium Lauryl Sulfate/ Acetic Acid  
 = 1000/500/10) = 55/45  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV265nm

Sample: 1; Camostat Mesilate (1.02µg)  
 2; Butyl *p*-Hydroxybenzoate (0.54µg)




AP-0528

## ● Carbazochrome Sodium Sulfonate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Ethanol/ 10.4mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>  
 = 75/925 (pH7 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV360nm

Sample: 1; Carbazochrome Sodium Sulfonate (1.0µg)  
 2; Carbazochrome (1.0µg)




AP-0535

## ● Carbidopa

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Ethanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 5/95  
 (pH2.7 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Methyldopa (10µg)  
 2; Carbidopa (10µg)




AP-0538

## ● Quinidine Sulfate and Quinine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 0.088% Acetic Acid,  
 0.16% Methanesulfonic Acid,  
 0.22% Diethylamine = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm

Sample: 1; Cinchonidine (1.0µg)  
 2; Quinidine Sulfate (10µg)  
 3; Quinine Hydrochloride (10µg)  
 4; Hydroquinidine Hydrochloride (10µg)  
 5; Dihydroquinine




AP-0828

## ● Potassium Clavulanate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 B; Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 (pH4.0 with H<sub>3</sub>PO<sub>4</sub>) = 50/50  
 B conc. 0% (0-4 min), 0→100% (4→15 min),  
 100% (15-25 min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm

Sample: 1; Potassium Clavulanate  
 2; Sodium Amoxicillin



AP-0810


## 1) Drugs

## ● Clarithromycin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 67mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample: 1; Clarithromycin (2.8µg)




AP-0578

## ● Griseofulvin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Griseofulvin (2.2µg)  
 2; Butyl *p*-Hydroxybenzoate (1.6µg)




AP-0676

## ● Clindamycin Phosphate

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 22/78  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV210nm

Sample: 1; Lincomycin Hydrochloride (4.0µg)  
 2; Clindamycin Phosphate (4.0µg)  
 3; Methyl *p*-Hydroxybenzoate (0.3µg)




AP-0584

## ● Clindamycin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH7.5 with KOH) = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV210nm

Sample: 1; Clindamycin Hydrochloride (15µg)




AP-0581

## ● Cloxacillin Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>  
 = 25/75 (pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Guaiacol Glycerol Ether (5.0µg)  
 2; Cloxacillin Sodium Monohydrate (10µg)




AP-0585

## ● Clofibrate

**COSMOSIL Application Data**

Column: 5CN-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Hexane/2-Propanol/Acetic Acid  
 = 1970/30/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm

Sample: 1; Clofibrate (2.0µg)  
 2; 4-Chlorophenol (0.12µg)  
 3; 4-Ethoxyphenol (0.12µg)



AP-0588


## 1) Drugs

## ● Clofibrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.1%Acetic Acid = 60/40  
 Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV275nm

Sample: 1; 4-Chlorophenol (4.0µg)  
 2; Ibuprofen (60µg)  
 3; Clofibrate (10µg)




NACALAI TESQUE, INC  
 AP-0590

## ● Clobetasol Propionate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 47.5/10/42.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV240nm

Sample: 1; Clobetasol Propionate (1.0µg)  
 2; Beclomethasone Dipropionate (2.0µg)




NACALAI TESQUE, INC  
 AP-1098

## ● Cloperastine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 100mmol/l KH<sub>2</sub>PO<sub>4</sub>, 0.16%Perchloric acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV222nm

Sample: 1; Cloperastine Hydrochloride (0.2µg)  
 2; Benzophenone (0.4µg)  
 3; 4-Chlorobenzophenone (0.2µg)




NACALAI TESQUE, INC  
 AP-0593

## ● Chlorpheniramine Maleate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 74.5mmol/l (NH<sub>4</sub>)H<sub>2</sub>PO<sub>4</sub>, 14.7mmol/l H<sub>3</sub>PO<sub>4</sub> = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV225nm

Sample: 1; Maleic Acid (Impurity)  
 2; Chlorpheniramine Maleate (0.06µg)




NACALAI TESQUE, INC  
 AP-0568

## ● Chlorpropamide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 1%Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm

Sample: 1; Chlorpropamide (1.0µg)




NACALAI TESQUE, INC  
 AP-0570

## ● Chlormadinone Acetate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 65/35  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV236nm

Sample: 1; Butyl *p*-Hydroxybenzoate (0.2µg)  
 2; 17 $\alpha$ -Acetoxypregesterone (0.2µg)  
 3; Chlormadinone Acetate (0.8µg)



NACALAI TESQUE, INC  
 AP-0563


## 1) Drugs

## ● Ketoprofen

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O/ 0.5mol KH<sub>2</sub>PO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 43/55/2  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV233nm

Sample: 1; Ketoprofen (0.4µg)



NACALAI TESQUE, INC  
 AP-0709


## ● Cortisone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Acetonitrile/ H<sub>2</sub>O = 30/70  
 B; Acetonitrile/ H<sub>2</sub>O = 70/30  
 B conc. 10% (0-5min)  
 10→90% (5→25min)  
 90% (25-30min)

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Cortisone (3.6µg)  
 2; Cortisone Acetate (3.75µg)  
 3; Hydrocortisone Acetate (3.6µg)




NACALAI TESQUE, INC  
 AP-0595

## ● Salicylic Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 1.67% Acetic Acid = 40/60  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV270nm

Sample: 1; *p*-Hydroxybenzoic Acid (0.025µg)  
 2; Phenol (0.10µg)  
 3; Salicylic Acid (0.50µg)




NACALAI TESQUE, INC  
 AP-1123

## ● Santonin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Santonin (0.50µg)  
 2; Ethyl *p*-Hydroxybenzoate (0.24µg)




NACALAI TESQUE, INC  
 AP-0849

## ● Cyanocobalamins

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 70mmol/l Na<sub>2</sub>HPO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 53/147  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV361nm

Sample: 1; Cyanocobalamin  
 2; Cyanocobalamin Derivative




NACALAI TESQUE, INC  
 AP-1101

## ● Digitoxin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Digoxin (1.0µg)  
 2; Gitoxin (0.5µg)  
 3; Digitoxin (0.5µg)  
 4; Acenaphthene (0.03µg)



NACALAI TESQUE, INC  
 AP-0620


## 1) Drugs

## ● Cyclosporin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.22% $H_3PO_4$  /  
*tert*-Butyl Methyl Ether = 50/45/5  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample: 1; Cyclosporin(Cyclosporin) (24 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0574

## ● Diclofenac Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.12%Acetic Acid = 70/30  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample:  
 1; Ethyl *p*-Hydroxybenzoate (0.14 $\mu$ g)  
 2; Propyl *p*-Hydroxybenzoate (0.2 $\mu$ g)  
 3; Diclofenac Sodium (0.4 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0613

## ● Diclofenamide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l Na<sub>2</sub>HPO<sub>4</sub>,  
 52mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample:  
 1; Diclofenamide (10 $\mu$ g)  
 2; Butyl *p*-Hydroxybenzoate (1.2 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0615

## ● Disulfiram

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV210nm

Sample: 1; Benzophenone (0.05 $\mu$ g)  
 2; Disulfiram (0.05 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0625

## ● Zidovudine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265nm

Sample: 1; Thymine (0.08 $\mu$ g)  
 2; Zidovudine (2.0 $\mu$ g)




NACALAI TESQUE, INC  
 AP-1130

## ● Dipyrindamole

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 7.5mmol/l KH<sub>2</sub>PO<sub>4</sub> = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; Dipyrindamole (2.8 $\mu$ g)  
 2; *p*-Terphenyl (1.2 $\mu$ g)



NACALAI TESQUE, INC  
 AP-0624


## 1) Drugs

## ● Josamycin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.8mol Sodium Perchlorate  
 (pH2.5 with HCl) = 40/60  
 Flow rate: 2.0 ml/min  
 Temperature: 40°C  
 Detection: UV231nm

Sample: 1; Josamycin (10 $\mu$ g)



NACALAI TESQUE, INC  
 AP-0706


## ● Diltiazem Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/Acetonitrile/  
 118mmol/l CH<sub>3</sub>COONa,  
 6.5mmol/l *d*-Camphorsulfonic Acid  
 = 25/25/50(pH6.68)

Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV240nm

Sample: 1; Diltiazem Hydrochloride (3.0 $\mu$ g)  
 2; Phenyl Benzoate (2.0 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0621

## ● Cetirizine Hydrochloride

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l Sulfuric Acid = 94/6  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Cetirizine Dihydrochloride (0.5 $\mu$ g)  
 2; 4-Dimethylaminoantipyrine (0.9 $\mu$ g)




NACALAI TESQUE, INC  
 AP-1096

## ● Cilastatin Sodium Salt

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A; 0.1% $H_3PO_4$   
 B; Acetonitrile/ 0.1% $H_3PO_4$  = 30/70  
 B conc. 15→30%(0→30min), 100%(30-40min)  
 Flow rate: 2.0 ml/min  
 Temperature: 50°C  
 Detection: UV210nm

Sample:  
 1; Cilastatin Sodium Salt (1.0 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0576

## ● Sulbactam Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra *n*-Butylammonium  
 Hydroxide(pH4.0 with H<sub>3</sub>PO<sub>4</sub>) = 25/75  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV220nm

Sample: 1; Sulbactam Sodium (10 $\mu$ g)  
 2; Ethyl *p*-Hydroxybenzoate (7.0 $\mu$ g)




NACALAI TESQUE, INC  
 AP-0853

## ● Cetirizine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 1.7mmol/l Sodium *I*-Heptanesulfonate  
 = 42/58 (pH3.0 with 0.5mol/l Sulfuric Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Cetirizine Dihydrochloride (2.0 $\mu$ g)  
 2; Propyl *p*-Hydroxybenzoate (4.0 $\mu$ g)



NACALAI TESQUE, INC  
 AP-1097


## 1) Drugs

## ● Cefaclor

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.4 with H<sub>3</sub>PO<sub>4</sub>) = 6/94  
 Flow rate: 2.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Cefaclor (2.0µg)  
 2; 4-Aminoacetophenone (2.9µg)



AP-0541


## ● Cefaclor

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A: 50mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH4.0 with H<sub>3</sub>PO<sub>4</sub>)  
 B: Acetonitrile/Buffer A = 45/55  
 B conc. 5→25%(0-30min), 25→100%(30-45min),  
 100%(45-55min)

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Cefaclor (1.0µg)



AP-0542

## ● Cefadroxil

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 10mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/340  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV262nm

Sample: 1; Cefadroxil (1.3µg)



AP-0545


## ● Cephalothin Sodium Salt

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/Ethanol/ 158mmol/l CH<sub>3</sub>COONa,  
 0.076%Acetic Acid(pH5.9 with NaOH)  
 =15/7/78

Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Similar compound (0.25µg)  
 2; Cephalothin Sodium Salt (0.25µg)



AP-0548


## ● Cefsulodin Sodium Salt

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: Acetonitrile/1%(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> = 3/97  
 B: Acetonitrile/1%(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> = 8/92  
 B conc. 0%(0-14min), 100%(14-30min)

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Isonicotinamide (0.2µg)  
 2; Cefsulodin Sodium Salt Hydrate (0.2µg)



AP-0553


## ● Ceftazidime

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>PO<sub>4</sub> (pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 13/87

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Ceftazidime (2.5µg)  
 2; Acetanilide (2.5µg)



AP-0555


## 1) Drugs

## ● Ceftazidime

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>PO<sub>4</sub> = 30/70(pH7.0 with NH<sub>3</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Ceftazidime (0.5µg)  
 2; Pyridine (0.5µg)




AP-0557

## ● Daunorubicin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 38/62 (pH2.2 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Daunorubicin Hydrochloride (5.0µg)  
 2; 2-Naphthalenesulfonic Acid (10µg)



AP-0602


## ● Ceftriaxone Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 18.1mmol/l tetra-*n*-heptylammonium bromide-Acetonitrile / 4.1mmol/l Na<sub>2</sub>HPO<sub>4</sub>, 2.6mmol/l KH<sub>2</sub>PO<sub>4</sub>, 0.88mmol/l Citric Acid, 1.8mmol/l NaOH = 45/55

Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Ceftriaxone Sodium (0.5µg)  
 2; Diethyl Terephthalate (0.9µg)




AP-0560

## ● Dantrolene Sodium

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Hexane/Acetic Acid/Ethanol = 90/10/9  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV300nm

Sample: 1; Theophylline (1.0µg)  
 2; Dantrolene Sodium (0.06µg)




AP-0601

## ● Thiamine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/Acetonitrile/ 5mmol/l Sodium *I*-Octanesulfonate, 1% H<sub>3</sub>PO<sub>4</sub> = 24/16/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Thiamine Hydrochloride (4.0µg)  
 2; Methyl Benzoate (20µg)




AP-0857

## ● Thiopental Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 30/70  
 Flow rate: 1.5 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Isopropyl *p*-Hydroxybenzoate (1.0µg)  
 2; Propyl *p*-Hydroxybenzoate (1.0µg)  
 3; Thiopental Sodium Isomer (4.8µg)  
 4; Thiopental Sodium (4.8µg)



AP-0859




## 1) Drugs

## • Ticarcillin Sodium

**COSMOSIL Application Data**

Column: 5TMS-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tetra-*n*-Butylammonium Bromide, 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH3.0 with H<sub>3</sub>PO<sub>4</sub>) / Acetic Acid = 225/1000/2.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; *o*-Toluic Acid (2.5µg)  
 2; Ticarcillin Sodium (7.6µg)




## • Tipepidine Hibenazate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran/ 1%Ammonium Acetate = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV254nm

Sample: 1; Hibenazic Acid  
 2; Tipepidine  
 3; Propyl *p*-Hydroxybenzoate (0.6µg)




## • Tipepidine Hibenazate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.2%Ammonium Acetate = 65/35  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Hibenazic Acid  
 2; Tipepidine  
 3; Xanthene (0.8µg)




## • Tipepidine Hibenazate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: (Acetonitrile/ 2-Propanol = 2/1) / 0.1%Sodium Lauryl Sulfate, 0.2%H<sub>3</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Hibenazic Acid  
 2; Tipepidine  
 3; Dibucaine Hydrochloride (1.6µg)




## • Timolol Maleate

**COSMOSIL Application Data**

Column: SPE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/Methanol/ 5mmol/l Sodium *l*-Hexanesulfonate, 0.3%Triethylamine, 0.4%Formic Acid(pH3.0 with Formic Acid) = 5/25/70  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; Maleic Acid  
 2; Timolol Maleate (7.5µg)




## • Dexamethasone

**COSMOSIL Application Data**

Column: SPE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 21 mmol/l Ammonium Formate (pH3.6 with Formic Acid) = 33/67  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Dexamethasone (0.18µg)




## 1) Drugs

## • Deferoxamine Mesilate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 2-Propanol/ 10.5mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>, 1.05mmol/l EDTA, 5.6mmol/l Sodium *l*-Heptanesulfonate (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 10/90  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Deferoxamine Mesilate (6.4µg)  
 2; Methyl *p*-Hydroxybenzoate (1.6µg)




## • Doxycycline Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ *N,N*-Dimethyloctylamine/ 100mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 550/3/450 (pH8.0 with NaOH)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Doxycycline Hydrochloride (10µg)




## • Doxifluridine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Doxifluridine (0.5µg)  
 2; Caffeine (1.0µg)




## • Doxorubicin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10.4mmol/l Sodium Lauryl Sulfate, 0.14%H<sub>3</sub>PO<sub>4</sub>=50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Doxorubicin Hydrochloride (5.0µg)




## • Tocopherols

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 98/2  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV292nm

Sample: 1;  $\alpha$ -Tocopherol (20µg)  
 2; Tocopherol Acetate (20µg)




## • Tocopherol Nicotinate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV264nm

Sample: 1; Nicotinic Acid (0.75µg)  
 2; Tocopherol (12.5µg)  
 3; Tocopherol Nicotinate (2.5µg)




## 1) Drugs

## • Todralazine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5.4mmol/l Sodium *I*-Heptanesulfonate = 40/60 (pH3.5 with Acetic Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm

Sample: 1; Potassium Biphthalate (0.4µg)  
 2; Todralazine Hydrochloride (0.4µg)




AP-0879

## • Dopamine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 200mmol/l Na<sub>2</sub>HPO<sub>4</sub> (pH3.0 with Citric Acid)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV280nm

Sample: 1; Uracil (0.5µg)  
 2; Dopamine Hydrochloride (1.0µg)




AP-0631

## • Dobutamine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l Tartrate buffer (pH3.0) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Dobutamine Hydrochloride (11µg)  
 2; Salicylamide (8.5µg)




AP-0628

## • Tolazamide

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Hexane/ Water-saturated Hexane/ Tetrahydrofuran/ Ethanol/ Acetic Acid = 475/ 475/ 20/ 15/ 9  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Tolbutamide (7.5µg)  
 2; Tolazamide (30µg)




AP-0882

## • Tranexamic Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 152mmol/l NaH<sub>2</sub>PO<sub>4</sub>, 8mmol/l Sodium Lauryl Sulfate, 0.83%Triethylamine(pH2.5 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV220nm

Sample: 1; Tranexamic Acid (6.0µg)  
 2; 4-(Aminomethyl) benzoic Acid(0.06µg)




AP-0888

## • Triamcinolone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 25/75  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Triamcinolone (1.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.4µg)



AP-0889


## 1) Drugs

## • Triamcinolone Acetonide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV240nm

Sample: 1; Prednisolone (0.4µg)  
 2; Triamcinolone Acetonide (0.8µg)




AP-0892

## • Trimetazidine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Methanol/ 14.2mmol/l Sodium *I*-Heptanesulfonate (pH 3.0 with 10% H<sub>3</sub>PO<sub>4</sub>) = 2/3  
 B; Methanol  
 B conc. 5%→25% 50min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: 1; *I*-(2,3,4-Trimethoxybenzyl)piperazine Dihydrochloride (40µg)




AP-1125

## • Trimetazidine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm

Sample: 1; Trimetazidine Hydrochloride (0.15µg)  
 2; *p*-Hydroxybenzoic Acid (0.18µg)




AP-0895

## • Tolnaftate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Diphenyl Phthalate (18µg)  
 2; Tolnaftate (4.0µg)




AP-0885

## • Domperidone

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 15.6mmol/l K<sub>2</sub>HPO<sub>4</sub> (pH3.5 with 20mmol/l H<sub>3</sub>PO<sub>4</sub>) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV287nm

Sample: 1; Domperidone (0.1µg)  
 2; Ethyl *p*-hydroxybenzoate (0.2µg)




AP-1102

## • Naphazoline and Chlorpheniramine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.2% Sodium Lauryl Sulfate, 0.1% H<sub>3</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; *o*-Ethoxybenzamide (5.1µg)  
 2; Naphazoline Nitrate (2.2µg)  
 3; Chlorpheniramine Maleate (4.1µg)



AP-0954


## 1) Drugs

## • Nabumetone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ Tetrahydrofuran/ 0.1% Acetic Acid = 28/12/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Nabumetone (10.0µg)




NACALAI TESQUE, INC  
 AP-1113

## • Nalidixic Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 40mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH2.8 with H<sub>3</sub>PO<sub>4</sub>) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV260nm

Sample: 1; Methyl *p*-Hydroxybenzoate (0.13µg)  
 2; Nalidixic Acid (0.005µg)




NACALAI TESQUE, INC  
 AP-0770

## • Nicardipine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.05%Perchloric Acid = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Nicardipine Hydrochloride (0.4µg)  
 2; Nifedipine (0.4µg)




NACALAI TESQUE, INC  
 AP-0775

## • Nicardipine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 10mmol/l KH<sub>2</sub>PO<sub>4</sub> = 68/32  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Nicardipine Hydrochloride (0.4µg)  
 2; Di-*n*-butyl Phthalate (1.3µg)




NACALAI TESQUE, INC  
 AP-0778

## • Nicotinic Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 5mmol/l Sodium *I*-Octanesulfonate-Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub>(pH3.0) = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV260nm

Sample: 1; Nicotinic Acid (1.0µg)  
 2; Caffeine (1.0µg)




NACALAI TESQUE, INC  
 AP-0788

## • Nicotinamide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Heptanesulfonate = 30/70  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Nicotinic Acid (4.0µg)  
 2; Nicotinamide (0.8µg)



NACALAI TESQUE, INC  
 AP-0784


## 1) Drugs

## • Nicergolin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub> (pH7.0 with Triethylamine) = 35/30/35  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV288nm

Sample: 1; Nicergoline (20µg)




NACALAI TESQUE, INC  
 AP-0781

## • Nitrendipine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/Tetrahydrofuran/H<sub>2</sub>O = 20/24/56  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Propyl *p*-Hydroxybenzoate (0.15µg)  
 2; Nitrendipine (0.50µg)




NACALAI TESQUE, INC  
 AP-0790

## • Neostigmine Methylsulfate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l NaH<sub>2</sub>PO<sub>4</sub> (pH3.0 with H<sub>3</sub>PO<sub>4</sub>), 5mmol/l Sodium *I*-Pentanesulfonate = 11/89  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV259nm

Sample: 1; Dimethylaminophenol (0.8µg)  
 2; Neostigmine Methylsulfate (5.0µg)




NACALAI TESQUE, INC  
 AP-0772

## • Baclofen

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.11%Acetic acid = 60/40  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV268nm

Sample: 1; Baclofen (5.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.063µg)




NACALAI TESQUE, INC  
 AP-0477

## • Haloperidol

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 3.5mmol/l Sodium Lauryl Sulfate-Methanol/ 10mmol/l *tri*-Sodium Citrate(pH3.5 with HCl) = 75/25  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV220nm

Sample: 1; *4*-(4-Chlorophenyl)-4-hydroxypiperidine (0.20µg)  
 2; Haloperidol (0.20µg)  
 3; Biphenyl (0.16µg)



NACALAI TESQUE, INC  
 AP-0680


## • Vancomycin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A: Acetonitrile/Tetrahydrofuran/ 0.2%Triethylamine (pH3.2 with H<sub>3</sub>PO<sub>4</sub>) = 7/1/92  
 B: Acetonitrile/Tetrahydrofuran/ 0.2%Triethylamine (pH3.2 with H<sub>3</sub>PO<sub>4</sub>) = 29/1/70  
 B conc. 0%(0-12min), 0→100%(12→20min), 100%(20-22min)

Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV280nm

Sample: 1; Similar compound 1  
 2; Vancomycin Hydrochloride (5.0µg)  
 3; Similar compound 2



NACALAI TESQUE, INC  
 AP-0901


## 1) Drugs

## ● Bisoprolol Fumarate

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 30mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH 2.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV225nm

Sample: 1; Fumaric Acid  
 2; Bisoprolol




AP-1090

## ● Hydrochlorothiazide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 (pH 3.0 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; 4-Amino-6-chlorobenzene-1,3-disulfonamide (1.0 µg)  
 2; Hydrochlorothiazide (3.0 µg)  
 3; 4-Aminoacetophenone (3.6 µg)



AP-0683


## 1) Drugs

## ● Piperacillin Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Acetic Acid,  
 25mmol/l Triethylamine = 21/79  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Ampicillin (10 µg)  
 2; Acetanilide (0.75 µg)  
 3; Piperacillin Sodium (2.5 µg)



AP-0799


## ● Piperacillin Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A: Acetonitrile/H<sub>2</sub>O/ 200mmol/l KH<sub>2</sub>PO<sub>4</sub> = 4/45/1  
 B: Acetonitrile/H<sub>2</sub>O/ 200mmol/l KH<sub>2</sub>PO<sub>4</sub> = 25/24/1  
 B conc. 0%(0-7 min), 0→17%(7→13 min),  
 17%(13-25 min), 17→80%(25→40min),  
 80%(40-44 min)

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Piperacillin Sodium (0.4 µg)




AP-0802

## ● Hydrocortisone

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Chloroform/Methanol/Acetic Acid  
 = 1000/20/1  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Prednisone (0.9 µg)  
 2; Hydrocortisone (2.0 µg)




AP-0685

## ● Hydrocortisone Succinate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 40mmol/l CH<sub>3</sub>COONa  
 (pH 4.0 with Acetic Acid) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Hydrocortisone (0.4 µg)  
 2; Hydrocortisone Succinate (1.0 µg)  
 3; Butyl p-Hydroxybenzoate (0.3 µg)




AP-0689

## ● Pimaricin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 0.1% Ammonium Acetate-  
 Methanol/Tetrahydrofuran/H<sub>2</sub>O = 44/2/47  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV303nm

Sample: 1; Pimaricin (0.02 µg)




AP-0796

## ● Vincristine Sulfate

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 1.67% Diethylamine  
 (pH 7.5 with H<sub>3</sub>PO<sub>4</sub>) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV297nm

Sample: 1; Vincristine Sulfate (10 µg)




AP-1128

## ● Hydrocortisone Acetate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Hydrocortisone (2.0 µg)  
 2; Hydrocortisone Acetate (4.0 µg)  
 3; Benzyl p-Hydroxybenzoate (2.0 µg)




AP-0686

## ● Hypromellose Phthalate (impurity)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Cyanoacetic Acid  
 = 15/85  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm

Sample: 1; Phthalic Acid (0.5 µg)




AP-0933

## ● Famotidine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/Methanol/  
 10mmol/l Sodium l-Heptanesulfonate  
 (pH 3.0 with Acetic Acid) = 19/3/78  
 Flow rate: 0.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Famotidine (0.5 µg)  
 2; Methyl p-Hydroxybenzoate (0.8 µg)




AP-0651

## ● Scopolamine Butyl Bromide

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 19mmol/l Sodium Lauryl Sulfate  
 = 68/37 (pH 3.6 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Scopolamine Hydrobromide (2.0 µg)  
 2; Scopolamine n-Butyl Bromide (2.0 µg)



AP-0852


## 1) Drugs

## • Bufexamac

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/Acetonitrile/  
 13.6mmol/l Sodium *I*-Octane Sulfonate  
 0.94%Acetic Acid, 1.9mmol/l EDTA  
 = 24/24/52  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm

Sample: 1; Bufexamac (2µg)  
 2; Diphenylimidazole (0.16µg)




AP-0522

## • Pravastatin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O/Acetic Acid/Triethylamine  
 = 500/500/1/1  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV238nm

Sample: 1; Ethyl *p*-Hydroxybenzoate (0.7µg)  
 2; Pravastatin Sodium (1.6µg)



AP-0811


## 1) Drugs

## • Fluocinolone Acetonide

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Water-saturated Chloroform/Methanol/  
 Acetic Acid = 200/3/2  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Triamcinolone Acetonide (3.0µg)  
 2; Fluocinolone Acetonide (3.0µg)




AP-0662

## • Fluorometholone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV254nm

Sample: 1; Fluorometholone (0.20µg)  
 2; Butyl *p*-Hydroxybenzoate (0.16µg)




AP-0666

## • Flavin Adenine Dinucleotide Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.2%KH<sub>2</sub>PO<sub>4</sub> = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV260nm

Sample: 1; Adenosine (1.0µg)  
 2; FAD (4.0µg)  
 3; FMN (4.0µg)  
 4; Riboflavin (1.0µg)




AP-0653

## • Fluoxymesterone

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: *n*-Butyl Chloride/Water-saturated *n*-Butyl  
 Chloride/Tetrahydrofuran/Methanol/  
 Acetic Acid = 95/95/14/7/6  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Fluoxymesterone (2.7µg)  
 2; Methyl Prednisolone (2.0µg)




AP-0669

## • Flurbiprofen

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.7%H<sub>3</sub>PO<sub>4</sub> = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Flurbiprofen (0.2µg)




AP-0671

## • Prednisolone Acetate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Prednisolone (1.5µg)  
 2; Prednisolone Acetate (1.0µg)  
 3; Cortisone Acetate (1.5µg)  
 4; Butyl *p*-Hydroxybenzoate (0.5µg)




AP-0814

## • Fluocinonide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Fluocinonide (4.0µg)  
 2; Propyl Benzoate (12.6µg)




AP-0663

## • Fluocinolone Acetonide

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Ethyl *p*-Hydroxybenzoate (0.6µg)  
 2; Isopropyl *p*-Hydroxybenzoate (1.0µg)  
 3; Fluocinolone Acetonide (4.0µg)  
 4; Propyl *p*-Hydroxybenzoate (1.0µg)




AP-0660

## • Procaine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Sodium *I*-Pentanesulfonate,  
 50mmol/l KH<sub>2</sub>PO<sub>4</sub>(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
 = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Procaine Hydrochloride (1.25µg)  
 2; Caffeine (1.25µg)




AP-0817

## • Procatamol Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Pentanesulfonate/  
 Acetic Acid = 23/76/1  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Procatamol Hydrochloride (0.06µg)



AP-0821




## 1) Drugs

## ● Propranolol Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 12mmol/l Sodium Lauryl Sulfate,  
 0.7mmol/l Tetra-*n*-butylammonium Phosphate/  
 Sulfuric Acid = 550/450/1(pH3.3 with NaOH)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV292nm

Sample:  
 1; Propranolol hydrochloride (0.08µg)



NACALAI TESQUE, INC


AP-0823

## ● Flopropione

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 1.2% H<sub>3</sub>PO<sub>4</sub> = 57/43  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV267nm

Sample: 1; Flopropione (0.40µg)  
 2; Ethyl 4-Hydroxybenzoate (0.26µg)



NACALAI TESQUE, INC


AP-0657

## ● Bromhexine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7.3mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH7.0 with NaOH) = 80/20  
 Flow rate: 2.0 ml/min  
 Temperature: 40°C  
 Detection: UV245nm

Sample: 1; Bamethane Sulfate (25µg)  
 2; Bromhexine Hydrochloride (1.25µg)



NACALAI TESQUE, INC


AP-0518

## ● Beclometasone Dipropionate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Beclometasone Dipropionate (1.6µg)  
 2; Testosterone Propionate (0.8µg)



NACALAI TESQUE, INC


AP-0479

## ● Bezafibrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Acetic Acid = 9/4  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; *p*-Chlorobenzoic Acid (0.55µg)  
 2; Bezafibrate (1.05µg)



NACALAI TESQUE, INC


AP-0510

## ● Bezafibrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 0.1% Acetic Acid = 9/4  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; *p*-Nitrophenol (0.32µg)  
 2; Bezafibrate (0.44µg)



NACALAI TESQUE, INC

AP-0513


## 1) Drugs

## ● Betahistine Mesilate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 8mmol/l Sodium Lauryl Sulfate-  
 Acetonitrile/ 2% Acetic Acid, 0.5% Diethylamine  
 = 37/63  
 Flow rate: 1.0 ml/min  
 Temperature: 35°C  
 Detection: UV261nm

Sample:  
 1; 2-Vinylpyridine (0.16µg)  
 2; Betahistine Mesilate (0.16µg)



NACALAI TESQUE, INC


AP-0488

## ● Betamethasone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV241nm

Sample: 1; Betamethasone (0.06µg)



NACALAI TESQUE, INC


AP-0497

## ● Betamethasone Sodium Phosphate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Tetra-*n*-butylammonium Bromide  
 8.9mmol/l Na<sub>2</sub>HPO<sub>4</sub>·50.7mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Betamethasone Sodium Phosphate (1.0µg)  
 2; Butyl *p*-Hydroxybenzoate (0.2µg)



NACALAI TESQUE, INC


AP-0494

## ● Betamethasone Valerate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 65/35  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Betamethasone Valerate (1.5µg)  
 2; Beclometasone Dipropionate (1.5µg)



NACALAI TESQUE, INC


AP-0503

## ● Berberine Chloride Hydrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 6mmol/l Sodium Lauryl Sulfate  
 25mmol/l KH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV345nm

Sample: 1; Palmatine Chloride (1.0µg)  
 2; Berberine Chloride (1.0µg)



NACALAI TESQUE, INC


AP-0487

## ● Benzylpenicillin Potassium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>  
 = 24/76 (pH8 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample:  
 1; Benzylpenicillin Potassium (4.0µg)  
 2; Methyl *p*-Hydroxybenzoate (0.06µg)



NACALAI TESQUE, INC


AP-0482

## 1) Drugs

## ● Calcium Folate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: [H<sub>2</sub>O/Acetonitrile/Methanol/  
 40%Tetra-*n*-butylammonium Hydroxide  
 = 760/200/8.6/9.4 (pH7.5 with NaH<sub>2</sub>PO<sub>4</sub>)  
 →1000(with H<sub>2</sub>O)]  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample: 1; Calcium Folate (2.8µg)  
 2; Folic Acid (0.8µg)




AP-0524

## ● Calcium Folate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 8mmol/l Na<sub>2</sub>HPO<sub>4</sub>/ Tetrabutylammonium  
 Hydroxide = 110/385/4 (pH 7.5 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 45°C  
 Detection: UV254nm  
 Sample: 1; Calcium Folate (2.0µg)  
 2; Folic Acid (2.0µg)




AP-1093

## ● Mitomycin C

**COSMOSIL Application Data**

Column: 5PE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l Ammonium Acetate,  
 0.025%Acetic Acid = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV365nm  
 Sample: 1; Mitomycin C (5.0µg)  
 2; Ethyl Vanillin (75µg)




AP-0768

## ● Mizoribine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.067% H<sub>3</sub>PO<sub>4</sub> aq.  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV279nm  
 Sample: 1; Mizoribine (1.0µg)




AP-1108

## ● Minocycline Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 2.8%Ammonium Oxalate/*N,N*-dimethylformamide/  
 100mmol/l EDTA = 11/5/4  
 (pH6.2 with Tetrabutylammonium Hydroxide)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; 4-*epi*-Minocycline  
 2; Minocycline Hydrochloride




AP-0767

## ● Mexiletine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 14.4mmol/l Sodium Lauryl Sulfate,  
 25mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 21/30  
 Flow rate: 2.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm  
 Sample: 1; Phenetylamine Hydrochloride (0.6µg)  
 2; Mexiletine Hydrochloride (1.0µg)




AP-0766

## 1) Drugs

## ● Mecobalamin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 100mmol/l Sodium *L*-Hexanesulfonate,  
 Acetonitrile/ 20mmol/l Phosphate  
 buffer(pH3.5) = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV266nm  
 Sample: 1; Cyanocobalamin (0.5µg)  
 2; Hydroxocobalamin Acetate (0.5µg)  
 3; Mecobalamin (1.0µg)




AP-0738

## ● Meticrane

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm  
 Sample: 1; Caffeine (0.2µg)  
 2; Meticrane (0.2µg)




AP-0754

## ● Meticrane

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV230nm  
 Sample: 1; Meticrane (0.4µg)  
 2; Methyl *p*-Hydroxybenzoate (0.4µg)




AP-0757

## ● Methotrexate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 200mmol/l Na<sub>2</sub>HPO<sub>4</sub>  
 (pH6.0 with 100mmol/l Citric Acid) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV302nm  
 Sample: 1; Folic Acid (1.0µg)  
 2; Methotrexate (1.0µg)




AP-0751

## ● Metformin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 4.5mmol/l Sodium Lauryl Sulfate,  
 0.04%H<sub>2</sub>PO<sub>4</sub> = 38/62  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV235nm  
 Sample: 1; Metformin Hydrochloride (0.50µg)  
 2; Isobutyl *p*-Hydroxybenzoate (0.95µg)




AP-0761

## ● Metronidazole

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV320nm  
 Sample: 1; Metronidazole (0.24µg)



AP-0763


## 1) Drugs

## ● Menatetrenone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV270nm

Sample: 1; Menatetrenone (0.6µg)  
 2; Phytanadione (0.6µg)




NACALAI TESQUE, INC  
 AP-0742

## ● Mepivacaine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 10mmol/l Sodium Lauryl Sulfate-Acetonitrile/  
 20mmol/l Phosphate buffer(pH3.0) = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Mepivacaine Hydrochloride (2.0µg)  
 2; Benzophenone (0.04µg)



NACALAI TESQUE, INC  
 AP-0746


## 1) Drugs

## ● Latamoxef Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 100mmol/l Ammonium  
 Acetate = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Latamoxef Sodium (isomer 1)  
 2; Latamoxef Sodium (isomer 2)




NACALAI TESQUE, INC  
 AP-0716

## ● Liothyronine Sodium

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Propyl *p*-Hydroxybenzoate (0.14µg)  
 2; Liothyronine Sodium (0.08µg)




NACALAI TESQUE, INC  
 AP-0728

## ● Meropenem

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 0.1% Triethylamine  
 (pH5.0 with H<sub>3</sub>PO<sub>4</sub>) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Meropenem Trihydrate (2.5µg)  
 2; Benzyl Alcohol (1.5µg)



NACALAI TESQUE, INC  
 AP-0749

## ● Ubidecarenone

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ Ethanol = 50/50  
 Flow rate: 0.5 ml/min  
 Temperature: 35°C  
 Detection: UV275nm

Sample: 1; Ubidecarenone (2.5µg)




NACALAI TESQUE, INC  
 AP-0900

## ● Lisinopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 B; Acetonitrile/ 25mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 40/60  
 B conc. 10→50%(0→10min), 50%(10-25min)  
 Flow rate: 1.5 ml/min  
 Temperature: 60°C  
 Detection: UV215nm

Sample: 1; Lisinopril (0.45µg)  
 2; Caffeine (0.015µg)




NACALAI TESQUE, INC  
 AP-0732

## ● Lisinopril

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 25mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 60°C  
 Detection: UV215nm

Sample: 1; Lisinopril (2.0µg)  
 2; Caffeine (0.5µg)




NACALAI TESQUE, INC  
 AP-0733

## ● Iodine, Salicylic Acid and Phenol

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Methanol/ 100mmol/l Phosphate  
 buffer(pH7.0) = 25/75  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm

Sample: 1; Benzoic Acid (2.0µg)  
 2; Theophylline (0.5µg)  
 3; Salicylic Acid (2.0µg)  
 4; Phenol (2.0µg)




NACALAI TESQUE, INC  
 AP-0935

## ● Latamoxef Sodiums

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 51mmol/l KH<sub>2</sub>PO<sub>4</sub>,  
 9mmol/l Na<sub>2</sub>HPO<sub>4</sub>,  
 5mmol/l Tetra-*n*-Butylammonium Bromide  
 = 25/75  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; 5-Mercapto-*l*-methyltetrazole (1.65µg)  
 2; Latamoxef Sodium (isomer 1)  
 3; Latamoxef Sodium (isomer 2)  
 4; *m*-Cresol (7.5µg)




NACALAI TESQUE, INC  
 AP-0711

## ● Lidocaine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 10mmol/l Sodium Lauryl Sulfate-Acetonitrile/  
 20mmol/l Phosphate buffer(pH3.0) = 45/55  
 Flow rate: 1.5 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Lidocaine (8.5µg)  
 2; Benzophenone (0.25µg)




NACALAI TESQUE, INC  
 AP-0724

## ● Ritodrine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 71mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>,  
 7.8mmol/l Sodium *l*-Heptanesulfonate  
 = 30/70(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV220nm

Sample: 1; Ritodrine Hydrochloride  
 2; *threo*-Ritodrine Hydrochloride  
 3; by-product



NACALAI TESQUE, INC  
 AP-0835


## 1) Drugs

## ● Ritodrine Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 71mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>,  
 7.8mmol/l Sodium *J*-Heptanesulfonate  
 = 30/70(pH3.0 with H<sub>3</sub>PO<sub>4</sub>)  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV274nm

Sample:  
 1; Ritodrine Hydrochloride (0.30µg)  
 2; Methyl *p*-Hydroxybenzoate (0.05µg)




AP-0838

## ● Rifampicin

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 20mmol/l Citric Acid, 11mmol/l Sodium Perchlorate-  
 Acetonitrile/H<sub>2</sub>O/ 1mol/l KH<sub>2</sub>PO<sub>4</sub>,  
 55mmol/l H<sub>3</sub>PO<sub>4</sub>(pH3.1) = 7/11/2  
 Flow rate: 2.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample:  
 1; Butyl *p*-Hydroxybenzoate (0.15µg)  
 2; Rifampicin (0.80µg)




AP-0832

## ● Rifampicin

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 15mmol/l Sodium Perchlorate,  
 28mmol/l Citric Acid, 17mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV254nm

Sample: 1; Rifampicin (0.08µg)




AP-0833

## ● Lincomycin Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/Methanol/ 1.35%*H*<sub>3</sub>PO<sub>4</sub>  
 (pH6.0 with Ammonia) = 15/15/78  
 Flow rate: 1.0 ml/min  
 Temperature: 46°C  
 Detection: UV210nm

Sample: 1; Lincomycin Hydrochloride (20µg)




AP-0729

## ● Reserpine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l KH<sub>2</sub>PO<sub>4</sub>  
 (pH3.0 with H<sub>3</sub>PO<sub>4</sub>) = 45/55  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV268nm

Sample: 1; Reserpine (0.20µg)  
 2; Butyl *p*-Hydroxybenzoate (0.08µg)




AP-0830

## ● Levallorphan Tartrate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 7mmol/l Sodium Lauryl Sulfate,  
 0.1%*H*<sub>3</sub>PO<sub>4</sub>(pH3.0 with NaOH) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV280nm

Sample: 1; *iso*-Butyl *p*-Hydroxybenzoate (0.3µg)  
 2; Levallorphan Tartrate (1.7µg)



AP-0718


## 1) Drugs

## ● Levothyroxine

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 0.15%*H*<sub>3</sub>PO<sub>4</sub> = 67/33  
 Flow rate: 1.0 ml/min  
 Temperature: 25°C  
 Detection: UV230nm

Sample: 1; Levothyroxine Sodium (0.1µg)  
 2; Ethinylestradiol (0.2µg)




AP-0721

## ● Roxatidine Acetate Hydrochloride

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O/ Triethylamine/  
 Acetic Acid = 60/340/2/1  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV274nm

Sample: 1; Roxatidine Acetate Hydrochloride (4.2µg)




AP-0840

## ● Roxithromycin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.2mmol/l (NH<sub>4</sub>)<sub>2</sub>HPO<sub>4</sub>  
 (pH5.3 with NaOH) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV205nm

Sample: 1; Roxithromycin (1.0µg)




AP-0845

## ● Tricyclic Drugs

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate  
 buffer(pH7) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Carbamazepine  
 2; Desipramine Hydrochloride  
 3; Imipramine Hydrochloride




AP-0066

## ● Bronchodilators

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 1.0AUFS

Sample: 1; Uracil (0.5µg)  
 2; Theobromine (1.5µg)  
 3; Theophylline (2.0µg)  
 4; Caffeine (2.0µg)  
 5; Phenol (0.8µg)




AP-0076

## ● Antiarrhythmic Drugs

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate  
 buffer(pH7) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample: 1; Phenytoin  
 2; Ketamine Hydrochloride  
 3; Quinidine  
 4; Lidocaine



AP-0067


## 1) Drugs

## ● Antiarrhythmic Drugs

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2) = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Procainamide (0.5µg)  
 2; N-Acetylprocainamide (0.5µg)




NACALAI TESQUE, INC  
 AP-0069

## ● Antiepileptics

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Barbital (7.59µg)  
 2; Phenytoin (10.38µg)  
 3; Carbamazepine (1.02µg)



NACALAI TESQUE, INC  
 AP-0077


## 1) Drugs

## ● Antihyperlipidemic Drugs

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Hydralazine Hydrochloride (1.0µg)  
 2; Todalazine Hydrochloride (1.0µg)




NACALAI TESQUE, INC  
 AP-0080

## ● Profens

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l Acetic Acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 1.0AUFS

Sample: 1; Ketoprofen (1.60µg)  
 2; Ibuprofen (1.69µg)  
 3; Flurbiprofen (1.57µg)




NACALAI TESQUE, INC  
 AP-0081

## ● Analgesics

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphoric Acid = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; Aminoantipyrine (0.74µg)  
 2; Acetaminophen (2.97µg)  
 3; Antipyrine (0.85µg)  
 4; Acetylsalicylic Acid (Aspirine) (0.54µg)  
 5; Phenacetin (1.0µg)  
 6; Salicylic acid (1.0µg)




NACALAI TESQUE, INC  
 AP-0073

## ● Analgesics

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Ethyl Acetate/Hexane = 1/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Acetanilide (1.0µg)  
 2; Phenacetin (1.0µg)  
 3; Acetaminophen (1.0µg)




NACALAI TESQUE, INC  
 AP-0074

## ● Cardiac Glycosides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.32AUFS

Sample: 1; Digitoxigenin (2.5µg)  
 2; Digitoxin (5.0µg)




NACALAI TESQUE, INC  
 AP-0083

## ● Anticancer Drugs

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 80/20  
 Flow rate: 0.3 ml/min  
 Temperature: Room temperature  
 Detection: UV226nm

Sample: 1; cis-Platin (CDDP) (1.46µg)  
 2; Guanosine (0.50µg)




NACALAI TESQUE, INC  
 AP-0380

## ● Histamine H1-Receptor Blockers

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 10mmol/l SDS, 0.1% Phosphoric Acid = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Uracil  
 2; Diphenhydramine Hydrochloride




NACALAI TESQUE, INC  
 AP-0078

## ● Histamine H2-Receptor Blockers

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; Famotidine (1.3µg)  
 2; Cimetidine (33.5µg)  
 3; Ranitidine (1.6µg)



NACALAI TESQUE, INC  
 AP-0079

## ● Tetracyclines Antibiotics

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 20°C  
 Detection: UV254nm

Sample: 1; Oxytetracycline  
 2; Tetracycline  
 3; Chlortetracycline




NACALAI TESQUE, INC  
 AP-0085

## ● Penicillin Antibiotics

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH7) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm, 0.2AUFS

Sample: 1; Carbenicillin (1.5µg)  
 2; Ampicillin (3.0µg)  
 3; Methicillin Sodium Salt (1.5µg)  
 4; Penicillin G Potassium Salt (3.0µg)



NACALAI TESQUE, INC  
 AP-0086




## 1) Drugs

## • Quinolone Antimicrobials

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 55/45  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm, 0.16AUFS

Sample: 1; Ofloxacin (0.39µg)  
 2; Oxolinic Acid (0.08µg)  
 3; Flumequine (1.08µg)  
 4; Piromidic Acid (0.08µg)




NACALAI TESQUE, INC  
 AP-0087

## • Nitrofurant Antimicrobials

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm, 0.5AUFS

Sample: 1; Nitrofurantoin (1.0µg)  
 2; Nitrofurazone (1.0µg)




NACALAI TESQUE, INC  
 AP-0089

## • Streptomycin Sulfate

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Sodium l-Hexanesulfonate, 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV205nm, 0.2AUFS

Sample: Streptomycin Sulfate (5.0µg)




NACALAI TESQUE, INC  
 AP-0091

## • Sulfa Drugs

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS

Sample: 1; Sulfaisomidin (0.24µg)  
 2; Sulfathiazole (0.24µg)  
 3; Sulfamethazine (0.24µg)  
 4; Sulfamethoxy-pyridazine (0.24µg)  
 5; Sulfamethoxazole (0.24µg)  
 6; Sulfachloropyridazine (0.24µg)  
 7; Sulfamonomethoxine (0.24µg)  
 8; Sulfadimethoxine (0.24µg)




NACALAI TESQUE, INC  
 AP-0090

## • Allantoin

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Allantoin  
 2; Allantoic Acid



NACALAI TESQUE, INC  
 AP-1060


## 2) Crude Drugs

## • Bearberry Leaf

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 1mmol/l HCl = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; Arbutin (5.0µg)  
 2; Hydroquinone (5.0µg)  
 3; Gallic Acid (5.0µg)




NACALAI TESQUE, INC  
 AP-0906

## • Scutellaria Root

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.68% Phosphoric Acid = 28/72  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV277nm

Sample: 1; Baicalin (0.1µg)  
 2; Methyl p-Hydroxybenzoate (0.2µg)




NACALAI TESQUE, INC  
 AP-1134

## • Phellodendron Bark

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 5.9mmol/l Sodium Lauryl Sulfate, 25mmol/l KH<sub>2</sub>PO<sub>4</sub> = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV345nm

Sample: 1; Palmatine Chloride Hydrate (1.0µg)  
 2; Berberine Chloride (1.0µg)




NACALAI TESQUE, INC  
 AP-0964

## • Pueraria Root

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV250nm

Sample: 1; Puerarin (1.3µg)




NACALAI TESQUE, INC  
 AP-0969

## • Kamishoyosan Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.11% H<sub>3</sub>PO<sub>4</sub> = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: 1; Geniposide (1.2µg)




NACALAI TESQUE, INC  
 AP-0942

## • Glycyrrhiza

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 2% H<sub>3</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Glycyrrhizic Acid (5.0µg)  
 2; Propyl p-Hydroxybenzoate (1.0µg)



NACALAI TESQUE, INC  
 AP-0927


## 2) Crude Drugs

## • Dried Yeast

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 9.2mmol/l Sodium *l*-Octanesulfonate, 20mmol/l KH<sub>2</sub>PO<sub>4</sub>(pH3.5 with H<sub>3</sub>PO<sub>4</sub>) = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Thiamine Hydrochloride (0.10µg)  
 2; Phenacetin (0.06µg)




NACALAI TESQUE, INC  
 AP-0919

## • Keishibukuryougan Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 1/5  
 Flow rate: 0.8 ml/min  
 Temperature: 45°C  
 Detection: UV210nm

Sample: 1; *D*-(-)-Amygdalin (2.0µg)



NACALAI TESQUE, INC  
 AP-0944


## 2) Crude Drugs

## • Saireito Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l NaH<sub>2</sub>PO<sub>4</sub> = 3/5  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; Saikosaponin b2 (1.0µg)



NACALAI TESQUE, INC  
 AP-0978

## • Saireito Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.5% H<sub>3</sub>PO<sub>4</sub> = 24/76  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV277nm

Sample: 1; Baicalin (0.5µg)




NACALAI TESQUE, INC  
 AP-0981

## • Red Ginseng and Ginseng

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV203nm

Sample: 1; Ginsenoside Rg1 (2.5µg)  
 2; Ginsenoside Re (2.5µg)




NACALAI TESQUE, INC  
 AP-0922

## • Red Ginseng and Ginseng

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV203nm

Sample: 1; Ginsenoside Rb1 (2.5µg)  
 2; Ginsenoside Rc (2.5µg)




NACALAI TESQUE, INC  
 AP-0923

## • Powdered Gardenia Fruit

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV240nm

Sample: 1; Caffeine (0.175µg)  
 2; Geniposide (0.175µg)




NACALAI TESQUE, INC  
 AP-0966

## • Peony Root

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.12% H<sub>3</sub>PO<sub>4</sub> = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV232nm

Sample: 1; Albiflorin (1.0µg)  
 2; Paeoniflorin (1.0µg)




NACALAI TESQUE, INC  
 AP-0960

## • Magnolia Bark

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 2% Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV289nm

Sample: 1; Honokiol (1.0µg)  
 2; Magnolol (1.0µg)




NACALAI TESQUE, INC  
 AP-0948

## • Bupleurum Root

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV206nm

Sample: 1; Saikosaponin a (1.0µg)  
 2; Saikosaponin d (1.0µg)




NACALAI TESQUE, INC  
 AP-0910

## • Toad Venom

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.1% H<sub>3</sub>PO<sub>4</sub> = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV300nm

Sample: 1; Bufalin (0.4µg)  
 2; Cinobufagin (0.4µg)  
 3; Resibufogenin (0.6µg)  
 4; Indometacin (Indomethacin) (0.6µg)




NACALAI TESQUE, INC  
 AP-1000

## • Senna Leaf

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 5mmol/l Tetra-*n*-heptylammonium Bromide- Acetonitrile/ 0.1mol/l Acetic Acid, 1% Sodium Acetate(pH5.0) = 32/68  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV340nm

Sample: 1; Sennoside B (1.0µg)  
 2; Sennoside A (0.5µg)



NACALAI TESQUE, INC  
 AP-0993


## 2) Crude Drugs

## • Swertia Herb

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 10/90  
 Flow rate: 0.5 ml/min  
 Temperature: 50°C  
 Detection: UV238nm

Sample: 1; Theophylline (1.0µg)  
 2; Swertiamarin (1.0µg)




AP-0996

## • Rhubarb

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 1.25%Acetic Acid = 20/80  
 Flow rate: 0.5 ml/min  
 Temperature: 40°C  
 Detection: UV340nm

Sample: 1; Senoside A (2.0µg)  
 2; Naringin (2.0µg)



AP-0972


## 2) Crude Drugs

## • Moutan Bark

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 3.1%Acetic Acid = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV274nm

Sample: 1; Paeonol (0.5µg)  
 2; Butyl *p*-Hydroxybenzoate (2.2µg)




AP-0950

## • Hochuekkito Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 1.2%*H*<sub>3</sub>*P*O<sub>4</sub> = 18/82  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV285nm

Sample: 1; Naringin (0.1µg)  
 2; Hesperidin (0.1µg)




AP-0930

## • Daiokanzoto Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.04%*H*<sub>3</sub>*P*O<sub>4</sub> = 540/2460  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV340nm

Sample: 1; Senoside A (0.26µg)




AP-0916

## • Capsicum

**COSMOSIL Application Data**

Column: 5PE-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 0.1%*H*<sub>3</sub>*P*O<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV281nm

Sample: 1; Vanillynonamide (0.4µg)  
 2; Capsaicin (0.4µg)




AP-0913

## • Povidone

**COSMOSIL Application Data**

Column: 5C<sub>8</sub>-MS  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV254nm

Sample: 1; *l*-Vinyl-2-pyrrolidone (0.1µg)  
 2; Vinyl Acetate (5.0µg)




AP-0965

## • Nux Vomica

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l *KH*<sub>2</sub>*PO*<sub>4</sub>,  
 1.11%Triethylamine = 10/90(pH3.0 with *H*<sub>3</sub>*P*O<sub>4</sub>)  
 Flow rate: 0.5 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Strychnine (0.8µg)  
 2; Barbitol Sodium (1.1µg)




AP-0958

## • Ipecac

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Sodium *l*-Heptanesulfonate  
 (pH4.0 with Acetic Acid) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: UV283nm

Sample: 1; Cephaeline Hydrobromide (2.0µg)  
 2; Emetine Hydrochloride (0.5µg)




AP-0938

## • Belladonna Root

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l *KH*<sub>2</sub>*P*O<sub>4</sub>,  
 1%Triethylamine(pH3.5 with *H*<sub>3</sub>*P*O<sub>4</sub>)  
 = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Atropine Sulfate (2.0µg)  
 2; Brucine Dihydrate (0.5µg)



AP-0909

## • Ryokeijutukanto Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 0.13%*H*<sub>3</sub>*P*O<sub>4</sub> = 25/75  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV273nm

Sample: 1; *trans*-Cinnamic Acid (0.1µg)




AP-0976

## • Scopolia Rhizome and Scopolia Extract

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l *KH*<sub>2</sub>*P*O<sub>4</sub>,  
 0.1%Triethylamine(pH3.5 with *H*<sub>3</sub>*P*O<sub>4</sub>)  
 = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Scopolamine Hydrobromide (0.6µg)  
 2; Atropine Sulfate (2.0µg)  
 3; Brucine Dihydrate (0.4µg)




AP-0989

## 3) Natural Compounds

## • Carotenes

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Tetrahydrofuran/Methanol = 20/80  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV470nm  
 Sample: 1;  $\alpha$ -Carotene  
 2;  $\beta$ -Carotene



## • Flavanones

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH2.5) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; Naringin (0.4µg)  
 2; Naringenin (0.2µg)  
 3; Apigenin (0.2µg)



## • Saikosaponins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 45/55  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD, Gain=6  
 Sample: 1; Saikosaponin c (1.5µg)  
 2; Saikosaponin a (1.5µg)  
 3; Saikosaponin b<sub>2</sub> (1.5µg)  
 4; Saikosaponin b<sub>1</sub> (1.5µg)  
 5; Saikosaponin d (1.5µg)




## 3) Natural Compounds

## • Hydroxyflavones

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 20mmol/l Phosphate buffer(pH2.5) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm  
 Sample: 1; 7-Hydroxyflavone (0.2µg)  
 2; 6-Hydroxyflavone (0.1µg)  
 3; 3-Hydroxyflavone (0.5µg)  
 4; 5-Hydroxyflavone (0.1µg)



## • Coumarins

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS  
 Sample: 1; Coumarin (0.84µg)  
 2; Dihydrocoumarin (3.86µg)



## • Anthraquinone dyes

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS  
 Sample: 1; Alizalin (0.05µg)  
 2; Chrysazin (0.1µg)  
 3; Anthrarufin (0.3µg)  
 6; Amylbenzene (10µg)



## • Chlorophyll

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A: Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 60/40  
 B: Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 95/5  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV405nm  
 Sample: Sodium Copper Chlorophyll (2.0µg)



## • Hinokitiol

**COSMOSIL Application Data**


Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 1mmol/l EDTA, 20mmol/l Phosphoric Acid = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS  
 Sample: Hinokitiol (1.0µg)



## • Anthraquinone dyes

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS  
 Sample: 4; Purpurin (0.2µg)  
 5; Quinizarin (0.1µg)  
 6; Amylbenzene (10µg)




## 3) Natural Compounds

## ● Capsaicin

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: Capsaicin (2.0µg)




AP-1034

## ● Alkaloids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Papaverine Hydrochloride (0.17µg)  
 2; Aconitine (8.1µg)  
 3; Reserpine (1.7µg)




AP-0133

## ● Ginkgo Biloba

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Ginkgolide C (1.3µg)  
 2; Ginkgolide A (1.3µg)  
 3; Ginkgolide B (1.3µg)




AP-1063

## ● Zingiberis Rhizoma

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV225nm

Sample: 1; 6-Gingerol (0.31µg)  
 2; 6-Shogaol (0.13µg)




AP-1016

## ● Watermelon

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; L-Citrulline (7.5µg)  
 2; Malic Acid (3.8µg)



AP-1062


## 4) Pesticides

## ● Asulam

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Phosphate buffer(pH3) = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm

Sample: Asulam




AP-0144

## ● Chlorophenoxyacetic Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 0.1% Acetic Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV275nm, 0.02AUFS

Sample: 1; MCP (0.22µg)  
 2; MCPP (0.20µg)  
 3; MCPB (0.24µg)




AP-0147

## ● Pesticides used at Golf Course

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 50mmol/l Phosphate buffer(pH3) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Thiram  
 2; Iprodione  
 3; Bensulide




AP-0145

## ● Diphenyl Ether Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: 1; Fluorodifen (0.6µg)  
 2; Chlormethoxynil (1.0µg)  
 3; Nitrofen (1.0µg)  
 4; Oxyfluorfen (1.2µg)  
 5; CNP (1.0µg)




AP-0146

## ● Aniline Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: 1; DCMU (0.11µg)  
 2; DCPA (0.08µg)  
 3; Linuron (0.11µg)  
 4; MCC (0.18µg)




AP-0148

## ● Dithiocarbamate Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Maneb (60µg)  
 2; Thiram (10µg)



AP-0149




## 4) Pesticides

## • Carbamate Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample:  
 1; Methomyl (0.5µg)  
 2; Pirimicarb (0.5µg)  
 3; MTMC (2.0µg)  
 4; PHC (2.0µg)  
 5; MPMC (2.0µg)  
 6; NAC (1.0µg)  
 7; XMC (1.0µg)  
 8; Ethiofencarb (2.0µg)  
 9; Isoprocarb (2.0µg)  
 10; BPMC (2.0µg)




NACALAI TESQUE, INC  
 AP-0150

## • Carbamate Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.05AUFS

Sample: Carbosulfan (6.0µg)




NACALAI TESQUE, INC  
 AP-0151

## • Triazine and Urea Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample:  
 1; PAC (0.5µg)  
 2; Bromacil (0.5µg)  
 3; CAT (0.5µg)  
 4; Methabenzthiazuron (0.5µg)  
 5; Chlorotoluron (0.5µg)  
 6; Isoproturon (0.5µg)  
 7; Monolinuron (0.5µg)  
 8; Metabromuron (0.5µg)  
 9; Dimetufun (0.5µg)  
 10; Propazine (0.5µg)  
 11; Terbutylazine (0.5µg)  
 12; Linuron (0.5µg)  
 13; Chloroxuron (0.5µg)  
 14; Prometryn (0.5µg)  
 15; Terbutryn (0.5µg)  
 16; Chloro IPC (2.5µg)  
 17; Ethofumesate (5.0µg)




NACALAI TESQUE, INC  
 AP-0152

## • Triazine and Urea Herbicides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample:  
 1; Ethidimuron (1.0µg)  
 2; Metoxuron (1.0µg)  
 3; Cyanazine (1.0µg)  
 4; Methabenzthiazuron (1.0µg)  
 5; Chlorotoluron (1.0µg)  
 6; Atrazine (1.0µg)  
 7; Isoproturon (1.0µg)  
 8; Metabromuron (1.0µg)  
 9; Metazachlor (5.0µg)  
 10; Propazine (1.0µg)  
 11; Terbutylazine (1.0µg)  
 12; Linuron (1.0µg)



NACALAI TESQUE, INC  
 AP-0153


## 5) Food Additives

## • Natural Colorants (Chlorophyll)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphoric Acid = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.02AUFS

Sample: 1; Sodium Copper Chlorophyllin (30µg)  
 2; Sodium Iron Chlorophyllin (30µg)




NACALAI TESQUE, INC  
 AP-0095

## • Natural Colorants (Carotenoid)

**COSMOSIL Application Data**

Column: 5SL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Chloroform/Hexane = 1/9  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; β-Carotene (0.5µg)  
 2; all-trans-Retinol Acetate (1.5µg)




NACALAI TESQUE, INC  
 AP-0094

## • Synthetic Sweeteners

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile / 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 1.0AUFS

Sample: 1; Acesulfame K (1.0µg)  
 2; Saccharin (1.0µg)  
 3; Diketopiperazine (1.0µg)  
 4; Aspartame (1.0µg)




NACALAI TESQUE, INC  
 AP-0107

## • Synthetic Colorants

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Methanol/ 20mmol Phosphate buffer(pH7) = 10/90  
 B; Methanol/ 20mmol Phosphate buffer(pH7) = 80/20  
 B conc. 15→100% Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1; Tartrazine(黄色4号) (1.0µg)  
 2; Amaranth(赤色2号) (1.0µg)  
 3; Indigo Carmine(青色2号) (1.0µg)  
 4; New Coccine(赤色102号) (1.0µg)  
 5; Sunset Yellow FCF(黄色5号) (1.0µg)  
 6; Fast Green FCF(绿色3号) (3.0µg)  
 7; Brilliant Blue FCF(青色1号) (3.0µg)  
 8; Acid Red(赤色106号) (1.0µg)  
 9; Erythrosine(赤色3号) (1.0µg)  
 10; Phloxine(赤色104号) (1.0µg)  
 11; Rose Bengale(赤色105号) (1.0µg)




NACALAI TESQUE, INC  
 AP-0092

## • Synthetic Colorants

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; Methanol/ 20mmol Phosphate buffer(pH7) = 10/90  
 B; Methanol/ 20mmol Phosphate buffer(pH7) = 80/20  
 B conc. 15→100% Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1; Tartrazine(黄色4号) (1.0µg)  
 2; Amaranth(赤色2号) (1.0µg)  
 3; Indigo Carmine(青色2号) (1.0µg)  
 4; New Coccine(赤色102号) (1.0µg)  
 5; Sunset Yellow FCF(黄色5号) (1.0µg)  
 6; Fast Green FCF(绿色3号) (3.0µg)  
 7; Brilliant Blue FCF(青色1号) (3.0µg)  
 8; Acid Red(赤色106号) (1.0µg)  
 9; Erythrosine(赤色3号) (1.0µg)  
 10; Phloxine(赤色104号) (1.0µg)  
 11; Rose Bengale(赤色105号) (1.0µg)




NACALAI TESQUE, INC  
 AP-0093

## • Food Preservatives

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/Methanol/ 20mmol/l Phosphate buffer(pH4.0) = 20/10/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm

Sample: 1; Benzoic Acid (0.075µg)  
 2; Sorbic Acid (0.1µg)  
 3; Dehydroacetic Acid (0.175µg)



NACALAI TESQUE, INC  
 AP-0378


## 5) Food Additives

## ● Preservatives (Fungicides)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV230nm, 0.5AUFS

Sample: 1; Thiabendazole (0.3µg)  
 2; *o*-Phenylphenol (0.3µg)  
 3; Imazalil (2.1µg)  
 4; Diphenyl (0.3µg)




NACALAI TESQUE, INC  
 AP-0097

## ● Glycyrrhizic Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile / (Acetic Acid/H<sub>2</sub>O=1/15)  
 = 2/3  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample: 1; Glycyrrhizic Acid (5.0µg)  
 2; Propyl *p*-Hydroxybenzoate (1.0µg)




NACALAI TESQUE, INC  
 AP-0109

## ● Antioxidants

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.16AUFS

Sample: 1; 2-*tert*-Butylhydroquinone (0.56µg)  
 2; 3-*tert*-Butyl-4-hydroxyanisole (0.58µg)  
 3; 2,6-Di-*tert*-butyl-4-hydroxytoluene (1.85µg)




NACALAI TESQUE, INC  
 AP-0100

## ● Antioxidants

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Propyl Gallate (1.0µg)  
 2; Rthoxyquin (1.0µg)



NACALAI TESQUE, INC  
 AP-0101

## ● Repellents

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: Piperonyl Butoxide (10µg)




NACALAI TESQUE, INC  
 AP-0104

## ● Umami Seasonings (Nucleic Acids)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ 5mmol/l Tributylammonium bromide, 20mmol/l Phosphate buffer(pH7)  
 = 5/95  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm, 0.32AUFS

Sample: 1; Sodium 5'-Cytidylate (1.0µg)  
 2; Sodium 5'-Uridylate (1.0µg)  
 3; Sodium 5'-Guanylate (1.0µg)  
 4; Sodium 5'-Inosinate (1.0µg)



NACALAI TESQUE, INC  
 AP-0106


## 6) Vitamins

## ● Hydrosoluble Vitamins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile /  
 5mmol/l Sodium *I*-Hexanesulfonate,  
 20mmol/l Phosphoric Acid = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 0.16AUFS

Sample: 1; Vitamin C (0.23µg)  
 2; Niacin (0.057µg)  
 3; Nicotinamide (0.042µg)  
 4; Vitamin B<sub>6</sub> (0.040µg)  
 5; Riboflavin Phosphate Sodium Salt (0.19µg)  
 6; Vitamin B<sub>1</sub> Hydrochloride (0.19µg)  
 7; Folic Acid (0.084µg)  
 8; Vitamin B<sub>2</sub> (0.57µg)




NACALAI TESQUE, INC  
 AP-0059

## ● Hydrosoluble Vitamins

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Nicotinamide (0.125µg)  
 2; Pyridoxine(Vitamin B<sub>6</sub>) (0.25µg)  
 3; Riboflavin (Vitamin B<sub>2</sub>) (0.25µg)  
 4; Nicotinic Acid (0.125µg)  
 5; *D*-Pantothenic Acid (3.125µg)  
 6; *L*(+)-Ascorbic Acid (0.875µg)




NACALAI TESQUE, INC  
 AP-0315

## ● Water-soluble Vitamins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 20mmol/l Phosphate buffer(pH2.5)  
 B; Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 60/40  
 B conc. 0→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Vitamin B<sub>1</sub> Hydrochloride (0.17µg)  
 2; *L*(+)-Ascorbic Acid (0.33µg)  
 3; Nicotinic Acid (0.05µg)  
 4; Nicotinamide (0.05µg)  
 5; Pyridoxine Hydrochloride (0.27µg)  
 6; *D*-Pantothenic Acid Sodium Salt (2.01µg)  
 7; Vitamin B<sub>2</sub> (0.20µg)  
 8; Folic Acid (0.26µg)  
 9; *D*-Biotin [Vitamin H] (2.02µg)  
 10; Flavin Mononucleotide Sodium Salt (0.26µg)  
 11; Vitamin B<sub>2</sub> (Riboflavin) (0.13µg)




NACALAI TESQUE, INC  
 AP-1055

## ● Ascorbic Acids

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 100mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Sorbic Acid (0.05µg)  
 2; Isoascorbic Acid (0.3µg)  
 3; Ascorbic Acid (0.3µg)



NACALAI TESQUE, INC  
 AP-0318


## ● Energy Drink

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Sodium *I*-Hexanesulfonate,  
 20mmol/l Phosphoric Acid = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Citric Acid (10mg/ml)  
 2; Carnitine (20mg/ml)  
 3; Nicotinamide (0.2mg/ml)  
 4; Vitamin B<sub>6</sub> (0.2mg/ml)  
 5; Vitamin B<sub>1</sub> (0.2mg/ml)  
 6; Caffeine (0.2mg/ml)  
 7; Riboflavin Phosphate (0.2mg/ml)

Injection Vol. 1.0µl



NACALAI TESQUE, INC  
 AP-1048


## 6) Vitamins

## • Fruit Juice

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Ascorbic Acid (1.5µg)  
 2; Malic Acid (3.0µg)




NACALAI TESQUE, INC  
 AP-0313

## • Fat-soluble Vitamins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/Methanol = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.64AUFS

Sample: 1; Vitamin K<sub>3</sub> (1.71µg)  
 2; Vitamin A Acetate (0.44µg)  
 3; Vitamin D<sub>3</sub> (0.57µg)  
 4; Vitamin E (4.97µg)  
 5; Vitamin E Acetate (4.79µg)  
 6; Vitamin K<sub>1</sub> (1.68µg)




NACALAI TESQUE, INC  
 AP-0060

## • Vitamin D

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV265nm

Sample: 1; Vitamin D<sub>2</sub> (0.3µg)  
 2; Vitamin D<sub>3</sub> (0.1µg)




NACALAI TESQUE, INC  
 AP-1035

## • Vitamin E

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol  
 πNAP Methanol/ H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV295nm

Sample: 1; δ-Tocotrienol  
 2; γ-Tocotrienol  
 3; β-Tocotrienol  
 4; α-Tocotrienol  
 5; δ-Tocopherol  
 6; γ-Tocopherol  
 7; β-Tocopherol  
 8; α-Tocopherol



NACALAI TESQUE, INC  
 AP-1071


## 6) Vitamins

## • Tocopherols

**COSMOSIL Application Data**

Column: SSL-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: n-Hexane/2-Propanol = 99.5/0.5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.16AUFS

Sample: 1; α-Tocopherol (3.0µg)  
 2; β-Tocopherol (7.5µg)  
 3; γ-Tocopherol (7.5µg)  
 4; δ-Tocopherol (7.5µg)




NACALAI TESQUE, INC  
 AP-0061

## • Vitamin A

**COSMOSIL Application Data**

Column: SSL-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: n-Hexane/Chloroform = 9/1  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; β-Carotene (0.5µg)  
 2; Vitamin A (1.5µg)




NACALAI TESQUE, INC  
 AP-0062

## • Vitamin A Acid

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / 20mmol/l Phosphate buffer(pH2.5) = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV350nm

Sample: 1; 13-cis-Retinoic Acid (0.04µg)  
 2; all-trans-Retinoic Acid (0.04µg)



NACALAI TESQUE, INC  
 AP-1036


## 7) Metabolites

## • Androgens

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV235nm, 0.5AUFS

Sample: 1; 4-Androstene-3,17-dione (1.0µg)  
 2; Testosterone (1.0µg)  
 3; 17-Methyltestosterone (1.0µg)  
 4; Testosterone Propionate (1.0µg)




NACALAI TESQUE, INC  
 AP-0111

## • Catecholamines

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 20mmol/l NaH<sub>2</sub>PO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV270nm, 0.1AUFS

Sample: 1; Norepinephrine (1.0µg)  
 2; Epinephrine (1.2µg)  
 3; L-DOPA (1.2µg)  
 4; Dopamine (1.0µg)  
 5; 3,4-Dihydroxyphenylacetic Acid (1.0µg)  
 6; 3-Methoxytyramine (1.0µg)




NACALAI TESQUE, INC  
 AP-0120

## • Estrogens

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 35/65  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.2AUFS

Sample: 1; Estriol (1.5µg)  
 2; Estradiol (1.5µg)  
 3; Estrone (1.5µg)




NACALAI TESQUE, INC  
 AP-0114

## • Estradiols

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm

Sample: 1; 17-β-Estradiol (2.0µg)  
 2; 17-α-Estradiol (2.0µg)




NACALAI TESQUE, INC  
 AP-0243

## • Hippuric Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 15/85  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV225nm, 0.32AUFS

Sample: 1; Creatinine (0.14µg)  
 2; Mandelic Acid (1.89µg)  
 3; Hippuric Acid (0.61µg)  
 4; o-Methyl Hippuric Acid(0.52µg)




NACALAI TESQUE, INC  
 AP-0112

## • Methylhippuric Acids

**COSMOSIL Application Data**

Column: SPYE  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: N-(o-Toluoyl) glycine (3.06µg)  
 N-(m-Toluoyl) glycine (0.65µg)  
 N-(p-Toluoyl) glycine (0.33µg)



NACALAI TESQUE, INC  
 AP-0289


## 7) Metabolites

## • Urate Metabolites

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 20mmol/l Phosphoric Acid  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV280nm, 0.1AUFS

Sample: 1; Hypoxanthine (0.9µg)  
 2; Uric Acid (18.3µg)  
 3; Xanthine (9.0µg)  
 4; Allopurinol (1.8µg)




NACALAI TESQUE, INC  
 AP-0117

## • Prostaglandins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-40%Acetonitrile  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Prostaglandin I<sub>2</sub> (2.0µg)  
 2; Prostaglandin F<sub>2</sub> α (2.0µg)  
 3; Prostaglandin E<sub>2</sub> (2.0µg)  
 4; Prostaglandin D<sub>2</sub> (2.0µg)




NACALAI TESQUE, INC  
 AP-0247

## • Adrenal Cortical Hormones

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol/ H<sub>2</sub>O = 50/50  
 πNAP Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Triamcinolone  
 2; Cortisone  
 3; Prednisolone  
 4; Hydrocortisone  
 5; Betamethasone  
 6; Dexamethasone  
 7; Triamcinolone Acetonide  
 8; Cortisone Acetate  
 9; Prednisolone Acetate  
 10; Fluocinolone Acetonide  
 11; Fluorometholone




(each 0.23µg)  
 NACALAI TESQUE, INC  
 AP-1037

## • Adrenal Cortical Hormones

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV240nm

Sample: 1; Prednisolone (0.1µg)  
 2; Hydrocortisone (0.1µg)  
 3; Cortisone (0.1µg)  
 4; Betamethasone (0.1µg)  
 5; Dexamethasone (0.1µg)



NACALAI TESQUE, INC  
 AP-0233


## 8) Carbohydrates

## ● Oligosaccharides

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Sucrose (10µg)  
 2; Maltose (10µg)  
 3; Trehalose (10µg)  
 4; Raffinose (10µg)



NACALAI TESQUE, INC


AP-0335

## ● Cyclodextrins

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 65/35  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; α-Cyclodextrin (10µg)  
 2; β-Cyclodextrin (10µg)  
 3; γ-Cyclodextrin (10µg)



NACALAI TESQUE, INC


AP-0336

## ● Xylooligosaccharides

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Xylooligosaccharides (50µg)  
 (Xylobiose, Xylotriose etc.)



NACALAI TESQUE, INC


AP-0321

## ● Isomaltooligosaccharides

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Isomaltooligosaccharides (50µg)  
 (Isomaltose, Isomaltotriose, Panose etc.)



NACALAI TESQUE, INC


AP-0322

## ● Anticarious Foods Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; meso-Erythritol (10µg)  
 2; Xylitol (10µg)  
 3; Palatinose (10µg)  
 4; Maltitol (10µg)



NACALAI TESQUE, INC


AP-0323

## ● Infusion Solution Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 85/15  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: RI

Sample: 1; Xylitol (10µg)  
 2; Fructose (10µg)  
 3; Glucose (10µg)



NACALAI TESQUE, INC

AP-0328

## 8) Carbohydrates


## ● Energy Drink

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 85/15  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: ELSD

Sample: 1; Fructose (1.0mg/ml)  
 2; Glucose (1.0mg/ml)  
 3; Glucitol(Sorbitol) (1.0mg/ml)  
 4; Sucrose (1.0mg/ml)  
 5; Inositol (1.0mg/ml)

Injection Vol. 5.0µl



NACALAI TESQUE, INC


AP-1049

## ● Chocolate Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Sucrose (10µg)  
 2; Lactose (10µg)



NACALAI TESQUE, INC


AP-0324

## ● Gum Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Xylitol (10µg)  
 2; Mannitol (10µg)  
 3; Maltitol (10µg)  
 4; Palatinol (10µg)  
 5; Maltotriitol (10µg)  
 6; Maltotetraitol (10µg)



NACALAI TESQUE, INC


AP-0325

## ● Cold Beverage Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Fructose (10µg)  
 2; Glucose (10µg)  
 3; Sucrose (10µg)  
 4; Maltose (10µg)



NACALAI TESQUE, INC


AP-0326

## ● Sports Drink Components

**COSMOSIL Application Data**

Column: Sugar-D  
 Column size: (4.6mmI.D.-250mm) × 2  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 75/25  
 Flow rate: 1.0 ml/min  
 Temperature: 50°C  
 Detection: RI

Sample: 1; meso-Erythritol (10µg)  
 2; Fructose (10µg)  
 3; Glucitol (10µg)  
 4; Glucose (10µg)  
 5; Sucrose (10µg)  
 6; Maltose (10µg)  
 7; Maltotriitol (10µg)  
 8; Maltotriose (10µg)



NACALAI TESQUE, INC

AP-0327




## 8) Carbohydrates

## ● PA-Glucose Oligomer

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 20mmol/l Acetate buffer(pH3.3)  
 B: 20mmol/l Acetate buffer(pH3.3)+0.5%Butanol  
 B(0→100%) 45min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: FLS at Ex.320nm  
 Em. 400nm

Sample:  
 PA-Glucose Oligomer(DP=3-22)



NACALAI TESQUE, INC


AP-0379

## ● Phosphorylated Sugars

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate  
 buffer(pH7.0) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample:  
 1; D-Fructose-6-phosphate (50µg)  
 2; D-Glucose-6-phosphate (50µg)



NACALAI TESQUE, INC


AP-0317

## ● Phosphorylated Sugars

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 20mmol/l Phosphate  
 buffer(pH7.0) = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample:  
 1; Glucose (10µg)  
 2; α-D-Glucose-1-phosphate (50µg)  
 3; D-Glucose-6-phosphate (50µg)



NACALAI TESQUE, INC


AP-0314

## ● Hyaluronic Acid

**COSMOSIL Application Data**

Column: CNT-1000  
 Column size: 7.5mmI.D.-300mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7),100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Hyaluronic Acid (100µg)



NACALAI TESQUE, INC


AP-1081

## ● Enzyme digests of Chondroitin Sulfate A

**COSMOSIL Application Data**

Column: Cholesterol  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 1mmol/l Tetrabutylammonium bisulfate/H<sub>2</sub>O  
 B: 1mmol/l Tetrabutylammonium bisulfate  
 -Acetonitrile/H<sub>2</sub>O=67/33  
 B conc. 20%→65%(7min)→65%(12min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: Chondroitin Sulfate A  
 Chondroitinase AC-II digested



NACALAI TESQUE, INC


AP-1082

## ● Enzyme digests of Chondroitin Sulfate A

**COSMOSIL Application Data**

Column: Cholesterol  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 1mmol/l Tetrabutylammonium bisulfate/H<sub>2</sub>O  
 B: 1mmol/l Tetrabutylammonium bisulfate  
 -Acetonitrile/H<sub>2</sub>O=67/33  
 B conc. 20%→65%(7min)→65%(12min)  
 Flow rate: 1.0 ml/min  
 Temperature: 40°C  
 Detection: UV240nm

Sample: Chondroitin Sulfate C  
 Chondroitinase AC-II digested



NACALAI TESQUE, INC

AP-1083


## 9) Lipids

## ● Fatty Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Oleic Acid (3.0µg)  
 2; Elaidic Acid (3.0µg)  
 3; Stearic Acid (3.0µg)



NACALAI TESQUE, INC


AP-1038

## ● Fatty Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; cis-Vaccenic Acid (3.0µg)  
 2; trans-Vaccenic Acid (3.0µg)



NACALAI TESQUE, INC


AP-1039

## ● Methylated Fatty Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Methyl Linolenate  
 Methyl Linoleate  
 Methyl Oleate  
 Methyl Margarate  
 Methyl Stearate  
 Methyl n-Nonadecanoate  
 (each 10µg)



NACALAI TESQUE, INC

AP-1040


## 9) Lipids

## ● Fatty Acid Derivatives

**COSMOSIL Application Data**

Column: Cholester  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Linolenic Acid *p*-Bromophenacyl Ester  
 2;  $\gamma$ -Linolenic Acid *p*-Bromophenacyl Ester



## ● Phosphatides

**COSMOSIL Application Data**

Column: SSL-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.05%TFA-Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: *L*- $\alpha$ -Phosphatidyl Choline (10 $\mu$ g)




## ● Phosphatides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.05%TFA-10%THF/90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: *L*- $\alpha$ -Phosphatidyl Choline (2.0 $\mu$ g)




## ● Phosphatides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.05%TFA-10%THF/90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: *L*- $\alpha$ -Phosphatidyl-L-serine




## ● Phosphatides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.05%TFA-10%THF/90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: *L*- $\alpha$ -Phosphatidylethanolamine (8.0 $\mu$ g)




## ● Phosphatides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: 0.05%TFA-10%THF/90%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: *L*- $\alpha$ -Phosphatidylinositol Sodium Salt



## 10) Nucleic Acid Related Substances


## ● Nucleobases and Nucleosides

**COSMOSIL Application Data**

Column:  $\pi$ NAP  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 20mmol/l Phosphate buffer(pH7.0)  
 B; Methanol/ 20mmol/l Phosphate buffer(pH7.0) = 10/90  
 B conc. 0→100% 15min Linear gradient

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm

Sample: 1; Cytosine (0.16 $\mu$ g)  
 2; Uracil (0.06 $\mu$ g)  
 3; Cytidine (0.32 $\mu$ g)  
 4; Uridine (0.24 $\mu$ g)  
 5; Guanine (0.16 $\mu$ g)  
 6; Thymine (0.16 $\mu$ g)  
 7; Adenine (0.08 $\mu$ g)  
 8; Guanosine (0.24 $\mu$ g)  
 9; Thymidine (0.24 $\mu$ g)  
 10; Adenosine (0.16 $\mu$ g)




## ● Nucleobases, Nucleosides and Nucleotides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm

Sample: 1; 5'-ATP (0.2 $\mu$ g)  
 2; 5'-ADP (0.2 $\mu$ g)  
 3; 5'-AMP (0.2 $\mu$ g)  
 4; Adenine (0.05 $\mu$ g)  
 5; Adenosine (0.1 $\mu$ g)




## ● Nucleotides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Methanol/ 5mmol/l Tetra-*n*-butylammonium Phosphate, 20mmol/l Phosphate buffer (pH2.5) = 20/80

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm

Sample: 1; 5'-AMP (0.33 $\mu$ g)  
 2; 5'-ADP (0.33 $\mu$ g)  
 3; 5'-ATP (0.33 $\mu$ g)




## ● Nucleotides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.5AUFS

Sample: 1; *d*CMP (1.0 $\mu$ g)  
 2; *d*GMP (1.2 $\mu$ g)  
 3; *d*AMP (1.0 $\mu$ g)




## ● DNA Oligomer

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 100mmol/l Ammonium Acetate  
 B; Acetonitrile  
 B conc. 9→11% 15min Linear gradient

Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm

Sample: pd(T)<sub>12-18</sub>




## ● DNA Ladder

**COSMOSIL Application Data**

Column: CNT-1000 + CNT-2000  
 Column size: 7.5mmI.D.-300mm  $\times$  2  
 Mobile phase: 20mmol/l Phosphate buffer(pH7), 100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV260nm

Sample: DNA (100b.p. Ladder)




11) Amino Acids, Peptides and Proteins

• Dansyl Amino Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; Acetonitrile/ 20mmol Phosphate buffer(pH7.0) =10/90  
 B; Acetonitrile/ 20mmol Phosphate buffer(pH7.0) =40/60  
 B conc. 0→100% 30min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: Ex.365nm Em.530nm RANGE 2

Sample: Dansyl Amino Acids (1.0µg each)



NACALAI TESQUE, INC  
AP-0003


• PTH-Amino Acids

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; Acetonitrile/ 20mmol/l Acetate buffer(pH4.8) =10/90  
 B; Acetonitrile/ 20mmol/l Acetate buffer(pH4.8) =50/50  
 Step wise gradient

Flow rate: 1.0 ml/min  
 Temperature: 60°C  
 Detection: UV270nm, 0.16AUFS  
 Sample: PTH-Asp (0.4mg/ml)  
 PTH-His (0.3mg/ml)  
 PTH-Arg (0.5mg/ml)  
 PTH-Ser (0.3mg/ml)  
 PTH-Trp (0.25mg/ml)  
 PTH-Lys (0.35mg/ml)  
 Others (0.2mg/ml)

Injection 2.0µl



NACALAI TESQUE, INC  
AP-0004


11) Amino Acids, Peptides and Proteins

• Amino-Acid Derivatives (Diastereomer)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Boc-(L)Phe-(L)Phe-OMe (7.5µg)  
 2; Boc-(D)Phe-(L)Phe-OMe (7.5µg)




NACALAI TESQUE, INC  
AP-0241

• The umami of Vegetables

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 5mmol/l Sodium L-Hexanesulfonate, 20mmol/l Phosphoric Acid  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Aspartic Acid (10µg)  
 2; Glutamic Acid (10µg)  
 3; 4-Amino-n-butyric Acid [GABA] (20µg)



NACALAI TESQUE, INC  
AP-1015


• Energy Drink

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; Taurine (10mg/ml)

Standard (10mg/ml)




NACALAI TESQUE, INC  
AP-1050

• Peptides

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 10→40% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Oxytocin (0.8µg)  
 2; Angiotensin II(Human) (0.8µg)  
 3; Angiotensin I(Human) (0.8µg)  
 4; Substance P (0.8µg)




NACALAI TESQUE, INC  
AP-0350

• Peptides

**COSMOSIL Application Data**

Column: Cholesterol  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 10→40% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: 1; Oxytocin (0.8µg)  
 2; Angiotensin II(Human) (0.8µg)  
 3; Angiotensin I(Human) (0.8µg)  
 4; Substance P (0.8µg)



NACALAI TESQUE, INC  
AP-0238

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index


## 11) Amino Acids, Peptides and Proteins

## ● Peptides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: A; 0.05%TFA-10%Acetonitrile  
 B; 0.05%TFA-30%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample:  
 1; Angiotensin II, [Sar<sup>1</sup>,Thr<sup>8</sup>] (0.49µg)  
 2; Angiotensin II, [Sar<sup>1</sup>,Ala<sup>8</sup>] (0.49µg)  
 3; Angiotensin II, Des-Asp<sup>1</sup>-[Ile<sup>6</sup>] (0.49µg)  
 4; Angiotensin II, [Sar<sup>1</sup>,Ile<sup>6</sup>] (0.49µg)  
 5; Angiotensin II, [Asn<sup>1</sup>,Val<sup>2</sup>] (0.49µg)  
 6; Angiotensin II, [Val<sup>2</sup>] (0.49µg)  
 7; Angiotensin II (Human) (0.49µg)




NACALAI TESQUE, INC  
 AP-1005

## ● Peptides

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mmI.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate  
 buffer(pH7.0) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample:  
 3; Angiotensin II, Des-Asp<sup>1</sup>-[Ile<sup>6</sup>] (0.5µg)  
 4; Angiotensin II, [Sar<sup>1</sup>,Ile<sup>6</sup>] (0.5µg)  
 5; Angiotensin II, [Asn<sup>1</sup>,Val<sup>2</sup>] (0.5µg)  
 6; Angiotensin II, [Val<sup>2</sup>] (0.5µg)




NACALAI TESQUE, INC  
 AP-0316

## ● Peptides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-22%Acetonitrile  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.16AUFS

Sample:  
 1; Bradykinin (1.0µg)  
 2; Angiotensin II (1.0µg)  
 3; Neurotensin (1.0µg)  
 4; Bombesin (1.0µg)  
 5; Substance P (1.0µg)




NACALAI TESQUE, INC  
 AP-0009

## ● Peptides

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: 0.05%TFA-22%Acetonitrile  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.16AUFS

Sample:  
 1; Met-Enkephalin (0.5µg)  
 2; [Ala<sup>2</sup>]-Met-Enkephalin (0.5µg)  
 3; Leu-Enkephalin (0.5µg)  
 4; [Ala<sup>2</sup>]-Leu-Enkephalin (0.5µg)




NACALAI TESQUE, INC  
 AP-0005

## ● Peptide Mappings

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 5mmol/l TFA-H<sub>2</sub>O  
 B; 5mmol/l TFA-60%Acetonitrile  
 B conc. 0→100% 30min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm, 0.32AUFS

Sample: Cytochrome C  
 Lysyl Endopeptidase digested  
 (2.5µl)




NACALAI TESQUE, INC  
 AP-0002

## ● Semi-purified Myosin

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Myosin (20µg)



NACALAI TESQUE, INC  
 AP-0346

## 11) Amino Acids, Peptides and Proteins

## ● Milk Protein


**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 20→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Milk

Sample Preparation:  
 •Ultracentrifuged at 90,000 g for 1 hr.  
 •Clear supernatant solution was injected.

Injection vol.: 1.0µl



NACALAI TESQUE, INC  
 AP-1079

## ● Soymilk Protein


**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A; 0.05%TFA-H<sub>2</sub>O  
 B; 0.05%TFA-Acetonitrile  
 B conc. 20→80% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Soybean Milk

Sample Preparation:  
 •Ultracentrifuged at 90,000 g for 1 hr.  
 •Clear supernatant solution was injected.

Injection vol.: 1.0µl




NACALAI TESQUE, INC  
 AP-1080

## ● Bacteria-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Choline Oxidase (6.0µg)  
 2; α-Amylase (3.0µg)  
 3; Glucose Oxidase (6.0µg)  
 4; Thermolysin (9.0µg)




NACALAI TESQUE, INC  
 AP-0337

## ● Bacteria-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Actinase E (13.4µg)  
 2; Alcohol Dehydrogenase (6.6µg)




NACALAI TESQUE, INC  
 AP-0338

## ● Human-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Transferrin (1.5µg)  
 2; Albumin(Human) (2.0µg)  
 3; Carbonic Anhydrase (1.0µg)




NACALAI TESQUE, INC  
 AP-0339

## ● Bovine-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample:  
 1; Albumin(Bovine) (2.0µg)  
 2; L-Glutamic Dehydrogenase (1.5µg)  
 3; Carbonic Anhydrase (1.5µg)



NACALAI TESQUE, INC  
 AP-0340




## 11) Amino Acids, Peptides and Proteins

## ● Bovine-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Fibrinogen (4.0µg)  
 2; Catalase (2.0µg)  
 3; Thyroglobulin (6.0µg)



NACALAI TESQUE, INC


AP-0341

## ● Bovine Spleen-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Ribonuclease A (1.0µg)  
 2; α-Chymotrypsinogen A (1.0µg)  
 3; Deoxyribonuclease I (3.0µg)



NACALAI TESQUE, INC


AP-0342

## ● Egg-derived Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Lysozyme (1.0µg)  
 2; Conalbumin (2.0µg)  
 3; Albumin(Ovalbumin)(1.8µg)



NACALAI TESQUE, INC


AP-0343

## ● Other Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Cytochrome C (1.0µg)  
 2; Albumin (Goat) (2.0µg)  
 3; Myoglobin (1.0µg)  
 4; Concanavalin A (0.5µg)



NACALAI TESQUE, INC


AP-0344

## ● Other Proteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: 1; Albumin(Rabbit) (3.2µg)  
 2; Peroxidase (3.2µg)



NACALAI TESQUE, INC


AP-0345

## ● Pyruvate Kinase

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Pyruvate Kinase (4.0µg)



NACALAI TESQUE, INC

AP-0348


## 11) Amino Acids, Peptides and Proteins

## ● Semi-purified Diaphorase

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A:0.05%TFA-20%Acetonitrile  
 B:0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 10min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.256AUFS

Sample: Diaphorase (6.0µg)



NACALAI TESQUE, INC


AP-0347

## ● Glycoproteins

**COSMOSIL Application Data**

Column: Protein-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 0.05%TFA-20%Acetonitrile  
 B: 0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm

Sample: Concanavalin A (1.0µg)



NACALAI TESQUE, INC


AP-0352

## ● Proteins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 0.05%TFA-20%Acetonitrile  
 B: 0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)



NACALAI TESQUE, INC


AP-0017

## ● Proteins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 0.05%TFA-20%Acetonitrile  
 B: 0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)



NACALAI TESQUE, INC


AP-0018

## ● Proteins

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 0.05%TFA-20%Acetonitrile  
 B: 0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)



NACALAI TESQUE, INC


AP-0019

## ● Proteins

**COSMOSIL Application Data**

Column: 5Ph-AR-300  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: 0.05%TFA-20%Acetonitrile  
 B: 0.05%TFA-60%Acetonitrile  
 B conc. 0→100% 20min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.64AUFS

Sample: 1; Insulin (1.5µg)  
 2; Trypsinogen (6.0µg)  
 3; Transferrin (4.0µg)  
 4; Trypsin Inhibitor, Soybean (5.0µg)  
 5; α-Chymotrypsinogen A (4.0µg)  
 6; Carbonic Anhydrase (3.0µg)



NACALAI TESQUE, INC

AP-0020




## 11) Amino Acids, Peptides and Proteins

## • Human Serum

**COSMOSIL Application Data**

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Human Serum




NACALAI TESQUE, INC  
 AP-0382

## • Human Serum

**COSMOSIL Application Data**

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Human Serum




NACALAI TESQUE, INC  
 AP-0383

## • Bovine Serum

**COSMOSIL Application Data**

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Bovine Serum




NACALAI TESQUE, INC  
 AP-0386

## • Bovine Serum

**COSMOSIL Application Data**

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Bovine Serum




NACALAI TESQUE, INC  
 AP-0387

## • Egg White

**COSMOSIL Application Data**

Column: 5Diol-120-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Egg White




NACALAI TESQUE, INC  
 AP-0384

## • Egg White

**COSMOSIL Application Data**

Column: 5Diol-300-II  
 Column size: 7.5mmI.D.-600mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH 7.0)  
 +100mmol/l Na<sub>2</sub>SO<sub>4</sub>  
 Flow rate: 1.0 ml/min  
 Temperature: Room temperature  
 Detection: UV280nm

Sample: Egg White



NACALAI TESQUE, INC  
 AP-0385


## 12) The others

## • Methoxyphenols

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *o*-Methoxyphenol (3.3µg)  
*m*-Methoxyphenol (3.3µg)  
*p*-Methoxyphenol (3.3µg)




NACALAI TESQUE, INC  
 AP-1041

## • Terphenyls

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *o*-Terphenyl (0.15µg)  
*m*-Terphenyl (0.05µg)  
*p*-Terphenyl (0.075µg)




NACALAI TESQUE, INC  
 AP-1042

## • Chlorophenols

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: *o*-Chlorophenol (2.0µg)  
*m*-Chlorophenol (2.0µg)  
*p*-Chlorophenol (4.0µg)



NACALAI TESQUE, INC  
 AP-1043


## 12) The others

## • Stilbenes

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; *cis*-Stilbene (0.10µg)  
 2; *trans*-Stilbene (0.20µg)




AP-1045

## • Halogenated Benzenes

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Toluene (3.3µg)  
 2; Chlorobenzene (3.3µg)  
 3; Bromobenzene (3.3µg)



NACALAI TESQUE, INC


AP-1046

## • Fluorinated Benzenes

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II Methanol/ H<sub>2</sub>O = 50/50  
 πNAP Methanol/ H<sub>2</sub>O = 50/50  
 NPE Methanol/ H<sub>2</sub>O = 40/60  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Benzene (6.0µg)  
 2; Fluorobenzene (2.0µg)  
 3; 1,4-Difluorobenzene (2.0µg)  
 4; 1,3-Difluorobenzene (2.0µg)  
 5; 1,2-Difluorobenzene (2.0µg)



NACALAI TESQUE, INC

AP-1054


## 12) The others

## • Pyridine and Phenol

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Pyridine (0.1µg)  
 2; Phenol (0.1µg)



NACALAI TESQUE, INC


AP-0239

## • Pyridine and Phenol

**COSMOSIL Application Data**

Column: Cholester  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Pyridine (0.1µg)  
 2; Phenol (0.1µg)



NACALAI TESQUE, INC


AP-0240

## • Uracil and Uridine

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Uracil (0.1µg)  
 2; Uridine (0.2µg)



NACALAI TESQUE, INC


AP-0299

## • Glycerol

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Trimethylene Glycol (20µg)  
 2; Glycerol (20µg)



NACALAI TESQUE, INC


AP-0301

## • Ethylene Glycol

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Trimethylene Glycol (20µg)  
 2; Ethylene Glycol (20µg)



NACALAI TESQUE, INC


AP-0303

## • Diethylene Glycol

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile / H<sub>2</sub>O = 95/5  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: 1; Diethylene Glycol (20µg)  
 2; Glycerol (20µg)



NACALAI TESQUE, INC

AP-0310


## 12) The others

## • Oxalic Acid

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 10mmol/l Phosphate buffer(pH7.0) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: 1; Oxamic Acid (0.2µg)  
 2; Oxalic Acid (1.0µg)




NACALAI TESQUE, INC  
 AP-0307

## • Hydrophilic Compounds (Ionicity)

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/ 50mmol/l Ammonium Acetate = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: ELSD

Sample: 1; meso-Erythritol (3.0µg)  
 2; Tris(hydroxymethyl)aminomethane (3.0µg)  
 3; Glyceric Acid (6.0µg)




NACALAI TESQUE, INC  
 AP-0305

## • Fluorine Compounds

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 0.1%TFA-95%Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: Perfluorotetradecanoic Acid (50µg)




NACALAI TESQUE, INC  
 AP-0381

## • Urea

**COSMOSIL Application Data**

Column: HILIC  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV210nm

Sample: Urea (20µg)




NACALAI TESQUE, INC  
 AP-1004

## • Acid Compounds

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ 0.1% Phosphoric Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Uracil (0.05µg)  
 2; Benzoic Acid (1.0µg)  
 3; o-Toluic Acid (1.0µg)  
 4; p-Ethylbenzoic Acid (0.2µg)  
 5; Benzene




NACALAI TESQUE, INC  
 AP-0164

## • Basic Compounds

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH3) = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.2AUFS

Sample: 1; 2-Ethylpyridine (0.4µg)  
 2; N,N-Dimethylaniline (0.6µg)



NACALAI TESQUE, INC  
 AP-0166


## 12) The others

## • Monosubstituted Benzenes (20 samples)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.16AUFS

Sample: 1; Benzamide (0.49µg) 14; Chlorobenzene (7.05µg)  
 2; Aniline (0.4µg) 15; Toluene (5.84µg)  
 3; Phenol (0.67µg) 16; Bromobenzene (15.37µg)  
 4; Benzonitrile (0.83µg) 17; Iodobenzene (3.66µg)  
 5; Acetophenone (0.04µg) 18; Ethylbenzene (6.87µg)  
 6; Styrene oxide (1.1µg) 19; n-Propylbenzene (14.1µg)  
 7; Nitrobenzene (0.06µg) 20; n-Butylbenzene (15.93µg)  
 8; Methylbenzoate (0.62µg)  
 9; Anisole (0.79µg)  
 10; Fluorobenzene (0.65µg)  
 11; Benzene (0.79µg)  
 12; N,N-Dimethylamine (0.15µg)  
 13; Thioanisole (0.12µg)




NACALAI TESQUE, INC  
 AP-0154

## • Furans

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 30/70  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV220nm, 0.2AUFS

Sample: 1; Furfuryl alcohol (0.13µg)  
 2; Furfural (0.25µg)  
 3; Furan (0.23µg)




NACALAI TESQUE, INC  
 AP-0155

## • Phenolphthalein

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.1AUFS

Sample: Phenolphthalein (0.6µg)




NACALAI TESQUE, INC  
 AP-0156

## • Anilines

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile / 0.1% Phosphoric Acid = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Aniline  
 2; p-Nitroaniline  
 3; 2,4-Dinitroaniline




NACALAI TESQUE, INC  
 AP-0160

## • Esters and others

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Acetonitrile/ H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.32AUFS

Sample: 1; Uracil  
 2; Pyridine  
 3; Phenol  
 4; Ethyl p-Hydroxybenzoate  
 5; Propyl p-Hydroxybenzoate  
 6; Methyl Salicylate  
 7; Toluene  
 8; Ethyl Salicylate




NACALAI TESQUE, INC  
 AP-0163

## • Aminophenols

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-PAQ  
 Column size: 4.6mm I.D.-250mm  
 Mobile phase: 20mmol/l Phosphate buffer(pH7.0)  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: p-Aminophenol (0.33µg)  
 m-Aminophenol (1.00µg)  
 o-Aminophenol (1.67µg)



NACALAI TESQUE, INC  
 AP-1059


## 12) The others

## ● Polyaromatic Compounds

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-AR-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: A: Methanol / H<sub>2</sub>O = 70/30  
 B: Methanol  
 B conc. 0→100%  
 3→15min Linear gradient  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm, 0.64AUFS

Sample:  
 1; Naphthalene  
 2; Acenaphthylene  
 3; Acenaphthene  
 4; Fluorene  
 5; Phenanthrene  
 6; Anthracene  
 7; Fluoranthene  
 8; Pyrene  
 9; Benz[*a*]anthracene  
 10; Chrysene  
 11; Benz[*b*]fluoranthene  
 12; Benz[*a*]pyrene  
 13; Benz[*a*]pyrene  
 14; Dibenz[*a,h*]anthracene  
 15; Benz[*g,h,i*]perylene  
 16; Indeno[1,2,3-*c,d*]pyrene




NACALAI TESQUE, INC  
 AP-0161

## ● Benzylpyridines

**COSMOSIL Application Data**

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l KH<sub>2</sub>PO<sub>4</sub> = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; 2-Benzylpyridine  
 2; 3-Benzylpyridine  
 3; 4-Benzylpyridine




NACALAI TESQUE, INC  
 AP-0182

## ● Carvylacetate

**COSMOSIL Application Data**

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV215nm

Sample: Carvylacetate (*cis*, *trans*)





NACALAI TESQUE, INC  
 AP-0184

## ● Diastereomers

**COSMOSIL Application Data**

Column: 5PYE  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; *Threo* form  
 2; *Erythro* form




NACALAI TESQUE, INC  
 AP-0186

## ● Surfactant

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran /Methanol = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Polyoxyethylene hexadecyl ether (500µg)




NACALAI TESQUE, INC  
 AP-0354

## ● Surfactant

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Tetrahydrofuran /Methanol = 10/90  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: RI

Sample: Polyoxyethylene Lauryl Ether (500µg)



NACALAI TESQUE, INC  
 AP-0364


## 12) The others

## ● Aromatic Compounds

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Diphenylmethane (5.80µg)  
 2; Fluorene (0.13µg)




NACALAI TESQUE, INC  
 AP-0358

## ● Aromatic Compounds

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Acetonitrile/H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Dibenzosuberone (0.25µg)  
 2; Dibenzosuberone (0.10µg)




NACALAI TESQUE, INC  
 AP-0360

## ● Aromatic Compounds

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Propylbenzene (27.9µg)  
 2; Allylbenzene (0.45µg)




NACALAI TESQUE, INC  
 AP-0362

## ● Aromatic Compounds

**COSMOSIL Application Data**

Column: 5PBB-R  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; 1,1'-Binaphthyl (0.60µg)  
 2; Perylene (0.72µg)




NACALAI TESQUE, INC  
 AP-0356

## ● Coloring agent (CBB)

**COSMOSIL Application Data**

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mmI.D.-150mm  
 Mobile phase: Methanol/ 20mmol/l Phosphate buffer(pH2.5) = 70/30  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: Coomassie Brilliant Blue G-250 (3.0µg)



NACALAI TESQUE, INC  
 AP-1061



## 4. Reference list

No.	Column	TITLE	AUTHOR	JOURNAL	YEAR	VOL. (ISSUE)	PAGE
1	Cholester	Improved separation of furocoumarins of essential oils by supercritical fluid chromatography	C. Desmorteux, M. Rothaupt, C. West and E. Lesellier	Journal of Chromatography A	2009	1216(42)	7088-7095
2	Cholester	Preparative Separation of cis- and trans-Isomers of Unsaturated Fatty Acid Methyl Esters Contained in Edible Oils by Reversed-Phase High-Performance Liquid Chromatography.	Tsuzuki W, Ushida K.	Lipids.	2009	44(4)	373-379
3	Cholester 5C <sub>18</sub> -AR-II	Isolation and Structure of a Galactocerebroside from the Sea Cucumber <i>Bohadschia argus</i>	Yuriko Ikeda, Masanori Inagaki, Koji Yamada, Xiao Wen Zhang, Bo Zhang, Tomofumi Miyamoto and Ryuichi Higuchi	CHEMICAL & PHARMACEUTICAL BULLETIN	2009	57(3)	315-317
4	Cholester	Influence of Cold Hardening on Chlorophyll and Carotenoid in <i>Chlorella vulgaris</i>	Watanabe Yuta, Yamada Naotaka, Machida Takeshi, Honjoh Ken-ichi, Kuwano Eiichi	Journal of the Faculty of Agriculture, Kyushu University	2009	54(1)	195-200
5	Cholester	Pyripyropenes, Fungal Sesquiterpenes Conjugated with $\alpha$ -Pyrone and Pyridine Moieties, Exhibits Anti-angiogenic Activity against Human Umbilical Vein Endothelial Cells	Asami Hayashi, Masayoshi Arai, Mayumi Fujita and Motomasa Kobayashi	Biological & Pharmaceutical Bulletin	2009	32(7)	1261-1265
6	Cholester	Practical Electrochemical Iodination of Aromatic Compounds	Kazuhide Kataoka, Yuji Hagiwara, Koji Midorikawa, Seiji Suga and Jun-ichi Yoshida	Org. Process Res. Dev.	2008	12(6)	1130-1136
7	Cholester	De novo synthesis of (Z)- and (E)-7-hexadecenyloctanoic acids by a selective lignin-degrading fungus, <i>Ceriporiopsis subvermispora</i>	Hiroshi Nishimura, Saeko Tsuda, Hito Shimizu, Yasunori Ohashi, Takahito Watanabe, Yoichi Honda and Takashi Watanabe	Phytochemistry	2008	69(14)	2593-2602
8	Cholester	Kadsuracoccinic Acids A-C, Ring-A seco-Lanostane Triterpenes from <i>Kadsura coccinea</i> and Their Effects on Embryonic Cell Division of <i>Xenopus laevis</i>	Heran Li, Liyan Wang, Syohei Miyata and Susumu Kitanaka	J. Nat. Prod.	2008	71(4)	739-741
9	Cholester	Peptide Thioester Synthesis via an Auxiliary-Mediated N-S Acyl Shift Reaction in Solution	Ken'ichiro Nakamura, Hiroaki Mori, Toru Kawakami, Hironobu Hojo, Yoshiaki Nakahara and Saburo Aimoto	International Journal of Peptide Research and Therapeutics	2007	13(1-2)	191-202
10	Cholester	Metachromins L-Q, new sesquiterpenoid quinones with an amino acid residue from sponge <i>Spongia</i> sp.	Yohei Takahashi, Takaaki Kubota, Jane Fromont and Jun'ichi Kobayashi	Tetrahedron	2007	63(36)	8770-8773
11	Cholester	Metachromins R-T, New Sesquiterpenoids from Marine Sponge <i>Spongia</i> sp.	Yohei Takahashi, Mika Yamada, Takaaki Kubota, Jane Fromont and Jun'ichi Kobayashi	CHEMICAL & PHARMACEUTICAL BULLETIN	2007	55(No.12)	1731-1733
12	$\pi$ NAP PYE PBB	Possibility of predicting separations in supercritical fluid chromatography with the solvation parameter model	C. West, J. Ogden and E. Lesellier	Journal of Chromatography A	2009	1216 (29)	5600-5607
13	PYE	Possibility of predicting separations in supercritical fluid chromatography with the solvation parameter model	C. West, J. Ogden and E. Lesellier	Journal of Chromatography A	2009	1216 (29)	5600-5607
14	PYE	Calculations and assignments of endohedral helium-3 chemical shifts of open-cage fullerenes and higher fullerenes	Guan-Wu Wang and Ping Wu	Theoretical Chemistry Accounts	2009	123 (5-6)	375-381
15	PYE	Effects of $\pi$ - $\pi$ Interactions on the Separation of PAHs on Phenyl-Type Stationary Phases	Paul G. Stevenson; Sindy Kayillo; Gary R. Dennis; R. Andrew Shalliker	Journal of Liquid Chromatography & Related Technologies	2008	31(3)	324-347
16	PYE	Organochlorine contaminants in endangered Steller sea lion pups ( <i>Eumetopias jubatus</i> ) from western Alaska and the Russian Far East	Matthew J. Myers, Gina M. Ylitalo, Margaret M. Krahn, Daryle Boyd, Don Calkins, Vladimir Burkanov and Shannon Atkinson	Science of The Total Environment	2008	396 (1)	60-69
17	PYE	Spatial extent, magnitude, and patterns of persistent organochlorine pollutants in Pacific herring ( <i>Clupea pallasii</i> ) populations in the Puget Sound (USA) and Strait of Georgia (Canada)	James E. West, Sandra M. O'Neill, and Gina M. Ylitalo	Science of The Total Environment	2008	394 (2-3)	369-378
18	PYE	Polychlorinated dibenzo-p-dioxins and dibenzofurans (PCDD/Fs) in marine mussels from French coasts: Levels, patterns and temporal trends from 1981 to 2005	C. Munsch, N. Guiot, K. Héas-Moisan, C. Tixier and J. Tronczyński	Chemosphere	2008	73 (6)	945-953
19	PYE	Analysis of chlorothalonil and degradation products in soil and water by GC/MS and LC/MS	Alicia Chaves, Damian Shea, David Danehower	Chemosphere	2008	71(4)	629-638
20	PYE	Polychlorinated Naphthalenes and Other Dioxin-Like Compounds in Elbe River Sediments	W Brack, L Bláha, J.P Giesy, M Grote, M Moeder, S Schrader and M Hecker	Environmental Toxicology and Chemistry	2008	27 (3)	519-528

No.	Column	TITLE	AUTHOR	JOURNAL	YEAR	VOL. (ISSUE)	PAGE
21	PYE	Extraction and high-performance liquid chromatographic analysis of C <sub>60</sub> , C <sub>70</sub> , and [6,6]-phenyl C <sub>61</sub> -butyric acid methyl ester in synthetic and natural waters	Dermont Bouchard, and Xin Ma	Journal of Chromatography A	2008	1203 (2)	153-159
22	PYE	Orthogonal screening system of columns for supercritical fluid chromatography	C. West, and E. Lesellier	Journal of Chromatography A	2008	1203 (1)	105-113
23	PYE	A unified classification of stationary phases for packed column supercritical fluid chromatography	C. West, E. Lesellier	Journal of Chromatography A	2008	1191(1-2)	21-39
24	PYE	Automated fractionation procedure for polycyclic aromatic compounds in sediment extracts on three coupled normal-phase high-performance liquid chromatography columns	Urte Lübcke-von Varel, Georg Streck, Werner Brack	Journal of Chromatography A	2008	1185 (1)	31-42
25	PYE	Characterisation of stationary phases in supercritical fluid chromatography with the solvation parameter model: V. Elaboration of a reduced set of test solutes for rapid evaluation	C. West, E. Lesellier	Journal of Chromatography A	2007	1169 (1-2)	205-219
26	PYE	Fractionation of chlorinated and brominated persistent organic pollutants in several food samples by pyrenyl-silica liquid chromatography prior to GC-MS determination	Belén Gómara, Carmen García-Ruiz, María José González, María Luisa Marina	Analytica Chimica Acta	2006	565 (2)	208-213
27	PYE	Rapid Method for Determination of Dioxin-Like Polychlorinated Biphenyls and Other Congeners in Marine Sediments Using Sonic Extraction and Photodiode Array Detection	J. Buzitis, G. M. Ylitalo and M. M. Krahn	Archives of Environmental Contamination and Toxicology	2006	51 (3)	337-346
28	PYE	Occurrence of Polybrominated Biphenyls, Polybrominated Dibenzo-p-dioxins, and Polybrominated Dibenzofurans as Impurities in Commercial Polybrominated Diphenyl Ether Mixtures	Nobuyasu Hanari, Kurunthachalam Kannan, Yuichi Miyake, Tsuyoshi Okazawa, Prasad Rao S. Kodavanti, Kenneth M. Aldous, and Nobuyoshi Yamashita	Environ. Sci. Technol.	2006	40 (14)	4400-4405
29	NPE 5C <sub>18</sub> -MS	Synthesis of cyclic bis(3'-5')diguanilyc acid (c-di-GMP) analogs	Mamoru Hyodo, Yumi Sato and Yoshihiro Hayakawa	Tetrahedron	2006	62(13)	3089-3094
30	HILIC	Development and validation of a reversed-phase high-performance liquid chromatographic method for quantification of peptide dendrimers in human skin permeation experiments	S. Mutalik, A.K. Hewavitharana, P.N. Shaw, Y.G. Anissimov, M.S. Roberts and H.S. Parekh	Journal of Chromatography B	2009	877(29)	3556-3562
31	HILIC	Direct Evidence for Efficient Transport and Minimal Metabolism of L-Cephalexin by Oligopeptide Transporter 1 in Budded Baculovirus Fraction	Keisuke Mitsuoka, Ikumi Tamai, Yasushi Morohashi, Yoshiyuki Kubo, Ryoichi Saitoh, Akira Tsuji and Yukio Kato	Biological & Pharmaceutical Bulletin	2009	32(8)	1459-1461
32	HILIC 5C <sub>18</sub> -MS-II	Structures of Acetylated Oleanane-type Triterpene Saponins, Rarasaponins IV, V, and VI, and Anti-hyperlipidemic Constituents from the Pericarps of <i>Sapindus rarak</i>	Yasunobu ASAO, Toshio MORIKAWA, Yuanyuan XIE, Masaki OKAMOTO, Makoto HAMAOKA, Hisashi MATSUDA, Osamu MURAOOKA, Dan YUAN, and Masayuki YOSHIKAWA	Chem. Pharm. Bull.	2009	57(2)	198-203
33	HILIC	Simultaneous measurement of diazolidinyl urea, urea, and allantoin in cosmetic samples by hydrophilic interaction chromatography	Takahiro Doi, Keiji Kajimura, Satoshi Takatori, Naoki Fukui, Shuzo Taguchi and Shozo Iwagami	Journal of Chromatography B	2009	877(10)	1005-1010
34	HILIC	Convergent synthesis of oligomers of triazole-linked DNA analogue (TLDNA) in solution phase	Tomoko Fujino, Naomi Yamazaki and Hiroyuki Isobe	Tetrahedron Letters	2009	50(28)	4101-4103
35	HILIC	Determination of para-aminohippuric acid (PAH) in human plasma and urine by liquid chromatography-tandem mass spectrometry	Phey Yen Han, P. Nicholas Shaw and Carl M.J. Kirkpatrick	Journal of Chromatography B	2009	877(27)	3215-3220
36	HILIC 5C <sub>18</sub> -MS-II	Medicinal Flowers. XXI. Structures of Perennisaponins A, B, C, D, E, and F, Acylated Oleanane-Type Triterpene Oligoglycosides, from the Flowers of <i>Bellis perennis</i>	Masayuki Yoshikawa, Xuezheng Li, Eriko Nishida, Seikou Nakamura, Hisashi Matsuda, Osamu Muraoka and Toshio Morikawa	CHEMICAL & PHARMACEUTICAL BULLETIN	2008	56(No.4)	559-568
37	HILIC	Separation efficiencies in hydrophilic interaction chromatography	Tohru Ikegami, Kouki Tomomatsu, Hirotaka Takubo, Kanta Horie and Nobuo Tanaka	Journal of Chromatography A	2008	1184(1-2)	474-503
38	HILIC	Triazole-Linked Analogue of Deoxyribonucleic Acid (TLDNA): Design, Synthesis, and Double-Strand Formation with Natural DNA	Hiroyuki Isobe, Tomoko Fujino, Naomi Yamazaki, Marine Guillot-Nieckowski and Eiichi Nakamura	ORGANIC LETTERS	2008	10(17)	3729-3732
39	HILIC	Tetrodotoxin poisoning evidenced by solid-phase extraction combining with liquid chromatography-tandem mass spectrometry	Hsiao-Chin Jen, Shin-Jung Lin, Yung-Hsiang Tsai, Chun-Hsiang Chen, Zu-Chun Lin and Deng-Fwu Hwang	Journal of Chromatography B	2008	871(1)	95-100
40	CNT-300 CNT-1000 CNT-2000	Chromatographic Length-Separation and Photoluminescence Study on DNA-Wrapped Single-Wall and Double-Wall Carbon Nanotubes	Yuki Asada, Toshiaki Sugai, Ryo Kitaura and Hisanori Shinohara	Journal of nanomaterials		in press	



# IV TECHNICAL NOTE

1. Selectivity of packing materials in reversed phase liquid chromatography ...	154
2. Preparation of mobile phase for HPLC .....	161
3. Sample pretreatment for HPLC .....	164
4. Baseline noise in gradient elution .....	168
5. Troubleshooting for increased pressure .....	169
6. Effect of guard column .....	172
7. Troubleshooting for normal phase chromatography .....	173
8. Inner diameter of column (scale down and scale up) .....	175












## 1. Selectivity of packing materials in reversed phase liquid chromatography

Reversed phase chromatography is the most commonly used method of HPLC, because of the high theoretical plate number, excellent separation characteristics, reproducibility, and ease of use. Columns packed with octadecyl group bonded type silica gel (C<sub>18</sub>, ODS) are the most widely used reversed phase chromatography. However, C<sub>18</sub> columns provide insufficient separation for compounds similar in hydrophobicity because the main separation mechanism of C<sub>18</sub> column is based on hydrophobic interaction. It may improve separation of compounds with similar hydrophobicity by using longer columns, changing mobile phases or changing temperature. However, in many cases, it is probably most effective to use different packing materials which retain compounds based on a secondary interaction in addition to hydrophobic interaction.

At Nacalai Tesque, we offer a variety of COSMOSIL reversed phase packing materials. Summary of these packing materials and their respective retention mechanism are in Table 1. Retention of compounds in each stationary phase depends on summation of the interactions. Therefore, comprehension of each interaction leads to selection of an appropriate column.


Table 1. Stationary phase and interaction of packing materials

	C <sub>18</sub> -MS-II	C <sub>18</sub> -AR-II	C <sub>8</sub> -MS	PE-MS	πNAP	PYE	NPE	PBB-R	Cholester
Silica gel	High purity porous spherical silica								
Average particle size	3 · 5 · 15 μm			5 μm					
Average pore size	approx. 120 Å								
Specific surface area	approx. 300 m <sup>2</sup> /g								
Stationary phase									
Types	Monomeric	Polymeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric	Monomeric
Interaction	Hydrophobic interaction	Hydrophobic interaction	Hydrophobic interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction π-π interaction	Hydrophobic interaction Dispersion interaction	Hydrophobic interaction Shape selectivity
End capping	Near-perfect treatment								
Carbon content	approx. 16%	approx. 17%	approx. 10%	approx. 10%	approx. 11%	approx. 18%	approx. 9%	approx. 8%	approx. 20%

## 1) Selectivity for polar functional group

### Selectivity

Selectivity for polar functional group is evaluated based on the separation of benzene, nitrobenzene, which has a nitro group, and anisole, which has a methoxy group. The chromatograms below show separation of the three compounds on four COSMOSIL columns: C<sub>18</sub>-MS-II, PE-MS, πNAP and PYE. Elution order on the C<sub>18</sub> column is as following: nitrobenzene, anisole and benzene. Elution orders on the aromatic columns are reversed. Separation on the C<sub>18</sub> column is based on hydrophobic interaction only. On the other hand, the packing materials on the other three columns have aromatic rings and reverse the elution order by π-π interaction.





The graph of selectivity for polar functional group is shown below. Among nine COSMOSIL columns, PYE and NPE columns have the highest selectivity factors for polar groups. As to mobile phases, methanol is more effective than acetonitrile for separation using π-π interaction.

### Selectivity for polar functional group

Column: 4.6mm I.D.-150mm  
 Mobile phase: Methanol / H<sub>2</sub>O = 50/50  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1; Nitrobenzene (0.13 μg)  
 2; Anisole (1.5 μg)  
 3; Benzene (4.0 μg)



### Application data



#### • Separation of toluenitrile position isomers

Toluenitriles have three position isomers. It is difficult to separate ortho and para isomers by C<sub>18</sub> or PE-MS column because of lack of poor π-π interaction. On the other hand, the isomers are well separated on PYE or NPE column which has strong π-π interaction.

### COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Mobile phase: 5C<sub>18</sub>-MS-II, 5PE-MS, 5NPE: Methanol/ H<sub>2</sub>O = 40/60  
 5PYE: Methanol/ H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: o-Toluenitrile (2.0 μg)  
 m-Toluenitrile (2.0 μg)  
 p-Toluenitrile (1.0 μg)



## 2) Selectivity for dipole

## Selectivity


Selectivity for dipole is evaluated based on the separation of 1,5-dinitronaphthalene and 1,8-dinitronaphthalene.

Dinitronaphthalenes (peak 1 and 2) were strongly retained on PYE and NPE because of  $\pi$ - $\pi$  interaction compared with dimethylnaphthalenes. However, there is a slight difference between these two columns. While 1,5-dinitronaphthalene (peak 2) was preferentially retained on PYE,

1,8-dinitronaphthalene (peak 1) was retained longer on NPE.


The results with NPE indicate the presence of strong dipole-dipole interaction. The two nitro group dipoles in

1,8-dinitronaphthalene are aligned for a much greater dipolar coupling with the bonded nitrophenyl group in NPE than 1,5-dinitronaphthalene.



## Selectivity for dipole

Column size: 4.6mm I.D.-150mm  
 Mobile phase: C<sub>18</sub>-MS-II Methanol / H<sub>2</sub>O = 80/20  
 NPE Methanol / H<sub>2</sub>O = 70/30  
 PYE Methanol / H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample:  
 1; 1,8-Dinitronaphthalene (1,8-DNN) (0.21  $\mu$ g)  
 2; 1,5-Dinitronaphthalene (1,5-DNN) (0.11  $\mu$ g)  
 3; Naphthalene (0.25  $\mu$ g)  
 4; 1-Methylnaphthalene (0.35  $\mu$ g)  
 5; 1,5-Dimethylnaphthalene (0.42  $\mu$ g)



## Application data


## ● Separation of phthalonitrile position isomers

Phthalonitriles have three position isomers. NPE or PYE completely separates these compounds due to  $\pi$ - $\pi$  interaction.

Furthermore, NPE strongly retains o-phthalonitrile due to dipole-dipole interaction.

## COSMOSIL Application Data

Column size: 4.6mm I.D.-150mm  
 Mobile phase: SC<sub>18</sub>-MS-II Methanol / H<sub>2</sub>O = 30/70  
 5NPE Methanol / H<sub>2</sub>O = 40/60  
 5PYE Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm  
 Sample:  
 o-; Phthalonitrile (0.3  $\mu$ g)  
 m-; Isophthalonitrile (3.0  $\mu$ g)  
 p-; Terephthalonitrile (0.15  $\mu$ g)




## 3) Selectivity for polyaromatic compounds

## Selectivity

Selectivity for polyaromatic compounds is evaluated based on the separation of benzene, naphthalene and anthracene.



The elution orders in all columns are the same : benzene, naphthalene and anthracene. Retention increases in all columns with increasing number of aromatic rings. In addition, highly dispersive packing materials such as PBB and PYE show much stronger retention for polyaromatic compounds due to dispersion interaction.



## Selectivity for polyaromatic compounds

Column size: 4.6mm I.D.-150mm  
 Mobile phase: SC<sub>18</sub>-MS-II Methanol/H<sub>2</sub>O = 90/10  
 5PBB-R Methanol/H<sub>2</sub>O = 90/10  
 5PYE Methanol/H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample:  
 1; Benzene (1.67  $\mu$ g)  
 2; Naphthalene (0.11  $\mu$ g)  
 3; Anthracene (0.0063  $\mu$ g)



## Application data



## ● Separation of dibenzosuberone and dibenzosuberone

C<sub>18</sub> retains dibenzosuberone (peak 1) longer than dibenzosuberone (peak 2). On the other hand, PBB-R and PYE retain dibenzosuberone (peak 2), which has a  $\pi$ -electron conjugated system, longer than dibenzosuberone (peak 1).

## COSMOSIL Application Data

Column size: 4.6mm I.D.-150mm  
 Mobile phase: SC<sub>18</sub>-MS-II Methanol/H<sub>2</sub>O = 80/20  
 5PBB-R Methanol  
 5PYE Methanol/H<sub>2</sub>O = 90/10  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample:  
 1; Dibenzosuberone (0.1  $\mu$ g)  
 2; Dibenzosuberone (0.025  $\mu$ g)



## 4) Selectivity for molecular shape

## Selectivity



Selectivity for molecular shape is evaluated based on the separation of chrysene and benz [a] anthracene. The isomers of two polyaromatic hydrocarbons, which consist of four benzene rings, are difficult to separate because of the similar hydrophobicity or aromaticity. However, PYE and Cholester columns, which recognize molecular shape, enable them to separate chrysene and benz [a] anthracene.



## Selectivity for molecular shape

Column: 4.6mm I.D.-150mm  
 Column size: 5C<sub>18</sub>-MS-II, 5PYE  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Cholester  
 Methanol  
 5PE-MS  
 Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Benz[a]anthracene (0.04μg)  
 2: Chrysene (0.04μg)



## Application data



## • Separation of diastereomers (threo- and erythro-)

C<sub>18</sub> cannot separate the threo and erythro forms. On the other hand, PYE retains the planar erythro form longer than the threo form.

## COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Column size: 5C<sub>18</sub>-MS-II, 5PYE  
 Mobile phase: Methanol / H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm


Sample: 1: Threo form  
 2: Erythro form



## 5) Selectivity for halides

## Selectivity



Selectivity for halide is evaluated based on the separation of chlorobenzene and bromobenzene. PBB-R shows the highest selectivity factor due to dispersion interaction of the five bromine atoms.



## Selectivity for halide

Column: 4.6mm I.D.-150mm  
 Column size: 5C<sub>18</sub>-MS-II, 5PBB-R, 5PYE  
 Mobile phase: Methanol / H<sub>2</sub>O = 60/40  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Toluene (3.3μg)  
 2: Chlorobenzene (3.3μg)  
 3: Bromobenzene (3.3μg)



## Application data



## • Separation of halogen exchange reaction products

PYE and PBB-R retain dispersed iodine atom longer than bromine atom. As a result, PYE and PBB-R can separate the complicated bromine and iodine compounds that C<sub>18</sub> cannot separate.

## COSMOSIL Application Data

Column: 4.6mm I.D.-150mm  
 Column size: 5C<sub>18</sub>-MS-II, 5PBB-R, 5PYE  
 Mobile phase: Methanol / H<sub>2</sub>O = 90/10  
 Methanol  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: R=I or Br




Sample courtesy of Dr. H. Yamamoto, RIKEN, Condensed Molecular Materials Laboratory

## 6) Selectivity for hydrophobicity

### Selectivity


Selectivity for hydrophobicity is evaluated based on the separation of alkylbenzenes. Two C<sub>18</sub> and Cholester show similar high selectivity for hydrophobicity. Other columns show less hydrophobic selectivity than C<sub>18</sub>.



### Selectivity for hydrophobicity

Column: 5C<sub>18</sub>-MS-II  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase: Methanol/H<sub>2</sub>O = 80/20  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Benzene (1.67µg)  
 2: Toluene (1.67µg)  
 3: Ethylbenzene (1.67µg)  
 4: Propylbenzene (1.67µg)  
 5: Butylbenzene (1.67µg)  
 6: Amylbenzene (1.67µg)




Lower concentration of organic solvent in mobile phase leads to much retention in reversed phase chromatography. In case of NPE, when methanol concentration is reduced to 60%, the retention times increase to those similar to C<sub>18</sub> with 80% methanol.

### Control of retention time

Column: 5NPE  
 Column size: 4.6mm I.D.-150mm  
 Mobile phase:  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C  
 Detection: UV254nm

Sample: 1: Benzene (1.67µg)  
 2: Toluene (1.67µg)  
 3: Ethylbenzene (1.67µg)  
 4: Propylbenzene (1.67µg)  
 5: Butylbenzene (1.67µg)  
 6: Amylbenzene (1.67µg)



## 2. Preparation of mobile phase for HPLC

### 1) Organic solvent / aqueous mixed mobile phase

1)-1. Preparation of methanol : water = 70 : 30 (v/v) 1L

- ① Measure 700 ml of methanol in a measuring cylinder.
- ② Measure 300 ml of distilled water in a measuring cylinder.
- ③ Mix ① and ② thoroughly and degas.


The better approach is to prepare the mobile phase gravimetrically rather than volumetrically. Following is example of preparation.

Composition table for mobile phase 1L using methanol and water			Composition table for mobile phase 1L using acetonitrile and water		
Methanol / Water	Methanol (g)	Distilled water (g)	Acetonitrile / Water	Acetonitrile (g)	Distilled water (g)
90/10 (v/v)	711.9	99.8	90/10 (v/v)	707.4	99.8
80/20 (v/v)	632.8	199.6	80/20 (v/v)	628.8	199.6
70/30 (v/v)	553.7	299.5	70/30 (v/v)	550.2	299.5
60/40 (v/v)	474.6	399.3	60/40 (v/v)	471.6	399.3
50/50 (v/v)	395.5	499.1	50/50 (v/v)	393.0	499.1
40/60 (v/v)	316.4	598.9	40/60 (v/v)	314.4	598.9
30/70 (v/v)	237.3	698.7	30/70 (v/v)	235.8	698.7
20/80 (v/v)	158.2	798.6	20/80 (v/v)	157.2	798.6
10/90 (v/v)	79.1	898.4	10/90 (v/v)	78.6	898.4

Caution : Methanol and acetonitrile are hazardous substances, do not use for medical purpose. Always process in a laboratory hood and wear an eye protection and a mask.

Influence of organic solvent composition in mobile phase on the retention time.


1% difference in the composition significantly changes the retention.



Special attention should be paid to measure correct amount of organic solvent as retention time is significantly changed by 1% different composition.

Column: 5C<sub>18</sub>-MS-446mm I.D. x 150mm  
 Flow rate: 1.0ml/min  
 Detection: 254nm 0.16AUFS  
 Temperature: 30°C  
 Sample: 1. Uracil  
 2. Naphthalene


### Differences between acetonitrile and methanol in reversed phase liquid chromatography



Column: COSMOSIL 5C<sub>18</sub>-MS-II  
 Column size: 4.6 mm I.D. x 150 mm  
 Flow rate: 1.0 ml/min  
 Temperature: 30°C

Back pressure differs depending on species of organic solvents and mixing ratio. Back pressure of acetonitrile is lower than that of methanol at the same concentration.





## 2) Organic solvent/ aqueous mixed mobile phase

### 2)-1. Preparation of methanol : 20 mmol/l phosphate buffer (pH2.5) = 80 : 20 (v/v) 1L

- Preparation of 20 mmol/l phosphate buffer (pH2.5)
  - Prepare 20 mmol/l sodium dihydrogenphosphate aqueous solution.
  - Prepare 20 mmol/l phosphoric acid aqueous solution.
  - Adjust the pH to 2.5 by mixing ① with ②.
  - Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.31 g of sodium dihydrogenphosphate and 1.05 g of phosphoric acid in distilled water to make 1 L solution.
- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).
- Confirm that the solution is pH2.5.
  - Preparation of methanol : 20 mmol/l phosphate buffer (pH2.5) = 80 : 20 1L
    - Measure 800 ml of methanol in a measuring cylinder.
    - Measure 200 ml of 20 mmol/l phosphate buffer (pH2.5) in a measuring cylinder.
    - Mix ① and ② thoroughly and degas.

### 2)-2. Preparation of methanol : 20 mmol/l phosphate buffer (pH7.0) = 80 : 20 (v/v) 1L

- Preparation of 20 mmol/l phosphate buffer (pH7.0)
  - Prepare 20 mmol/l sodium dihydrogenphosphate aqueous solution.
  - Prepare 20 mmol/l di-sodium hydrogenphosphate aqueous solution.
  - Adjust the pH7.0 by mixing ① with ②.
  - Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.14 g of sodium dihydrogenphosphate and 1.49 g of di-sodium hydrogenphosphate in distilled water to make 1L solution.
- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).
- Confirm that the solution is pH7.0.

- Preparation of Methanol : 20 mmol/l phosphate buffer (pH7.0) = 80 : 20 1 L

- Measure 800 ml of methanol in a measuring cylinder.
- Measure 200 ml of 20 mmol/l phosphate buffer (pH7.0) in a measuring cylinder.
- Mix ① and ② thoroughly and degas.

The better approach is to prepare the mobile phase gravimetrically rather than volumetrically. Following is example of preparation.

Methanol : 20 mmol/l Phosphate buffer	Methanol (g)	20 mmol/l Phosphate buffer (pH2.5) (g)	20 mmol/l Phosphate buffer (pH7.0) (g)
90 / 10 (v/v)	711.9	99.8	99.9
80 / 20 (v/v)	632.8	199.6	199.8
70 / 30 (v/v)	553.7	299.4	299.7
60 / 40 (v/v)	474.6	399.2	399.6
50 / 50 (v/v)	395.5	499.0	499.5
40 / 60 (v/v)	316.4	598.8	599.4
30 / 70 (v/v)	237.3	698.6	699.3
20 / 80 (v/v)	158.2	798.4	799.2
10 / 90 (v/v)	79.1	898.2	899.1

Caution : Methanol and acetonitrile are hazardous substances, do not use for medical purpose. Always process in a laboratory hood and wear an eye protection and a mask.

## 3) Preparation of ion pair reagent containing mobile phase

### 3)-1. Preparation of 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l phosphate buffer (pH2.5)

- Prepare 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l sodium dihydrogenphosphate aqueous solution.
- Prepare 5 mmol/l sodium 1-butanefulfonate containing 20 mmol/l phosphoric acid aqueous solution
- Adjust the pH to 2.5 by mixing ① with ②.
- Filter ③ under reduced pressure to remove insoluble substance that may deteriorate pump-seal and clog columns (0.45 μm or smaller pore size is recommended).

#### (Easy method)

- Dissolve 1.31 g of sodium dihydrogenphosphate, 1.05 g of phosphoric acid and 0.80 g of sodium 1-butanefulfonate in distilled water to make 1L solution.
- Filter the solution under reduced pressure to remove insoluble substance (0.45 μm or smaller pore size is recommended).
- Confirm that the solution is pH2.5.



※ 0.5 M Sodium 1-butanefulfonate aqueous solution is also available from Nacalai Tesque.

### 3. Sample pretreatment for HPLC

Pretreatment before HPLC analysis is often required for samples of low concentration or samples containing analytical contaminants. It improves reproducibility and sensitivity in analysis, and protects HPLC columns. The preparation methods are different according to the each sample. The followings are examples of different pretreatments.

#### 1) Filtration


Filtration is a common method used for separating solids from liquids. It extends a column's life by minimizing column damages from solid contaminants such as particles, sediments and colloid substances. It also improves reproducibility of analytical data. We offer both syringe-type and spin-type filters for sample filtration.

	■ Syringe filter	■ Centrifugal filter
Product	Cosmonice filter	Cosmospin filter
Configuration		
Usage	Easy to use. Just attach a filter on top of a syringe.	Easy to use by centrifugation.
Type	<ul style="list-style-type: none"> <li>• W (aqueous system)</li> <li>• S (solvent system)</li> </ul>	<ul style="list-style-type: none"> <li>• Pore diameter : 0.2 µm</li> <li>• Pore diameter : 0.45 µm</li> </ul>
Required equipment	Syringe • Sample bottle	Centrifuge
Page	refer to page 74	refer to page 74

#### Cosmonice filter

How to use :

- ① Fill a syringe with the sample you want to filter.
- ② Attach a Cosmonice filter to the syringe.
- ③ Push the syringe plunger to filter the sample.
- ④ Analyze the filtered sample by HPLC.




#### Cosmospin filter

Components : • Sample reservoir  
• Sample collection tube


How to use :

- ① Insert a Cosmospin sample reservoir into a Cosmospin sample collection tube.
- ② Add a sample into the Cosmospin sample reservoir.
- ③ Close the sample collection tube cap and centrifuge.
- ④ Remove the sample reservoir and collect the filtered sample in the sample collection tube.
- ⑤ Analyze the filtered sample by HPLC.



#### 2) Protein precipitation

Protein precipitation is commonly used to remove proteins in samples for downstream analysis. For example, when analyzing drug concentration in blood samples, proteins have to be removed first. Otherwise, proteins may be adsorbed in columns and interfere with the analysis. Common methods for protein precipitation include salting out, isoelectric point precipitation and precipitation with organic solvents. The following shows a general procedure for protein precipitation with organic solvents.




Procedure for protein precipitation :

#### 3) Ultrafiltration


Ultrafiltration is a method to concentrate proteins or other macromolecules through a semipermeable membrane with defined pores. Ultrafiltration is applicable for sample desalting, concentrating proteins from dilute solution such as urine samples, or deproteinizing samples with high protein concentration (i.e. blood serum or plasma). Following is a general procedure for ultrafiltration.

Procedure for ultrafiltration :




#### 4) Solvent extraction method

Solvent extraction is a method to separate compounds due to their unequal solubility in two immiscible liquid phases, usually water and an organic solvent. The method is used to concentrate highly hydrophobic compounds, and consequently increase analytical sensitivity. A buffer solution is added to sample to optimize the pH and target substance is then extracted by an organic solvent such as ether and chloroform. However, when target substance is combined with proteins, solvent extraction may not work well.




#### 5) Ion exchange

Pretreatment by ion-exchange resin may be effective for samples that the solvent extraction method cannot be adapted due to its emulsification. A preliminary experiment may be required for the selection of resin and experimental conditions. For example, a negatively charged compound is strongly adsorbed on an anion-exchange resin such as DEAE cellulose resin. Therefore, the target compound is collected by increasing salt concentration of buffer solution or adjusting pH of elution buffer after washing off other weakly adsorbed undesired substances.



## 4. Baseline noise in gradient elution

In gradient analysis, incomplete mixing of mobile phases or impurities in water of mobile phase can cause baseline noise. In the former case, it can be improved by using a proper mixer before injector (Baseline 1→2). In the latter case, it can be improved by using a pre-column. Impurities in water are adsorbed on the pre-column (Baseline 2→3). COSMOSIL 5C<sub>18</sub>-AR-II 4.6 mm I.D. x 10 mm or 10 mm I.D. x 20 mm as a pre-column.




Column COSMOSIL 5C<sub>18</sub>-AR-300 4.6 mm I.D. x 150mm  
 Pre-column COSMOSIL 5C<sub>18</sub>-AR-II 4.6 mm I.D. x 10mm  
 Mobile phase A: 0.1% TFA containing water  
 B: 0.1% TFA containing 95% acetonitrile  
 B: 0% → 100%/30 min liner gradient  
 Flow rate 1.0 ml/min  
 Temperature 30°C  
 Detection UV 220nm

## 5. Troubleshooting for increased pressure

Repeated analysis may increase back pressure. Continuous use of HPLC columns under high pressure can cause deterioration and overload of the equipment. Therefore, it is important to monitor column back pressure regularly and solve the problem timely.

The back pressure increase can be due to clogging of a column or clogging of the equipment. First of all, identify the clogging site.

Schematic diagram of HPLC



Remove analytical column first and connect the plumbing from the back of HPLC system directly to the detector. Measure the pressure of flowing mobile phase without an HPLC column. Generally equipment should hardly generate any pressure. If there is significant flow pressure, disconnect the system components one by one to identify the clogged component(s). Possible causes and solutions of clogged equipment are discussed in section II below.

If the flow pressure without a column is normal, then pressure increase is due to clogging of a column. In this case, one needs to determine the causes and whether it is time to replace the column. Possible causes and solutions of clogged column are discussed in greater details in sections I and III.

Symptom	Possible Cause
Pressure increase rapidly in short-term use	→ Flow pressure without a column is 0-0.3 MPa → Clogging of column Refer to section I
	Yes Flow pressure with a column is 0-0.3 MPa or higher → Clogging of equipment Refer to section II
Pressure increase gradually in long-term use	→ Deterioration of column due to long-term use refer to section III

### I. Solution in case an HPLC column is clogged in short-term use.

Select the possible cause of clogging according to the following flow chart.

STEP 1	• Salt deposition • Use mobile phase of high concentration organic solvent right after using buffer	→ YES	Cause 1
↓ NO			
STEP 2	• Forget to filter mobile phase • Sample is not dissolved enough	→ YES	Cause 2
↓ NO			
STEP 3	• Analyzing samples which tend to adsorb to a column (i.e. protein samples) • Sample deposition in column	→ YES	Cause 3



**Cause 1** Salt is deposited on a column.  
**Solution :** Wash columns for 30 minutes at half of the analytical flow rate with 10% organic solvent (methanol or acetonitrile) in water to dissolve deposited salt. If the situation is not improved, wash with 100% water under the same condition.  
**Prevention :** To switch to high concentration organic solvent after using a buffer, first wash a column with a mobile phase not containing salt (with the same concentration of organic solvent as the buffer), then switch to the mobile phase of higher organic solvent concentration. Example : To change mobile phase from 10/90 (v/v) acetonitrile/20mmol/l phosphate buffer (pH2.5) to 90/10 (v/v) acetonitrile/water, first wash for 15 minutes with 10/90 (v/v) acetonitrile/water, and then switch to 90/10 (v/v) acetonitrile/water.

**Cause 2** Column filter is clogged by sample or impurities.  
**Solution :** Connect the column in reverse direction, and then wash the column for 30 minutes at half of the usual analytical flow rate with the mobile phase used for analysis. If the situation is not improved, change the end fitting in the front of column. (We can replace end fittings with a paid service fee.)  
**Prevention :** We recommend filtering sample and/or mobile phase. For more information, please refer to page 164, TECHNICAL NOTE 3. Sample pretreatment for HPLC 1) filtration.

**Cause 3** Sample may be adsorbed to packing material or deposited in a column.  
**Solution :** Wash for 30 minutes at half of analytical flow rates with a solvent which adsorbed substances are dissolved in. The followings are how to wash each type of columns.  
 [Reversed phase columns]  
 a) When an adsorbed substance is not protein, wash with methanol or tetrahydrofuran.  
 b) When an adsorbed substance is protein, wash with 50-70% of acetonitrile/water (containing 0.1% of trifluoroacetic acid). However proteins may be deposited in high concentration of organic solvent depending on varieties.  
 [COSMOSIL Sugar-D/NH<sub>2</sub>/HILIC columns]  
 Wash with 50/50 (v/v) acetonitrile/water for NH<sub>2</sub>-MS and 100% water for Sugar-D and HILIC columns.  
 [COSMOSIL SL-II]  
 Wash with methanol, tetrahydrofuran or ethanol.

**Prevention :** Choose appropriate pretreatment for each sample. For more information, please refer to page 164, TECHNICAL NOTE 3. Sample pretreatment for HPLC 1) filtration. We also recommend using guard column. For more information for guard columns, please refer to page 172, TECHNICAL NOTE 6. Effect of guard columns.

**Caution :**

- When wash columns, do not connect column exit and let the solution through.
- Long term of washing may deteriorate the performance of columns.
- Do not use strongly alkaline solution (more than pH 7.5) or strongly acidic solution (less than pH 1.5) for silica gel base packing material.
- Store columns with manufacturer recommended storage solvent after washing, When the situation is not improved, replace the column.

## II. Solutions in case pressure is too high because of clogged equipment.

First, identify the specific clogging site by disconnecting the components in the system one by one and checking the flow pressure. The followings are possible common causes.

**Cause 1** Salt is deposited in plumbing.  
**Solution :** Flow water to the plumbing without connecting a column and any other equipment. Washing out the plumbing in a reversing connection is also an effective way. If the situation is not improved, replace it with a new one.

**Cause 2** Check-valve of pump is clogged by stain  
**Solution :** Wash the check-valve with a stain dissolving solvent. Take apart the washable part, soak it in the solvent, then clean in an ultrasonic cleaner.

**Cause 3** Manual injector is clogged with stain  
**Solution :** Wash with a stain dissolving solvent. Soak rotor seal and line filter in water and clean them in an ultrasonic cleaner. If the situation is not improved, replace the injector.  
**Prevention :** It extends the life time of an HPLC system to maintain regular wash of the system. Wash the system the same as wash an HPLC column. When the mobile phase contains salt, wash for 10-15 minutes with a mobile phase which has the same composition but not containing salt. For example, when using 50/50 (v/v) methanol/20mmol/l phosphate buffer, wash with 50/50 (v/v) methanol/water. When the mobile phase contains halogen, acid and/or base, wash for 10-15 minutes with mobile phase which has the same composition but not containing halogen, acid and/or base.

## III. Solutions in case a column is damaged from long term use

Every column will have to be replaced eventually. Performance of a column is expected to deteriorate slowly after long term use. One has to decide whether it is time to replace the column.

**Cause 1** Column deterioration result from long term use.  
**Solution :** Wash according to Section I, Cause 3.  
**Prevention :** Same as prevention in Section I, Cause 3. When the column condition is not improved, you could continue to use the column if peak shapes do not change and the maximum pressure is less than 20 MPa. However, we recommend replacing the column because it place extra burden on the equipment.

**Cause 2** Silica gel in the column may be cracked because of long term use.  
**Solution :** Replace the columns.




## 6. Effect of guard column

The use of guard columns to protect both analytical and preparative columns is highly recommended. COSMOSIL guard columns are packed with packing materials identical to that used in analytical and preparative columns. As a result, the COSMOSIL guard columns do not contribute to any decrease in the performance of the main column.


### Example of using guard columns

The following chromatograms show analysis examples using a COSMOSIL 5C<sub>18</sub>-MS-II analytical column (4.6 mm I.D. × 150 mm) and the same column connected with its guard column (4.6 mm I.D. × 10 mm). There is no change in separation characteristics since the packing material of the guard column is identical to that used in the main column.



### Inspection of columns

We recommend that the performance of a column be examined regularly. Deteriorated columns shall be replaced timely. If a deteriorated guard column continues to be used, the packed column will also deteriorate.



### Ordering information

Please refer to respective pages of each column.

Please refer to page 76 for the parts necessary to use guard columns or guard column cartridges.

## 7. Troubleshooting for normal phase chromatography

**Q 1** : How can I convert from reversed phase mode to normal phase mode or vice versus using the same HPLC equipment?

**A 1** : To convert from reversed phase mode to normal phase mode, or vice versus, flush the equipment with a solvent that is miscible with both the current mobile phase and the intended mobile phase. Connect the HPLC pump directly with the detector, and replace the solvents according to following instructions.

Solvent conversion from reversed phase to normal phase

① To convert a mobile phase without buffer solution in reversed phase to normal phase, replace solvents according to the following steps :

- 1) Flush the equipment with a solvent for reversed phase. For example, methanol/H<sub>2</sub>O (v/v=50/50).
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example, tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent for normal phase.

② To convert a mobile phase with buffer solution in reversed phase to normal phase, replace solvents according to the following steps :

- 1) Flush the equipment with a solvent with buffer solution for reversed phase.
- 2) Flush the equipment with a solvent with the composition same as 1) and without salt.
- 3) Flush the equipment with a solvent miscible to both mobile phases. For example, tetrahydrofuran, ethanol.
- 4) Flush the equipment with a solvent for normal phase.

Solvent conversion from normal phase to reversed phase

① To convert from normal phase to a mobile phase without buffer solution in reversed phase, replace solvents according to following steps :

- 1) Flush the equipment with a solvent for normal phase.
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent for reversed phase.

② To convert from normal phase to a mobile phase with buffer solution in reversed phase, replace solvents according to following steps :

- 1) Flush the equipment with a solvent for normal phase. For example, hexane/ethyl acetate.
- 2) Flush the equipment with a solvent miscible to both mobile phases. For example tetrahydrofuran, ethanol.
- 3) Flush the equipment with a solvent with the composition same as 4) and without salt.
- 4) Flush the equipment with a solvent for reversed phase.

**Q 2** : My flow rate is not stable. How can I troubleshoot?

**A 2** : Possible causes for unstable flow rate can be a malfunctioning check valve or air in a mobile phase, wash the check valve thoroughly by ultrasonic cleaner. Solvents with a low boiling point such as n-Hexan and n-Heptane generate air easily. To prevent air generation, degas the mobile phase sufficiently.

**Q 3** : In spite of using the same condition, the retention time is different. How can I solve the problem?

**A 3** : One possible cause is unstable flow rate. Please refer to Q2 section.

Another possible cause is variation of polar component in mobile phases. In normal phase chromatography, the retention time depends on the concentration of small amounts of very polar constituents in the mobile phase. This is especially true for water content in a mobile phase. In this case, always use fresh solvents in a mobile phase. If sample solvent includes water, change to a solvent without water or decrease injection volume. If a column contains water, remove water from the column by washing it with ethanol.

**Q 4** : How can I wash the COSMOSIL SL-II column?

**A 4** : The SL-II column can be washed with tetrahydrofuran, methanol, ethanol, methylene chloride, n-Hexane or n-Heptane.

**Q 5** : How can I store the COSMOSIL SL-II column?

**A 5** : Store the SL- II column with the shipping screw tighten in order to prevent the solvent in the column from volatilization. In case where a solvent containing halogens is used, replace the solvent in the column with a solvent without halogens such as n-Heptane before storing.

**Q 6** : I get peak tailing in my run. What can I do about it?

**A 6** : ① In case where the sample contains acidic compounds, add approx. 0.5% of acetic acid to the mobilephase.  
② In case where the sample contains basic compounds, add approx. 0.5% of triethylamine to the mobile phase.

**Q 7** : I get no peaks. How can I troubleshoot?

**A 7** : First, make sure that there is no problem with the system. If the problems are with the sample, sample solvent or mobile phase, try following check list.  
1. The analyte may not be eluted from a column because the retention of the analyte is too strong. In this case, use a stronger eluent (mobile phase).  
2. The sample contains chelating compounds or basic compounds. They may be adsorbed to the packing materials. In this case, add 0.1% - 1% acid (trifluoroacetic acid or acetic acid) to the mobile phase.

## 8. Inner diameter of column (scale down and scale up)

The figure below shows general parameters for 1.0 mm to 50 mm I.D. COSMOSIL columns : flow rate, equipment, inner diameter of pipe, application, surface ratio (compared with 4.6 mm I.D.) and particle size. It may help to scale up or down from the most commonly used 4.6 mm I.D. column.

Inner diameter (mm I.D.)	1.0	2.0	3.0	4.6	10	20	28	50
Flow rate (ml/min)	0.05	0.2	0.4	1.0	5.0	19	37	70
Detector cell · Injector	for Semi-micro		for Analytical			for Preparative		
Inner diameter of pipe (mm)	0.05	0.1	0.2-0.3			1.0		
Application	LC-MS Solvent saving		Solvent saving with standard system	Standard	Preparative (small scale)	Preparative (medium scale)	Preparative (large scale)	Preparative (super large scale)
Surface ratio with 4.6 mm I.D.	0.05	0.19	0.43	1.00	4.73	18.90	37.05	118.15
Particle size (µm)	3 or 5			5		15 or more		

### Scale down

When scaling down from the most commonly used analytical column (4.6 mm I.D.) to a semi-micro or 3.0 mm I.D. analytical HPLC column (of the same column length), sample loading dose is proportionate to the cross section of column. The 3.0 mm I.D. columns provide high sensitivity and solvent saving without the need to change the existing equipment settings. Semi-micro columns (2.0 mm I.D. and 1.0 mm I.D.) provide higher sensitivity and enable analysis of minor components, but one needs to change the piping of HPLC equipment, the injector and the detector cell for semimicro columns.

Column size	4.6 mm I.D. × 150 mm	3.0 mm I.D. × 150 mm	2.0 mm I.D. × 150 mm	1.0 mm I.D. × 150 mm
Chromatogram				
Flow rate (ml/min)	1.0	0.4	0.2	0.05
Pressure (MPa)	3.4	3.6	3.8	3.6
Injection volume(µl)	1.0	0.4	0.2	0.05
Detector Cell · Injector	for Analytical		for Semi-micro	
Detector sensitivity(AUFS)	0.08		0.04	
Inner diameter of pipe (mm)	0.25		0.10	0.05
	Column Mobile phase Flow rate Temperature Detection	COSMOSIL 5C <sub>18</sub> -MS- II Acetonitrile : Water = 70 : 30 1.0 ml/min 30°C UV 254 nm		Sample 1. Benzene 2. Toluene 3. Ethylbenzene 4. Propylbenzene 5. Butylbenzene 6. Amylbenzene

The 'V' logo is a white, bold, sans-serif letter 'V' centered within a blue rounded square.

# INDEX

Sample name .....	178
Packing material's name .....	188



Sample name			
I. COSMOSIL HPLC column	A	p-Aminophenol	147
		p-Aminophenol Hydrochloride	31, 80
		Acenaphthene	87, 148
		Acenaphthylene	148
		Acesulfame K	123
		Acetaminophen	31, 80, 110
		Acetanilide	90, 99, 110
		Acetic Acid	18
		Acetoaminophen	15
		Acetophenone	14, 30, 147
II. Liquid chromatography related product		4'-Acetoxyacetanilide	31, 80
		17 $\alpha$ -Acetoxyprogesterone	86
		Acetrizic Acid	81
		N-Acetylprocainamide	14, 110
		Acetylsalicylic Acid (Aspirine)	110
		Acid Red	123
		Aconitine	120
		Acrinol	80
		Acrylic Acid	73
		Actinase E	139
III. Application data		Adenine	18, 135
		Adenosine	18, 100, 135
		5'-ADP	135
		L-Adrenaline	18, 72
		[Ala2]-Leu-Enkephalin	138
		[Ala2]-Met-Enkephalin	138
		Albiflorin	115
		Albumin (Bovine)	44, 53, 139
		Albumin (Goat)	140
		Albumin (Human)	139
IV. Technical note		Albumin (Ovalbumin)	140
		Albumin (Rabbit)	140
		Alcohol Dehydrogenase	139
		Alizalin	119
		Allantoic Acid	112
		Allantoin	112
		Allopurinol	129
		Allylbenzene	149
		Amaranth	123
		Amidotrizic Acid	81
V. Index		Amikacin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative	81
		4-Aminoacetophenone	90, 98
		Aminoantipyrine	110
		4-Aminoantipyrine	15
		4-Aminobenzoic Acid	80
		4-Amino-n-butyric Acid [GABA]	137
		4-Amino-6-chlorobenzene-1,3-disulfonamide	98
		4-(Aminomethyl) benzoic Acid	94
		m-Aminophenol	147
		o-Aminophenol	147
	B	Baclofen	97
		Baicalein	23
		Baicalin	113, 115
		Bamethane Sulfate	102
		Barbital	110
		Barbital Sodium	117
		Beclometasone Dipropionate	86, 102, 103
		Bensulide	121
		Benzamide	147
		Benz[a]anthracene	148, 158
	C	Benzene	14, 21, 30, 31, 34, 144, 146, 147, 155, 157, 160
		Benzo[b]fluoranthene	148
		Benzoic Acid	36, 80, 106, 123, 146
		Benzonitrile	147
		Benzo[g,h,i]perylene	148
		Benzophenone	86, 88, 106, 107
		Benzo[a]pyrene	148
		Benz[a]pyrene	148
		Benzyl Alcohol	9, 14, 106
		Benzyl p-Hydroxybenzoate	98
		Benzylpenicillin Potassium	103
		2-Benzylpyridine	148
		3-Benzylpyridine	148
		4-Benzylpyridine	148
		Berberine Chloride	103, 113
		Berberine Hydrochloride	25
		Betahistine Mesilate	103
		Betamethasone	103, 129
		Betamethasone Sodium Phosphate	103
		Betamethasone Valerate	103
		Betamethasone 17-Valerate	172
		Bezafibrate	102
		Bibenzyl	22
		1,1'-Binaphthyl	149
		D-Biotin [Vitamin H]	125
		Biphenyl	97
		Bisoprolol	98
		Boc-(D)Phe-(L)Phe-OMe	137
		Boc-(L)Phe-(L)Phe-OMe	137
		Bombesin	138
		Bovine Serum	142
		BPMC	120
		Bradykinin	138
		Brilliant Blue FCF	123
		Bromacil	122
		Bromhexine Hydrochloride	102
		Bromobenzene	144, 147, 159
		Brucine Dihydrate	116, 117
		Bufalin	115
		Bufexamac	100
		Butylbenzene	12, 21, 34, 160
		n-Butylbenzene	147
		Butylbenzoate	12
		2-tert-Butylhydroquinone	124
		3-tert-Butyl-4-hydroxyanisol	124
		Butyl p-Hydroxybenzoate	16, 82, 84, 85, 86, 88, 98, 101, 103, 108, 117
		iso-Butyl p-Hydroxybenzoate	108
		Butyl Salicylate	16
		C60	54, 56, 58, 59
		C60O	56
		C70	54, 56, 58
		C76	56, 58
		C84	56, 58
		C86	57
		Caffeine	23, 31, 35, 93, 96, 101, 105, 107, 109, 115, 125
		Calcium Folate	104
		Camostat Mesilate	84
		Campesterol	26
		Capsaicin	27, 116, 120
		Captopril	83
		Carbamazepine	109, 110
		Carbazochrome	84
		Carbazochrome Sodium Sulfonate	84
		Carbazole	33
		Carbenicillin	111
		Carbidopa	84
		Carbonic Anhydrase	139, 141
		Carbosulfan	122
		Carnitine	125
		$\alpha$ -Carotene	118
		$\beta$ -Carotene	118, 123, 127
		Carvylacetate (cis, trans)	148
		CAT	122
		Catalase	140
		Catechin (C)	23, 35
		Catechin gallate(CG)	23, 35
		Cefaclor	90
		Cefadroxil	90
		Cefsulodin Sodium Salt Hydrate	90
		Ceftazidime	90, 91
		Ceftriaxone Sodium	91
		Cephaeline Hydrobromide	116
		Cephalothin Sodium Salt	90
		Cetirizine Dihydrochloride	89
		Chlomethoxnyl	121
		Chlormadinone Acetate	86
		Chlorobenzene	144, 147, 159
		4-Chlorobenzoic Acid	82
		p-Chlorobenzoic Acid	102
		4-Chlorobenzophenone	86
		Chloro IPC	122
		m-Chlorophenol	143
		o-Chlorophenol	143
		p-Chlorophenol(4-)	85, 86, 93
		4-(4-Chlorophenyl)-4-hydroxypiperidine	97
		Chlorotoluron	122
		Chloroxuron	122

I. COSMOSIL HPLC column	Chlorpheniramine Maleate	86, 95
	Chlorpropamide	86
	Chlortetracycline	111
	Cholesterol	26
	Choline Oxidase	139
	Chondroitinase AC-II digested	132
	Chondroitin Sulfate A	132
	Chondroitin Sulfate C	132
	Chrysazin	119
	Chrysene	148, 158
II. Liquid chromatography related product	Chrysin	23
	$\alpha$ -Chymotrypsinogen A	140, 141
	Ciclosporin (Cyclosporin)	88
	Cilastatin Sodium Salt	89
	Cimetidine	110
	Cinchonidine	84
	trans-Cinnamic Acid	117
	Cinnamyl Alcohol	36
	Cinobufagin	115
	Citric Acid	18, 125
III. Application data	L-Citrulline	120
	Clarithromycin	85
	Clindamycin Hydrochloride	85
	Clindamycin Phosphate	85
	Clobetasol Propionate	86
	Clofibrate	85, 86
	Cloperastine Hydrochloride	86
	Cloxacillin Sodium Monohydrate	85
	CNP	121
	Conalbumin	48, 140
IV. Technical note	Concanavalin A	140, 141
	Coomassie Brilliant Blue G-250	149
	Cortisone	24, 87, 129
	Cortisone Acetate	87, 101, 129
	Coumarin	119
	Creatinine	128
	m-Cresol	106
	Crotonic Acid	73
	Crude fullerenes	58
	Cyanazine	122
V. Index	Cyanocobalamin	87, 105
	Cyanocobalamin Derivative	87
	Cyanuric Acid	39
	$\alpha$ -Cyclodextrin	130
	$\beta$ -Cyclodextrin	130
	$\gamma$ -Cyclodextrin	130
	Cytidine	18, 135
	Cytochrome C	44, 138, 140
	Cytosine	18, 135
	D	
Daidzein	25	

Daidzin	25
Dansyl Amino Acids	136
Dantrolene Sodium	91
dATP	18
Daunorubicin Hydrochloride	83, 91
dCMP	135
DCMU	121
DCEPA	121
dCTP	18
Deferoxamine Mesilate	93
Dehydroacetic Acid	123
Deoxyribonuclease I	140
2'-Deoxyuridine	82
Desipramine Hydrochloride	109
Dexamethasone	92, 129
dGMP	135
dGTP	18
Diaphorase	141
Dibenz[a,h]anthracene	148
Dibenzosuberone	149, 157
Dibenzosuberone	149, 157
Dibucaine Hydrochloride	92
2,6-Di-tert-butyl-4-hydroxytoluene	124
Di-n-butyl Phthalate	96
Diclofenac Sodium	88
Diclofenamide	88
Diethylene Glycol	145
Diethyl Terephthalate	91
1,2-Difluorobenzene	144
1,3-Difluorobenzene	144
1,4-Difluorobenzene	144
Digitoxigenin	111
Digitoxin	87, 111
Digoxin	87
Dihydrocoumarin	119
Dihydroquinine	84
7,8-Dihydroxyflavone	23
3,4-Dihydroxyphenylacetic Acid	128
Diketopiperazine	123
Diltiazem Hydrochloride	89
Dimefuron	122
4-Dimethylaminoantipyrine	89
Dimethylaminophenol	97
N,N-Dimethylaniline	146, 147
N,N-Dimethylbenzylamine	9
1,5-Dimethylnaphthalene	28, 156
2,4-Dinitroaniline	147
1,2-Dinitrobenzene	83
1,5-Dinitronaphthalene(1,5-DNN)	28, 156
1,8-Dinitronaphthalene(1,8-DNN)	28, 156
Diphenhydramine Hydrochloride	110
Diphenyl	124
Diphenylimidazole	100

Diphenylmethane	149
Diphenyl Phthalate	95
Dipyridamole	88
Disulfiram	88
DNA	135
DNA-wrapped CNT	60
Dobutamine Hydrochloride	94
Domperidone	95
L-DOPA	18, 72, 128
Dopamine	128
Dopamine Hydrochloride	18, 72, 94
Doxifluridine	93
Doxorubicin Hydrochloride	83, 93
Doxycycline Hydrochloride	93
dTTP	18
E	
E.coli proteins	49
Egg White	142
Elaidic Acid	133
Emetine Hydrochloride	116
Epicatechin(EC)	23, 35
Epicatechin gallate(ECG)	23, 35
Epigallocatechin(EGC)	23, 35
Epigallocatechin gallate(EGCG)	23, 35
Epinephrine	128
Epirubicin Hydrochloride	83
meso-Erythritol	41, 130, 131, 146
Erythrosine	123
Estradiol	128
17- $\alpha$ -Estradiol	128
17- $\beta$ -Estradiol	128
Estradiol Benzoate	83
Estriol	83, 128
Estrone	128
Ethidimuron	122
Ethinylestradiol	109
Ethiofencarb	122
Ethofumesate	122
o-Ethoxybenzamide	95
4-Ethoxyphenol	85
Ethylbenzene	21, 34, 147, 160
p-Ethylbenzoic Acid	146
Ethylene Glycol	145
Ethyl p-Hydroxybenzoate(4-)	16, 84, 87, 88, 89, 95, 100, 102, 147
2-Ethylpyridine	146
Ethyl Salicylate	16, 147
Ethyl Vanillin	104

F	
FAD	100
Famotidine	99, 110
Fast Green FCF	123
Fibrinogen	140
Fisetin	23
Flavin Mononucleotide Sodium Salt	125
Flavone	23
Flopropione	102
Flumequine	112
Fluocinolone Acetonide	100, 101, 129
Fluocinonide	100
Fluoranthene	148
Fluorene	148, 149
Fluorobenzene	144, 147
Fluorodifen	121
Fluorometholone	101, 129
Fluoxymesterone	100
Flurbiprofen	101, 111
FMN	100
Folic Acid	104, 105, 125
Fructooligosaccharides	42
Fructose	41, 42, 130, 131
D-Fructose-6-phosphate	132
Fullerene toluene extract	55
Fumaric Acid	18, 98
Furan	147
Furfural	147
Furfuryl Alcohol	147
G	
Gabexate Mesilate	84
Galactose	41
Gallic Acid	16, 113
Gallocatechin(GC)	23, 35
Gallocatechin gallate(GCG)	23, 35
Gd@C82(I)	57
Geniposide	113, 115
Genistein	25
Genistin	25
Gentisic Acid	16
6-Gingerol	120
Ginkgolide A	120
Ginkgolide B	120
Ginkgolide C	120
Ginsenoside Rb1	114
Ginsenoside Rc	114
Ginsenoside Re	114
Ginsenoside Rg1	114
Gitoxin	87
Glucitol (Sorbitol)	41, 131



Glucose	40, 41, 42, 130, 131, 132
α-D-Glucose-1-phosphate	132
D-Glucose-6-phosphate	132
Glucose Oxidase	48, 139
β-Glucosidase	53
Glucuronic Acid	42
Glutamic Acid	137
L-Glutamic Dehydrogenase	139
Glyceric Acid	146
Glycerol	41, 145
Glycine	38
Glycitein	25
Glycitin	25
Glycolic Acid	18
Glycylglycine	38
Glycyrrhizic Acid	113, 124
Griseofulvin	85
Guaiacol Glycerol Ether	81, 85
Guanine	18, 135
Guanosine	18, 111, 135
<b>H</b>	
Haloperidol	97
Hemoglobin	53
Hemoglobin, Bovine	46
Hesperidin	117
Hibenzic Acid	92
Hinokitiol	119
Hippuric Acid	128
Histidine	73
Homovanillic Acid	18
Honokiol	114
Human Serum	142
Hyaluronic Acid	132
Hydralazine Hydrochloride	111
Hydrochlorothiazide	98
Hydrocortisone	24, 98, 129
Hydrocortisone Acetate	87, 98
Hydrocortisone Succinate	98
Hydroquinidine Hydrochloride	84
Hydroquinone	113
Hydroxocobalamin Acetate	105
p-Hydroxybenzoic Acid	87, 95
3-Hydroxyflavone	23, 119
5-Hydroxyflavone	23, 119
6-Hydroxyflavone	23, 119
7-Hydroxyflavone	23, 119
Hypoxanthine	129
<b>I</b>	
Ibuprofen	86, 111

Idoxuridine	82
Imazalil	124
Imipramine Hydrochloride	109
Indeno[1,2,3-c,d]pyrene	148
Indigo Carmine	123
Indometacin (Indomethacin)	82, 115
Inositol	41, 131
Insulin	151
Iodobenzene	147
5-Iodouracil	82
Ipratropium Bromide	82
Ipratropium Bromide Derivative	82
Iprodione	121
Isoamyl Benzoate	172
Isoascorbic Acid	125
Isobutyl p-Hydroxybenzoate	81, 105
Isoleucine	39
Isomaltooligosaccharides	130
Isoniazid	82
Isonicotinamide	90
Isonicotinic Acid	82
Isophthalonitrile	26, 156
Isoproc carb	122
Isopropyl p-Hydroxybenzoate	91, 100
Isoproturon	122
Isoxsuprine Hydrochloride	82
<b>J</b>	
Josamycin	89
<b>K</b>	
Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative	81
Ketamine Hydrochloride	109
Ketoprofen	87, 111
<b>L</b>	
Lactic Acid	18
β-Lactoglobulin	53
Lactose	43, 131
Latamoxef Sodium (isomer 1)	106, 107
Latamoxef Sodium (isomer 2)	106, 107
Leu-Enkephalin	138
Levallorphan Tartrate	108
Levothyroxine Sodium	109
Lidocaine	107, 109
Lincomycin Hydrochloride	85, 108
Linolenic Acid p-Bromophenacyl Ester	134
γ-Linolenic Acid p-Bromophenacyl Ester	134
Linuron	121, 122
Liothyronine Sodium	107

Lisinopril	107
Luteolin	23
Lypoxytase	51
Lysozyme	44, 140
Lysyl Endopeptidase digested	138
<b>M</b>	
Magnolol	114
Maleic Acid	18, 86, 92
Malic Acid	18, 120, 126
Maltitol	41, 130, 131
Maltoheptaose	42
Maltohexaose	42
Maltopentaose	42
Maltose	40, 41, 42, 43, 130, 131
Maltotetraitol	131
Maltotetraose	42
Maltotriose	42, 131
Mandelic Acid	128
Maneb	121
Mannitol	131
Mannosamine	42
Mannose	42
MCC	121
MCP	121
MCPB	121
MCPP	121
Mecobalamin	105
Melamine	39
Menatetrenone	106
Mepivacaine Hydrochloride	106
5-Mercapto-1-methyltetrazole	106
6-Mercaptopurine Hydrate	80
Meropenem Trihydrate	106
Met-Enkephalin	138
Metabromuron	122
Metazachlor	122
Metformin Hydrochloride	105
Methabenzthiazuron	122
Methicillin Sodium Salt	111
Methomyl	122
Methotrexate	105
DL-3-Methoxy-4-hydroxymandelic Acid	18
6-Methoxyflavone	23
m-Methoxyphenol	143
o-Methoxyphenol	143
p-Methoxyphenol	143
3-Methoxytyramine	128
Methyl Benzoate	14, 30, 83, 91, 147
N-Methylbenzylamine	9
Methyldopa	84
o-Methyl Hippuric Acid	128

Methyl p-Hydroxybenzoate	16, 85, 93, 94, 96, 97, 99, 103, 105, 108, 113
Methyl Linoleate	133
Methyl Linolenate	133
Methyl Margarate	133
1-Methylnaphthalene	28, 156
Methyl n-Nonadecanoate	133
Methyl Oleate	133
Methyl Prednisolone	100
Methyl Salicylate	16, 147
Methyl Stearate	133
17-Methyltestosterone	128
Meticrane	105
Metoxuron	122
Metronidazole	105
Mexiletine Hydrochloride	104
Milk	139
4-epi-Minocycline	104
Minocycline Hydrochloride	104
Mitomycin C	104
Mizoribine	104
Monolinuron	122
MPMC	122
MTMC	122
Myoglobin	44, 48, 51, 53, 140
Myosin	138
Myricetin	23
<b>N</b>	
Nabumetone	96
NAC	122
Nalidixic Acid	96
Naphazoline Nitrate	95
Naphthalene	28, 148, 156, 157
2-Naphthalenesulfonic Acid	91
2-Naphthalenesulfonic Acid Sodium	83
1-Naphthol	37
2-Naphthol	37
Naringenin	118
Naringin	116, 117, 118
Neostigmine Methylsulfate	97
Neurotensin	138
New Coccine	123
Niacin	125
Nicardipine Hydrochloride	96
Nicergoline	97
Nicotinamide	96, 125
Nicotinic Acid	93, 96, 125
Nifedipine	96
Nitrendipine	97
p-Nitroaniline	147
Nitrobenzene	147, 155

I. COSMOSIL HPLC column	II. Liquid chromatography related product	III. Application data	IV. Technical note	V. Index
p-Nitrobenzyl Alcohol	2-Phenylethyl-β-gentiobioside			19
Nitrofen	2-Phenylethyl-β-lactoside			19
Nitrofurantoin	2-Phenylethyl-β-maltoside			19
Nitrofurazone	2-Phenylethyl-β-melibioside			19
p-Nitrophenol	o-Phenylphenol			124
L-Noradrenaline	Phenytol			109, 110
Norepinephrine	Phloxine			123
	L-α-Phosphatidyl-L-serine			134
<b>O</b>	L-α-Phosphatidyl Choline			134
Ofloxacin	L-α-Phosphatidylethanolamine			134
Oleic Acid	L-α-Phosphatidylinositol Sodium Salt			134
Omeprazole	Phthalic Acid			16, 98
Ovalbumin	Phthalonitrile			26, 156
Oxalic Acid	Phytonadione			106
Oxamic Acid	Pimaricin			99
Oxine-copper	Piperacillin Sodium			99
Oxolinic Acid	Piperonyl Butoxide			124
Oxyfluorfen	Pirimicarb			122
Oxytetracycline	Piromidic Acid			112
Oxytocin	cis-Platin (CDDP)			111
Oxytocin Acetate Salt	Polyethylene Glycol			49
	Poly-L-Glutamic Acid Sodium Salt			45
<b>P</b>	Polyoxyethylene Hexadecyl Ether			148
PAC	Polyoxyethylene Lauryl Ether			148
Paeoniflorin	Potassium Biphthalate			94
Paeonol	Potassium Clavulanate			84
PA-Glucose Oligomer (DP=3-22)	Pravastatin Sodium			100
Palatinit	Prednisolone			24, 95, 101, 129
Palatinose	Prednisolone Acetate			101, 129
Palmitine Chloride	Prednisone			98
Palmitine Chloride Hydrate	Procainamide			14, 110
Palmitine Hydrochloride	Procaine Hydrochloride			101
D-Pantothenic Acid	Procaterol Hydrochloride			101
D-Pantothenic Acid Sodium Salt	Progesterone			83
Papaverine Hydrochloride	Prometryn			122
pd(T)12-18	Propazine			122
Penicillin G Potassium Salt	Propionic Acid			18
Perfluorotetradecanoic Acid	Propranolol Hydrochloride			102
Peroxidase	Propylbenzene			21, 34, 149, 160
Perylene	n-Propylbenzene			147
PHC	Propyl Benzoate			100
Phenacetin	Propyl Gallate			124
Phenanthrene	Propyl p-Hydroxybenzoate			16, 88, 89, 91, 92, 97, 100, 107, 113, 124, 147
Phenethylamine Hydrochloride	Propyl Salicylate			16
Phenol	Prostaglandin D2			129
Phenolphthalein	Prostaglandin E2			129
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	Prostaglandin F2α			129
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	Prostaglandin I2			129
Phenylalanine	Protocatechuic Acid			16
Phenyl Benzoate	PTH-Arg			136
2-Phenylethyl-β-cellobioside	PTH-Asp			136
	PTH-His			136
	PTH-Lys			136
	PTH-Ser			136
	PTH-Trp			136
	Puerarin			113
	Pullulan			49
	Purpurin			119
	Pyrene			148
	Pyridine			91, 145, 147
	Pyridoxine (Vitamin B6)			125
	Pyridoxine Hydrochloride			125
	Pyruvate Kinase			140
	<b>Q</b>			
	Quercetin			23
	Quinidine			109
	Quinidine Sulfate			84
	Quinine Hydrochloride			84
	Quinizarin			119
	<b>R</b>			
	Raffinose			41, 130
	Ranitidine			110
	Reserpine			108, 120
	Resibufogenin			115
	13-cis-Retinoic Acid			127
	all-trans-Retinoic Acid			127
	all-trans-Retinol Acetate			123
	Rhamnose			41, 42
	Riboflavin			100
	Riboflavin (Vitamin B2)			125
	Riboflavin Phosphate			125
	Riboflavin Phosphate Sodium Salt			125
	Ribonuclease A			44, 140
	Rifampicin			108
	Ritodrine Hydrochloride			107, 108
	threo-Ritodrine Hydrochloride			107
	Rose Bengale			123
	Roxatidine Acetate Hydrochloride			109
	Roxithromycin			109
	Rthoxyquin			124
	<b>S</b>			
	Saccharin			123
	Saikosaponin a			35, 114, 118
	Saikosaponin b1			35, 118
	Saikosaponin b2			35, 115, 118
	Saikosaponin c			35, 118
	Saikosaponin d			35, 114, 118
	Saikosaponin h			35
	Salicylamide			36, 94
	Salicylic Acid			36, 87, 106, 110
	Santonin			87
	Sc2@C76(I)			57
	Sc2@C78			57
	Sc2@C80(I)			57
	Scopolamine n-Butyl Bromide			99
	Scopolamine Hydrobromide			15, 99, 117
	Sennoside A			115, 116
	Sennoside B			115
	6-Shogaol			120
	Sitosterol			26
	Sodium Amoxicillin			84
	Sodium Copper Chlorophyll			119
	Sodium Copper Chlorophyllin			123
	Sodium 5'-Cytidylate			124
	Sodium 5'-Guanylate			124
	Sodium 5'-Inosinate			124
	Sodium Iron Chlorophyllin			123
	Sodium Polystyrene Sulfonate			49
	Sodium 5'-Uridylate			124
	Sorbic Acid			73, 123, 125
	Soybean Milk			139
	Stearic Acid			133
	Stigmasterol			26
	cis-Stilbene			144
	trans-Stilbene			144
	Streptomycin Sulfate			112
	Strychnine			117
	Styrene oxide			147
	Substance P			16, 137, 138
	Succinic Acid			18
	Sucrose			41, 43, 130, 131
	Sulbactam Sodium			89
	Sulfachloropyridazine			112
	Sulfadimethoxine			112
	Sulfaisomidin			112
	Sulfamethazine			112
	Sulfamethoxazole			112
	Sulfamethoxypyridazine			112
	Sulfamonomethoxine			112
	Sulfathiazole			82, 112
	Sunset Yellow FCF			123
	Swertiamarin			116
	<b>T</b>			
	Tartaric Acid			18
	Tartrazine			123
	Taurine			137
	Terbutryn			122
	Terbutylazine			122
	Terephthalonitrile			26, 156
	m-Terphenyl			143

o-Terphenyl	12, 21, 143
p-Terphenyl	88, 143
Testosterone	128
Testosterone Propionate	102, 128
Tetracycline	111
Theobromine	31, 109
Theophylline	31, 91, 106, 109, 116
Thermolysin	139
Thiabendazole	124
Thiamine Hydrochloride	91, 114
Thioanisole	147
Thiopental Sodium	91
Thiopental Sodium Isomer	91
Thiram	121
Thymidine	18, 135
Thymine	18, 88, 135
Thyroglobulin	48, 140
Ticarcillin Sodium	92
Tiglic Acid	73
Timolol Maleate	92
Tipepidine	92
Tocopherol	93
α-Tocopherol	27, 93, 127
β-Tocopherol	27, 127
γ-Tocopherol	27, 127
δ-Tocopherol	27, 127
Tocopherol Acetate	93
Tocopherol Nicotinate	93
α-Tocotrienol	126
β-Tocotrienol	126
γ-Tocotrienol	126
δ-Tocotrienol	126
Todalazine Hydrochloride	94, 111
Tolazamide	94
Tolbutamide	94
Tolnaftate	95
Toluene	14, 21, 30, 31, 34, 144, 147, 159, 160
o-Toluic Acid	92, 146
m-Tolunitrile	24, 155
o-Tolunitrile	24, 155
p-Tolunitrile	24, 155
N-(m-Toluoyl) glycine	128
N-(o-Toluoyl) glycine	128
N-(p-Toluoyl) glycine	128
Tranexamic Acid	94
Transferrin	51, 139, 142
Trehalose	43, 130
Triamcinolone	94, 129
Triamcinolone Acetonide	95, 101, 129
Trimetazidine Hydrochloride	95
1-(2,3,4-Trimethoxybenzyl) piperazine Dihydrochloride	95
Trimethylene Glycol	145
Triphenylene	12, 21

Tris (hydroxymethyl) aminomethane	146
Triton X-100	29
Trypsin Inhibitor	51
Trypsin Inhibitor, Soybean	141
Trypsinogen	141
Tryptophan	73
Tyrosine	73
<b>U</b>	
Ubidecarenone	106
Uracil	9, 18, 94, 109, 110, 135, 145, 146, 147
Urea	146
Uric Acid	129
Uridine	18, 135, 145
<b>V</b>	
cis-Vaccenic Acid	133
trans-Vaccenic Acid	133
Valerophenone	12
Valine	39
Vancomycin Hydrochloride	97
Vanillylnonanamide	116
N-Vanillylnonanamide	27
[Arg8]-Vasopressin	83
Vincristine Sulfate	99
Vinyl Acetate	117
2-Vinylpyridine	103
1-Vinyl-2-pyrrolidone	117
Vitamin A	127
Vitamin A Acetate	126
Vitamin B1	125
Vitamin B1 Hydrochloride	125
Vitamin B2	125
Vitamin B6	125
Vitamin B12	125
Vitamin C	125
Vitamin D2	126
Vitamin D3	126
Vitamin E	126
Vitamin E Acetate	126
Vitamin K1	126
Vitamin K3	126
<b>X</b>	
Xanthene	92
Xanthine	129
XMC	122
p-Xylene	36
Xylitol	41, 130, 131
Xylooligosaccharides	130

Xylose	42
<b>Z</b>	
Zidovudine	88

Packing material's name

C<sub>18</sub>-MS-II

Acesulfame K	123
Acetaminophen	31, 80
Acetanilide	99
Acetoaminophen	15
Acetophenone	14, 30, 147
4'-Acetoxyacetanilide	31, 80
17 $\alpha$ -Acetoxypregesterone	86
Acettrizoic Acid	81
N-Acetylprocainamide	14
Acid Red	123
Aconitine	120
Acrylic Acid	73
Adenosine	100
Alizalin	119
Amaranth	123
Amidotrizoic Acid	81
4-Aminoantipyrine	15
p-Aminophenol Hydrochloride	31, 80
Amlodipine Besilate	81
Ampicillin	99, 111
Ampicillin Sodium Salt	81
D-(-)-Amygdalin	114
Amylbenzene	12, 21, 34, 119, 160
Aniline	147
Anisole	147, 155
Anthracene	157
Anthrarufin	119
Antipyrine	15
Apigenin	118
L(+)-Ascorbic Acid	125
Aspartame	123
Asulam	121
Atrazine	122
Atropine Sulfate	15
Baicalein	23
Baicalin	115
Bamethane Sulfate	102
Barbital	110
Beclometasone Dipropionate	86, 102, 103
Bensulide	121
Benzamide	147
Benz[a]anthracene	158
Benzene	14, 21, 30, 31, 34, 144, 146, 147, 155, 157, 160
Benzoic Acid	106, 123, 146
Benzonitrile	147
Benzophenone	88
Benzyl Alcohol	9, 14

Benzyl p-Hydroxybenzoate	98
Benzympenicillin Potassium	103
Berberine Hydrochloride	25
Betahistine Mesilate	103
Betamethasone	103, 129
Betamethasone Sodium Phosphate	103
Betamethasone Valerate	103
Betamethasone 17-Valerate	172
D-Biotin [Vitamin H]	125
Boc-(D)Phe-(L)Phe-OMe	137
Boc-(L)Phe-(L)Phe-OMe	137
BPMC	122
Brilliant Blue FCF	123
Bromacil	122
Bromhexine Hydrochloride	102
Bromobenzene	144, 147, 159
Butylbenzene	12, 21, 34, 160
n-Butylbenzene	147
Butylbenzoate	12
2-tert-Butylhydroquinone	124
3-tert-Butyl-4-hydroxyanisol	124
Butyl p-Hydroxybenzoate	82, 84, 85, 86, 88, 98, 101, 103, 117
C60	54
C70	54
Caffeine	23, 31, 101, 105, 107, 115
Calcium Folate	104
Campesterol	26
Capsaicin	27, 120
Captopril	83
Carbamazepine	109, 110
Carbenicillin	111
Carbosulfan	122
$\alpha$ -Carotene	118
$\beta$ -Carotene	118
CAT	122
Catechin (C)	23
Catechin gallate(CG)	23
Cefaclor	90
Cefadroxil	90
Ceftazidime	91
Ceftriaxone Sodium	91
Cephaeline Hydrobromide	116
Cephalothin Sodium Salt	90
Chlormethoxynil	121
Chlormadinone Acetate	86
Chlorobenzene	144, 147, 159
4-Chlorobenzoic Acid	82
Chloro IPC	122
m-Chlorophenol	143
o-Chlorophenol	143

p-Chlorophenol	143
Chlorotoluron	122
Chloroxuron	122
Chlortetracycline	111
Cholesterol	26
Chrysazin	119
Chrysene	158
Chrysin	23
Cimetidine	110
Clarithromycin	85
Clindamycin Hydrochloride	85
Clobetasol Propionate	86
Cloxacillin Sodium Monohydrate	85
CNP	121
Coomassie Brilliant Blue G-250	149
Cortisone	24, 87, 129
Cortisone Acetate	87, 101, 129
Coumarin	119
Creatinine	128
m-Cresol	106
Crotonic Acid	73
Cyanazine	122
Daidzein	25
Daidzin	25
Dansyl Amino Acids	136
Daunorubicin Hydrochloride	91
DCMU	121
DCPA	121
Dehydroacetic Acid	123
Desipramine Hydrochloride	109
Dexamethasone	129
Dibenzosuberone	157
Dibenzosuberone	157
2,6-Di-tert-butyl-4-hydroxytoluene	124
Di-n-butyl Phthalate	96
Diclofenamide	88
Diethyl Terephthalate	91
1,2-Difluorobenzene	144
1,3-Difluorobenzene	144
1,4-Difluorobenzene	144
Dihydrocoumarin	119
7,8-Dihydroxyflavone	23
Diketopiperazine	123
Diltiazem Hydrochloride	89
Dimefuron	122
Dimethylaminophenol	97
N,N-Dimethylaniline	147
N,N-Dimethylbenzylamine	9
1,5-Dimethylnaphthalene	28, 156
2,4-Dinitroaniline	147
1,5-Dinitronaphthalene(1,5-DNN)	28, 156
1,8-Dinitronaphthalene(1,8-DNN)	28, 156
Diphenhydramine Hydrochloride	110

Diphenyl	124
Disulfiram	88
Dobutamine Hydrochloride	94
Dopamine Hydrochloride	94
Doxorubicin Hydrochloride	93
Doxycycline Hydrochloride	93
Elaidic Acid	133
Emetine Hydrochloride	116
Epicatechin(EC)	23
Epicatechin gallate(EGC)	23
Epigallocatechin(EGC)	23
Epigallocatechin gallate(EGCG)	23
Erythrosine	123
17- $\alpha$ -Estradiol	128
17- $\beta$ -Estradiol	128
Estradiol Benzoate	83
Estriol	83
Ethidimuron	122
Ethiofencarb	122
Ethofumesate	122
Ethylbenzene	21, 34, 147, 160
p-Ethylbenzoic Acid	146
Ethyl p-Hydroxybenzoate	84, 87, 89, 100, 147
Ethyl Salicylate	147
FAD	100
Famotidine	110
Fast Green FCF	123
Fisetin	23
Flavin Mononucleotide Sodium Salt	125
Flavone	23
Flumequine	112
Fluocinolone Acetonide	129
Fluocinonide	100
Fluorobenzene	144, 147
Fluorodifen	121
Fluorometholone	101, 129
Flurbiprofen	111
FMN	100
Folic Acid	104, 105, 125
Fullerene toluene extract	55
Furan	147
Furfural	147
Furfuryl Alcohol	147
Gabexate Mesilate	84
Gallocatechin(GC)	23
Gallocatechin gallate(GCG)	23
Geniposide	115
Genistein	25
Genistin	25
6-Gingerol	120
Ginkgolide A	120
Ginkgolide B	120
Ginkgolide C	120

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

Ginsenoside Rb1 ..... 114  
 Ginsenoside Rc ..... 114  
 Glycitein ..... 25  
 Glycitin ..... 25  
 Griseofulvin ..... 85  
 Guaiacol Glycerol Ether ..... 81, 85  
 Hibenzic Acid ..... 92  
 Hippuric Acid ..... 128  
 Histidine ..... 73  
 Hydralazine Hydrochloride ..... 111  
 Hydrocortisone ..... 24, 98, 129  
 Hydrocortisone Acetate ..... 87, 98  
 Hydrocortisone Succinate ..... 98  
 p-Hydroxybenzoic Acid ..... 95  
 3-Hydroxyflavone ..... 23, 119  
 5-Hydroxyflavone ..... 23, 119  
 6-Hydroxyflavone ..... 23, 119  
 7-Hydroxyflavone ..... 23, 119  
 Ibuprofen ..... 111  
 Imazalil ..... 124  
 Imipramine Hydrochloride ..... 109  
 Indigo Carmine ..... 123  
 Indometacin (Indomethacin) ..... 82  
 Iodobenzene ..... 147  
 Iprodione ..... 121  
 Isoamyl Benzoate ..... 172  
 Isobutyl p-Hydroxybenzoate ..... 81  
 Isoniazid ..... 82  
 Isonicotinic Acid ..... 82  
 Isophthalonitrile ..... 26, 156  
 Isoprocarb ..... 122  
 Isopropyl p-Hydroxybenzoate ..... 91  
 Isoproturon ..... 122  
 Isoxsuprine Hydrochloride ..... 82  
 Ketamine Hydrochloride ..... 109  
 Ketoprofen ..... 111  
 Latamoxef Sodium (isomer 1) ..... 106  
 Latamoxef Sodium (isomer 2) ..... 106  
 Lidocaine ..... 109  
 Linuron ..... 121, 122  
 Lisinopril ..... 107  
 Luteolin ..... 23  
 Mandelic Acid ..... 128  
 Maneb ..... 121  
 MCC ..... 121  
 MCP ..... 121  
 MCPB ..... 121  
 MCPP ..... 121  
 Menatetrenone ..... 106  
 5-Mercapto-1-methyltetrazole ..... 106  
 Metabromuron ..... 122  
 Metazachlor ..... 122  
 Methabenzthiazuron ..... 122

Methicillin Sodium Salt ..... 111  
 Methomyl ..... 122  
 Methotrexate ..... 105  
 6-Methoxyflavone ..... 23  
 m-Methoxyphenol ..... 143  
 o-Methoxyphenol ..... 143  
 p-Methoxyphenol ..... 143  
 Methyl Benzoate ..... 14, 30, 83, 147  
 N-Methylbenzylamine ..... 9  
 o-Methyl Hippuric Acid ..... 128  
 Methyl p-Hydroxybenzoate ..... 94, 103, 105  
 Methyl Linoleate ..... 133  
 Methyl Linolenate ..... 133  
 Methyl Margarate ..... 133  
 1-Methylnaphthalene ..... 28, 156  
 Methyl n-Nonadecanoate ..... 133  
 Methyl Oleate ..... 133  
 Methyl Salicylate ..... 147  
 Methyl Stearate ..... 133  
 Meticrane ..... 105  
 Metoxuron ..... 122  
 Metronidazole ..... 105  
 Monolinuron ..... 122  
 MPMC ..... 122  
 MTMC ..... 122  
 Myricetin ..... 23  
 NAC ..... 122  
 Naphthalene ..... 28, 156, 157  
 2-Naphthalenesulfonic Acid ..... 91  
 Naringenin ..... 118  
 Naringin ..... 116, 118  
 Neostigmine Methylsulfate ..... 97  
 New Coccine ..... 123  
 Niacin ..... 125  
 Nicardipine Hydrochloride ..... 96  
 Nicergoline ..... 97  
 Nicotinamide ..... 96, 125  
 Nicotinic Acid ..... 93, 96, 125  
 Nifedipine ..... 96  
 Nitrendipine ..... 96  
 p-Nitroaniline ..... 147  
 Nitrobenzene ..... 147, 155  
 Nitrofen ..... 121  
 Ofloxacin ..... 112  
 Oleic Acid ..... 133  
 Oxine-copper ..... 14  
 Oxolinic Acid ..... 112  
 Oxyfluorfen ..... 121  
 Oxytetracycline ..... 111  
 Oxytocin Acetate Salt ..... 83  
 PAC ..... 122  
 Paeonol ..... 8, 117  
 Palmatine Hydrochloride ..... 25

D-Pantothenic Acid Sodium Salt ..... 125  
 Papaverine Hydrochloride ..... 120  
 Penicillin G Potassium Salt ..... 111  
 PHC ..... 122  
 Phenacetin ..... 15  
 Phenol ..... 106, 145, 147  
 Phenolphthalein ..... 147  
 Phenylalanine ..... 73  
 Phenyl Benzoate ..... 89  
 o-Phenylphenol ..... 124  
 Phenytoin ..... 109, 110  
 Phloxine ..... 123  
 L- $\alpha$ -Phosphatidyl-L-serine ..... 134  
 L- $\alpha$ -Phosphatidyl Choline ..... 134  
 L- $\alpha$ -Phosphatidylethanolamine ..... 134  
 L- $\alpha$ -Phosphatidylinositol Sodium Salt ..... 134  
 Phthalonitrile ..... 26, 156  
 Phytonadione ..... 106  
 Pimaricin ..... 99  
 Piperacillin Sodium ..... 99  
 Pirimicarb ..... 122  
 Piromidic Acid ..... 112  
 Potassium Biphthalate ..... 94  
 Pravastatin Sodium ..... 100  
 Prednisolone ..... 24, 95, 101, 129  
 Prednisolone Acetate ..... 101, 129  
 Procaïnamide ..... 14  
 Procaine Hydrochloride ..... 101  
 Progesterone ..... 83  
 Prometryn ..... 122  
 Propazine ..... 122  
 Propranolol Hydrochloride ..... 102  
 Propylbenzene ..... 21, 34, 160  
 n-Propylbenzene ..... 147  
 Propyl Benzoate ..... 100  
 Propyl Gallate ..... 124  
 Propyl p-Hydroxybenzoate ..... 91, 92, 97, 147  
 Prostaglandin D2 ..... 129  
 Prostaglandin E2 ..... 129  
 Prostaglandin F2 $\alpha$  ..... 129  
 Prostaglandin I2 ..... 129  
 Puerarin ..... 113  
 Purpurin ..... 119  
 Pyridine ..... 91, 145, 147  
 Pyridoxine Hydrochloride ..... 125  
 Quercetin ..... 23  
 Quinidine ..... 109  
 Quinizarin ..... 119  
 Ranitidine ..... 110  
 Reserpine ..... 120  
 13-cis-Retinoic Acid ..... 127  
 all-trans-Retinoic Acid ..... 127  
 Riboflavin ..... 100

Riboflavin Phosphate Sodium Salt ..... 125  
 Rose Bengale ..... 123  
 Roxatidine Acetate Hydrochloride ..... 109  
 Rthoxyquin ..... 124  
 Saccharin ..... 123  
 Saikosaponin a ..... 114, 118  
 Saikosaponin b1 ..... 118  
 Saikosaponin b2 ..... 115, 118  
 Saikosaponin c ..... 118  
 Saikosaponin d ..... 114, 118  
 Salicylamide ..... 97  
 Salicylic Acid ..... 106  
 Santonin ..... 87  
 Scopolamine Hydrobromide ..... 15  
 Sennoside A ..... 115, 116  
 Sennoside B ..... 115  
 6-Shogaol ..... 120  
 Sitosterol ..... 26  
 Sodium Copper Chlorophyll ..... 119  
 Sorbic Acid ..... 73, 123  
 Stearic Acid ..... 133  
 Stigmasterol ..... 26  
 cis-Stilbene ..... 144  
 trans-Stilbene ..... 144  
 Streptomycin Sulfate ..... 112  
 Styrene oxide ..... 147  
 Sulbactam Sodium ..... 89  
 Sulfachloropyridazine ..... 112  
 Sulfadimethoxine ..... 112  
 Sulfaisomidin ..... 112  
 Sulfamethazine ..... 112  
 Sulfamethoxazole ..... 112  
 Sulfamethoxypyridazine ..... 112  
 Sulfamonomethoxine ..... 112  
 Sulfathiazole ..... 112  
 Sunset Yellow FCF ..... 123  
 Swertiamarin ..... 116  
 Tartrazine ..... 123  
 Terbutryn ..... 122  
 Terbutylazine ..... 122  
 Terephthalonitrile ..... 26, 156  
 m-Terphenyl ..... 143  
 o-Terphenyl ..... 12, 21, 143  
 p-Terphenyl ..... 143  
 Testosterone Propionate ..... 102  
 Tetracycline ..... 111  
 Theobromine ..... 31  
 Theophylline ..... 31, 106, 116  
 Thiabendazole ..... 124  
 Thioanisole ..... 147  
 Thiopental Sodium ..... 91  
 Thiopental Sodium Isomer ..... 91  
 Thiram ..... 121

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index



I. COSMOSIL HPLC column	II. Liquid chromatography related product	III. Application data	IV. Technical note	V. Index
Thymine ······ 88				
Tiglic Acid ······ 73				
Tipepidine ······ 92				
Tocopherol ······ 93				
α-Tocopherol ······ 27, 93				
β-Tocopherol ······ 27				
γ-Tocopherol ······ 27				
δ-Tocopherol ······ 27				
Tocopherol Acetate ······ 93				
Tocopherol Nicotinate ······ 93				
α-Tocotrienol ······ 126				
β-Tocotrienol ······ 126				
γ-Tocotrienol ······ 126				
δ-Tocotrienol ······ 126				
Todalazine Hydrochloride ······ 94, 111				
Toluene ······ 14, 21, 30, 31, 34, 144, 147, 159, 160				
o-Toluic Acid ······ 146				
m-Tolunitrile ······ 24, 155				
o-Tolunitrile ······ 24, 155				
p-Tolunitrile ······ 24, 155				
Triamcinolone ······ 94, 129				
Triamcinolone Acetonide ······ 94, 129				
Trimetazidine Hydrochloride ······ 95				
1-(2,3,4-Trimethoxybenzyl) piperazine Dihydrochloride ······ 95				
Triphenylene ······ 21				
Triton X-100 ······ 29				
Tryptophan ······ 73				
Tyrosine ······ 73				
Uracil ······ 9, 94, 110, 146, 147				
cis-Vaccenic Acid ······ 133				
trans-Vaccenic Acid ······ 133				
Valerophenone ······ 12				
Vancomycin Hydrochloride ······ 97				
N-Vanillylnonanamide ······ 27				
[Arg8]-Vasopressin ······ 83				
2-Vinylpyridine ······ 103				
Vitamin A Acetate ······ 126				
Vitamin B1 Hydrochloride ······ 125				
Vitamin B2 [Riboflavin] ······ 125				
Vitamin B6 ······ 125				
Vitamin B12 ······ 125				
Vitamin C ······ 125				
Vitamin D2 ······ 126				
Vitamin D3 ······ 126				
Vitamin E ······ 126				
Vitamin E Acetate ······ 126				
Vitamin K1 ······ 126				
Vitamin K3 ······ 126				
Xanthene ······ 92				
XMC ······ 122				
Zidovudine ······ 88				
	<b>C<sub>18</sub>-AR-II</b>			
	Acenaphthene ······ 148			
	Acenaphthylene ······ 148			
	Acetaminophen ······ 110			
	Acetanilide ······ 90			
	N-Acetylprocainamide ······ 110			
	Acetylsalicylic Acid (Aspirine) ······ 110			
	Acid Red ······ 123			
	Albiflorin ······ 115			
	Allopurinol ······ 129			
	Amaranth ······ 123			
	4-Aminoacetophenone ······ 98			
	Aminoantipyrine ······ 110			
	4-Aminobenzoic Acid ······ 80			
	4-Amino-6-chlorobenzene-1,3-disulfonamide ······ 98			
	Amylbenzene ······ 12			
	4-Androstene-3,17-dione ······ 128			
	Angiotensin I (Human) ······ 16			
	Angiotensin II (Human) ······ 16, 138			
	Angiotensin II, [Asn1,Val5] ······ 138			
	Angiotensin II, [Sar1,Ala8] ······ 138			
	Angiotensin II, [Sar1,Ile8] ······ 138			
	Angiotensin II, [Sar1,Thr8] ······ 138			
	Angiotensin II, [Val5] ······ 138			
	Angiotensin II, Des-Asp1-[Ile8] ······ 138			
	Anthracene ······ 22, 148			
	Antipyrine ······ 110			
	Atenolol ······ 81			
	Azathioprine ······ 80			
	Aztreonam ······ 80			
	Baclofen ······ 97			
	Baicalin ······ 113			
	Benz[a]anthracene ······ 148			
	Benzo[b]fluoranthene ······ 148			
	Benzoic Acid ······ 80			
	Benzo[g,h,i]perylene ······ 148			
	Benzophenone ······ 86, 106, 107			
	Benzo[a]pyrene ······ 148			
	Benz[a]pyrene ······ 148			
	Benzyl Alcohol ······ 106			
	Bezafibrate ······ 102			
	Bibenzyl ······ 22			
	Biphenyl ······ 97			
	Brilliant Blue FCF ······ 123			
	Bufalin ······ 115			
	Bufexamac ······ 100			
	Butylbenzene ······ 12			
	Butylbenzoate ······ 12			
	Butyl p-Hydroxybenzoate ······ 16, 84, 108			
	iso-Butyl p-Hydroxybenzoate ······ 108			
	Butyl Salicylate ······ 16			
	Caffeine ······ 93, 96, 109			
	Camostat Mesilate ······ 84			
	Carbazole ······ 33			
	Ceftazidime ······ 90			
	p-Chlorobenzoic Acid ······ 102			
	4-Chlorobenzophenone ······ 86			
	4-Chlorophenol ······ 86			
	4-(4-Chlorophenyl)-4-hydroxypiperidine ······ 97			
	Chlorpheniramine Maleate ······ 95			
	Chlorpropamide ······ 86			
	Chrysene ······ 148			
	Cilastatin Sodium Salt ······ 89			
	trans-Cinnamic Acid ······ 117			
	Cinobufagin ······ 115			
	Clofibrate ······ 86			
	Cloperastine Hydrochloride ······ 86			
	dAMP ······ 135			
	dCMP ······ 135			
	dGMP ······ 135			
	Cytochrome C ······ 138			
	Deferoxamine Mesilate ······ 93			
	Dibenz[a,h]anthracene ······ 148			
	Dibucaine Hydrochloride ······ 92			
	Diclofenac Sodium ······ 88			
	Digitoxigenin ······ 111			
	Digitoxin ······ 111			
	3,4-Dihydroxyphenylacetic Acid ······ 128			
	N,N-Dimethylaniline ······ 146			
	Diphenylimidazole ······ 100			
	L-DOPA ······ 128			
	Dopamine ······ 128			
	Doxifluridine ······ 93			
	Epinephrine ······ 128			
	Erythrosine ······ 123			
	Estradiol ······ 128			
	Estriol ······ 128			
	Estrone ······ 128			
	Ethinylestradiol ······ 109			
	o-Ethoxybenzamide ······ 95			
	Ethyl p-Hydroxybenzoate(4-) ······ 16, 88, 100, 102			
	2-Ethylpyridine ······ 146			
	Ethyl Salicylate ······ 16			
	Famotidine ······ 99			
	Fast Green FCF ······ 123			
	Flopropione ······ 102			
	Fluocinolone Acetonide ······ 100			
	Fluoranthene ······ 148			
	Fluorene ······ 148			
	Flurbiprofen ······ 101			
	Gallic Acid ······ 16			
	Geniposide ······ 113			
	Gentisic Acid ······ 16			
	Glycyrrhizic Acid ······ 113, 124			
	Haloperidol ······ 97			
	Hesperidin ······ 117			
	Hibenzic Acid ······ 92			
	Hinokitiol ······ 119			
	Honokiol ······ 114			
	Hydrochlorothiazide ······ 98			
	p-Hydroxybenzoic Acid ······ 87			
	Hypoxanthine ······ 129			
	Ibuprofen ······ 86			
	Indeno[1,2,3-c,d]pyrene ······ 148			
	Indigo Carmine ······ 123			
	Indometacin (Indomethacin) ······ 115			
	Isobutyl p-Hydroxybenzoate ······ 105			
	Isopropyl p-Hydroxybenzoate ······ 100			
	Josamycin ······ 89			
	Ketoprofen ······ 87			
	Levallorphan Tartrate ······ 108			
	Levothyroxine Sodium ······ 109			
	Lidocaine ······ 107			
	Lysyl Endopeptidase digested ······ 138			
	Magnolol ······ 114			
	Mepivacaine Hydrochloride ······ 106			
	6-Mercaptopurine Hydrate ······ 80			
	Meropenem Trihydrate ······ 106			
	Metformin Hydrochloride ······ 105			
	3-Methoxytyramine ······ 128			
	Methyl Benzoate ······ 91			
	Methyl p-Hydroxybenzoate ······ 16, 93, 96, 97, 99, 108, 113			
	Methyl Salicylate ······ 16			
	17-Methyltestosterone ······ 128			
	Nabumetone ······ 96			
	Nalidixic Acid ······ 96			
	Naphazoline Nitrate ······ 95			
	Naphthalene ······ 148			
	Naringin ······ 117			
	New Coccine ······ 123			
	Nicotinic Acid ······ 96			
	Nitrofurantoin ······ 112			
	Nitrofurazone ······ 112			
	p-Nitrophenol ······ 102			
	Norepinephrine ······ 128			
	Oxytocin ······ 16			
	Paeoniflorin ······ 115			
	Perfluorotetradecanoic Acid ······ 146			
	Phenacetin ······ 110			
	Phenanthrene ······ 148			
	Phenol ······ 87, 109			
	Phloxine ······ 123			
	Phthalic Acid ······ 16, 98			
	Piperonyl Butoxide ······ 124			
	Poly-L-Glutamic Acid Sodium Salt ······ 45			
	Procainamide ······ 110			
	Procaterol Hydrochloride ······ 101			
	Propyl p-Hydroxybenzoate ······ 16, 88, 100, 113, 124			

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

Propyl Salicylate ..... 16  
 Protocatechuic Acid ..... 16  
 PTH-Arg ..... 136  
 PTH-Asp ..... 136  
 PTH-His ..... 136  
 PTH-Lys ..... 136  
 PTH-Ser ..... 136  
 PTH-Trp ..... 136  
 Pyrene ..... 148  
 Reserpine ..... 108  
 Resibufogenin ..... 115  
 Ritodrine Hydrochloride ..... 107, 108  
 threo-Ritodrine Hydrochloride ..... 107  
 Rose Bengale ..... 123  
 Salicylic Acid ..... 87, 110  
 Sodium Copper Chlorophyllin ..... 123  
 Sodium 5'-Cytidylate ..... 124  
 Sodium 5'-Guanylate ..... 124  
 Sodium 5'-Inosinate ..... 124  
 Sodium Iron Chlorophyllin ..... 123  
 Sodium 5'-Uridylate ..... 124  
 Substance P ..... 16  
 Sunset Yellow FCF ..... 123  
 Tartrazine ..... 123  
 o-Terphenyl ..... 12  
 Testosterone ..... 128  
 Testosterone Propionate ..... 128  
 Theobromine ..... 109  
 Theophylline ..... 109  
 Thiamine Hydrochloride ..... 91  
 Tipepidine ..... 92  
 Uracil ..... 109  
 Uric Acid ..... 129  
 Valerophenone ..... 12  
 Xanthine ..... 129

**C<sub>18</sub>-PAQ**

Acenaphthene ..... 87  
 Acetic Acid ..... 18  
 Acrinol ..... 80  
 Adenine ..... 18, 135  
 Adenosine ..... 18, 135  
 5'-ADP ..... 135  
 L-Adrenaline ..... 18, 72  
 Amikacin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative ..... 81  
 4-Aminoacetophenone ..... 91  
 4-Amino-n-butyric Acid [GABA] ..... 137  
 4-(Aminomethyl) benzoic Acid ..... 94  
 m-Aminophenol ..... 147  
 o-Aminophenol ..... 147  
 p-Aminophenol ..... 147  
 Amoxicillin ..... 81

5'-AMP ..... 135  
 Amylbenzene ..... 12  
 Arbutin ..... 113  
 Ascorbic Acid ..... 19  
 Aspartic Acid ..... 137  
 5'-ATP ..... 135  
 Atropine Sulfate ..... 116, 117  
 Barbital Sodium ..... 117  
 Berberine Chloride ..... 103, 113  
 Brucine Dihydrate ..... 116, 117  
 Butylbenzene ..... 12  
 Butylbenzoate ..... 12  
 Caffeine ..... 107, 125  
 Carbazochrome ..... 84  
 Carbazochrome Sodium Sulfonate ..... 84  
 Carbidopa ..... 84  
 Carnitine ..... 125  
 Cefaclor ..... 90  
 Cefsulodin Sodium Salt Hydrate ..... 90  
 Chlorpheniramine Maleate ..... 86  
 Ciclosporin (Cyclosporin) ..... 88  
 Cinchonidine ..... 84  
 Citric Acid ..... 18, 125  
 Cyanocobalamin ..... 105  
 Cytidine ..... 18  
 Cytosine ..... 18  
 dATP ..... 18  
 dCTP ..... 18  
 2'-Deoxyuridine ..... 82  
 dGTP ..... 18  
 Digitoxin ..... 87  
 Digoxin ..... 87  
 Dihydroquinine ..... 84  
 Diphenyl Phthalate ..... 95  
 L-DOPA ..... 18, 72  
 Dopamine Hydrochloride ..... 18, 72  
 dTTP ..... 18  
 Fumaric Acid ..... 18  
 Gallic Acid ..... 113  
 Ginsenoside Re ..... 114  
 Ginsenoside Rg1 ..... 114  
 Gitoxin ..... 87  
 Glutamic Acid ..... 137  
 Glycine ..... 38  
 Glycolic Acid ..... 18  
 Glycylglycine ..... 38  
 Guanine ..... 18  
 Guanosine ..... 18  
 Homovanillic Acid ..... 18  
 Hydroquinidine Hydrochloride ..... 84  
 Hydroquinone ..... 113  
 Hydroxocobalamin Acetate ..... 105  
 Idoxuridine ..... 92

5-Iodouracil ..... 92  
 Isonicotinamide ..... 90  
 Kanamycin Sulfate 2,4,6-trinitrobenzenesulfonic Acid Derivative ..... 81  
 Lactic Acid ..... 18  
 Latamoxef Sodium (isomer 1) ..... 107  
 Latamoxef Sodium (isomer 2) ..... 107  
 Liothyronine Sodium ..... 107  
 Lisinopril ..... 107  
 Maleic Acid ..... 18, 86  
 Malic Acid ..... 18  
 Mecobalamin ..... 105  
 DL-3-Methoxy-4-hydroxymandelic Acid ..... 18  
 Methyl dopa ..... 84  
 Mizoribine ..... 104  
 Nicotinamide ..... 125  
 L-Noradrenaline ..... 18, 72  
 PA-Glucose Oligomer (DP=3-22) ..... 132  
 Palmatine Chloride ..... 103  
 Palmatine Chloride Hydrate ..... 113  
 Phenacetin ..... 114  
 2-Phenylethyl-β-cellobioside ..... 19  
 2-Phenylethyl-β-gentiobioside ..... 19  
 2-Phenylethyl-β-lactoside ..... 19  
 2-Phenylethyl-β-maltoside ..... 19  
 2-Phenylethyl-β-melibioside ..... 19  
 Potassium Clavulanate ..... 84  
 Propionic Acid ..... 18  
 Propyl p-Hydroxybenzoate ..... 107  
 Quinidine Sulfate ..... 84  
 Quinine Hydrochloride ..... 84  
 Riboflavin Phosphate ..... 125  
 Roxithromycin ..... 109  
 Scopolamine Hydrobromide ..... 117  
 Sennoside A ..... 116  
 Sodium Amoxicillin ..... 84  
 Strychnine ..... 117  
 Succinic Acid ..... 18  
 Sulfathiazole ..... 82  
 Tartaric Acid ..... 18  
 o-Terphenyl ..... 12  
 Thiamine Hydrochloride ..... 114  
 Thymidine ..... 18  
 Thymine ..... 18  
 Tolnaftate ..... 95  
 Tranexamic Acid ..... 94  
 Triphenylene ..... 12  
 Ubidecarenone ..... 106  
 Uracil ..... 18  
 Uridine ..... 18  
 Valerophenone ..... 12  
 Vitamin B1 ..... 125  
 Vitamin B6 ..... 125

**Cholester**

Amylbenzene ..... 21, 160  
 Angiotensin I (Human) ..... 137  
 Angiotensin II (Human) ..... 137  
 Anthracene ..... 22  
 Apigenin ..... 118  
 Baicalein ..... 23  
 Benz[a]anthracene ..... 158  
 Benzene ..... 21, 160  
 Bibenzyl ..... 22  
 Bromobenzene ..... 144  
 Butylbenzene ..... 21, 160  
 Caffeine ..... 23, 35  
 α-Carotene ..... 118  
 β-Carotene ..... 118  
 Catechin (C) ..... 23, 35  
 Catechin gallate(CG) ..... 23, 35  
 Chlorobenzene ..... 144  
 m-Chlorophenol ..... 143  
 o-Chlorophenol ..... 143  
 p-Chlorophenol ..... 143  
 Chondroitinase AC-II digested ..... 132  
 Chondroitin Sulfate A ..... 132  
 Chondroitin Sulfate C ..... 132  
 Chrysene ..... 158  
 Chrysin ..... 23  
 7,8-Dihydroxyflavone ..... 23  
 Elaidic Acid ..... 133  
 Epicatechin(EC) ..... 23, 35  
 Epicatechin gallate(EGC) ..... 23, 35  
 Epigallocatechin(EGC) ..... 23, 35  
 Epigallocatechin gallate(EGCG) ..... 23, 35  
 Ethylbenzene ..... 21, 160  
 Fisetin ..... 23  
 Flavone ..... 23  
 Gallocatechin(GC) ..... 23, 35  
 Gallocatechin gallate(GCG) ..... 23, 35  
 3-Hydroxyflavone ..... 23  
 5-Hydroxyflavone ..... 23  
 6-Hydroxyflavone ..... 23  
 7-Hydroxyflavone ..... 23  
 Linolenic Acid p-Bromophenacyl Ester ..... 134  
 γ-Linolenic Acid p-Bromophenacyl Ester ..... 134  
 Luteolin ..... 23  
 6-Methoxyflavone ..... 23  
 m-Methoxyphenol ..... 143  
 o-Methoxyphenol ..... 143  
 p-Methoxyphenol ..... 143  
 Methyl Linoleate ..... 133  
 Methyl Linolenate ..... 133  
 Methyl Margarate ..... 133  
 Methyl n-Nonadecanoate ..... 133

I. COSMOSIL HPLC column

II. Liquid chromatography related product

III. Application data

IV. Technical note

V. Index

I. COSMOSIL HPLC column	II. Liquid chromatography related product	III. Application data	IV. Technical note	V. Index
Methyl Oleate ..... 133	1,4-Difluorobenzene ..... 144	Campesterol ..... 26	[6,6]-Phenyl-C61 Butyric Acid Methyl Ester ..... 59	
Methyl Stearate ..... 133	Fluocinolone Acetonide ..... 129	Capsaicin ..... 27	Phthalonitrile ..... 156	
Myricetin ..... 23	Fluorobenzene ..... 144	Carvylacetate (cis, trans) ..... 148	Propylbenzene ..... 160	
Naringenin ..... 118	Fluorometholone ..... 129	Chlorobenzene ..... 159	Terephthalonitrile ..... 156	
Naringin ..... 118	Genistein ..... 25	Cholesterol ..... 26	Toluene ..... 160	
Oleic Acid ..... 133	Genistin ..... 25	Chrysene ..... 158	m-Tolunitrile ..... 155	
Oxytocin ..... 137	Glycitein ..... 25	Dibenzosuberone ..... 157	o-Tolunitrile ..... 155	
Phenol ..... 145	Glycitin ..... 25	Dibenzosuberone ..... 157	p-Tolunitrile ..... 155	
Propylbenzene ..... 21, 160	Guanine ..... 135	1,5-Dimethylnaphthalene ..... 28, 156		
Pyridine ..... 145	Guanosine ..... 135	1,5-Dinitronaphthalene(1,5-DNN) ..... 28, 156	<b>PBB-R</b>	
Quercetin ..... 23	Hydrocortisone ..... 24, 129	1,8-Dinitronaphthalene(1,8-DNN) ..... 28, 156	Allylbenzene ..... 149	
13-cis-Retinoic Acid ..... 127	Isophthalonitrile ..... 26	Ethylbenzene ..... 160	Amylbenzene ..... 160	
all-trans-Retinoic Acid ..... 127	Nitrobenzene ..... 155	Isophthalonitrile ..... 26, 156	Anthracene ..... 157	
Saikosaponin a ..... 35, 118	Palmatine Hydrochloride ..... 25	1-Methylnaphthalene ..... 28, 156	Benzene ..... 157, 160	
Saikosaponin b1 ..... 35, 118	Phthalonitrile ..... 26	Naphthalene ..... 28, 156, 157	1,1'-Binaphthyl ..... 149	
Saikosaponin b2 ..... 35, 118	Prednisolone ..... 24, 129	Nitrobenzene ..... 155	Bromobenzene ..... 159	
Saikosaponin c ..... 35, 118	Prednisolone Acetate ..... 129	Phthalonitrile ..... 26, 156	Butylbenzene ..... 160	
Saikosaponin d ..... 35, 118	Sitosterol ..... 26	Propylbenzene ..... 160	Chlorobenzene ..... 159	
Saikosaponin h ..... 35	Stigmasterol ..... 26	Sitosterol ..... 26	Dibenzosuberone ..... 149, 157	
Stearic Acid ..... 133	Terephthalonitrile ..... 26	Stigmasterol ..... 26	Dibenzosuberone ..... 149, 157	
cis-Stilbene ..... 144	Thymidine ..... 135	Terephthalonitrile ..... 26, 156	Diphenylmethane ..... 149	
trans-Stilbene ..... 144	Thymine ..... 135	$\alpha$ -Tocopherol ..... 27	Ethylbenzene ..... 160	
Substance P ..... 137	$\alpha$ -Tocopherol ..... 27	$\beta$ -Tocopherol ..... 27	Fluorene ..... 149	
m-Terphenyl ..... 143	$\beta$ -Tocopherol ..... 27	$\gamma$ -Tocopherol ..... 27	Naphthalene ..... 157	
o-Terphenyl ..... 21, 143	$\gamma$ -Tocopherol ..... 27	$\delta$ -Tocopherol ..... 27	Perylene ..... 149	
p-Terphenyl ..... 143	$\delta$ -Tocopherol ..... 27	Toluene ..... 159, 160	Polyoxyethylene Hexadecyl Ether ..... 148	
Toluene ..... 21, 144, 160	$\alpha$ -Tocotrienol ..... 126	m-Tolunitrile ..... 155	Polyoxyethylene Lauryl Ether ..... 148	
Triphenylene ..... 21	$\beta$ -Tocotrienol ..... 126	o-Tolunitrile ..... 155	Propylbenzene ..... 149, 160	
cis-Vaccenic Acid ..... 133	$\gamma$ -Tocotrienol ..... 126	p-Tolunitrile ..... 155	Toluene ..... 159, 160	
trans-Vaccenic Acid ..... 133	$\delta$ -Tocotrienol ..... 126	N-(m-Toluoyl) glycine ..... 128	Triton X-100 ..... 29	
Vitamin D2 ..... 126	m-Tolunitrile ..... 24	N-(o-Toluoyl) glycine ..... 128		
Vitamin D3 ..... 126	o-Tolunitrile ..... 24	N-(p-Toluoyl) glycine ..... 128		
	p-Tolunitrile ..... 24	N-Vanillylnonanamide ..... 27		
<b><math>\pi</math>NAP</b>	Triamcinolone ..... 129		<b>CN-MS</b>	
Adenine ..... 135	Triamcinolone Acetonide ..... 129		Acetaminophen ..... 31	
Adenosine ..... 135	Uracil ..... 135	<b>NPE</b>	Acetophenone ..... 30	
Anisole ..... 155	Uridine ..... 135	Amylbenzene ..... 160	4'-Acetoxyacetanilide ..... 31	
Benzene ..... 144, 155	N-Vanillylnonanamide ..... 27	Benzene ..... 144, 160	p-Aminophenol Hydrochloride ..... 31	
Berberine Hydrochloride ..... 25		Butylbenzene ..... 160	Benzene ..... 30, 31	
Betamethasone ..... 129	<b>PYE</b>	C60 ..... 54, 59	Caffeine ..... 31	
Campesterol ..... 26	Amylbenzene ..... 160	C70 ..... 54	4-Chlorophenol ..... 85	
Capsaicin ..... 27, 120	Anisole ..... 155	1,2-Difluorobenzene ..... 144	Clofibrate ..... 85	
Cholesterol ..... 26	Anthracene ..... 157	1,3-Difluorobenzene ..... 144	4-Ethoxyphenol ..... 85	
Cortisone ..... 24, 129	Benz[a]anthracene ..... 158	1,4-Difluorobenzene ..... 144	Methyl Benzoate ..... 30	
Cortisone Acetate ..... 129	Benzene ..... 155, 157, 160	1,5-Dimethylnaphthalene ..... 28, 156	Theobromine ..... 31	
Cytidine ..... 129	2-Benzylpyridine ..... 148	1,5-Dinitronaphthalene(1,5-DNN) ..... 28, 156	Theophylline ..... 31	
Cytosine ..... 129	3-Benzylpyridine ..... 148	1,8-Dinitronaphthalene(1,8-DNN) ..... 28, 156	Toluene ..... 30, 31	
Daidzein ..... 25	4-Benzylpyridine ..... 148	Ethylbenzene ..... 160		
Daidzin ..... 25	Bromobenzene ..... 159	Fluorobenzene ..... 144	<b>C<sub>22</sub>-AR-II</b>	
Dexamethasone ..... 129	Butylbenzene ..... 160	Isophthalonitrile ..... 156	Acetophenone ..... 30	
1,2-Difluorobenzene ..... 144	C60 ..... 54	1-Methylnaphthalene ..... 28, 156	Benzene ..... 30	
1,3-Difluorobenzene ..... 144	C70 ..... 54	Naphthalene ..... 28, 156	Methyl Benzoate ..... 30	
		[6,6]-Phenyl-C61 Butyric Acid Butyl Ester ..... 59		

Toluene	30
<b>C<sub>8</sub>-MS</b>	
Acetophenone	30
Benzene	30
Bisoprolol	98
Butyl p-Hydroxybenzoate	108
Cetirizine Dihydrochloride	89
Clindamycin Phosphate	85
Cyanocobalamin	87
Cyanocobalamin Derivative	87
1,2-Dinitrobenzene	83
Dipyridamole	88
Domperidone	95
Ethyl p-Hydroxybenzoate	95
Fumaric Acid	98
Ipratropium Bromide	82
Ipratropium Bromide Derivative	82
Lincomycin Hydrochloride	85, 108
Methyl Benzoate	30
Methyl p-Hydroxybenzoate	85
Mexiletine Hydrochloride	104
4-epi-Minocycline	104
Minocycline Hydrochloride	104
Omeprazole	83
Phenethylamine Hydrochloride	104
Propyl p-Hydroxybenzoate	89
Rifampicin	108
Scopolamine n-Butyl Bromide	99
Scopolamine Hydrobromide	99
p-Terphenyl	88
Toluene	30
Vincristine Sulfate	99
Vinyl Acetate	117
1-Vinyl-2-pyrrolidone	117

**C<sub>4</sub>-MS**

Acetophenone	30
Benzene	30
Methyl Benzoate	30
Toluene	30

**TMS-MS**

Acetophenone	30
Benzene	30
Daunorubicin Hydrochloride	83
Doxorubicin Hydrochloride	83
Epirubicin Hydrochloride	83
Methyl Benzoate	30
2-Naphthalenesulfonic Acid Sodium	83

Ticarcillin Sodium	92
Toluene	30
o-Toluic Acid	92

**PE-MS**

Acetophenone	30
Anisole	155
Benz[a]anthracene	158
Benzene	30, 155
Berberine Hydrochloride	25
Capsaicin	116, 120
Chrysene	158
Cortisone	24
Daidzein	25
Daidzin	25
Dexamethasone	92
Ethyl Vanillin	104
Genistein	25
Genistin	25
Glycitein	25
Glycitin	25
Hydrocortisone	24
Maleic Acid	92
Methyl Benzoate	30
Mitomycin C	104
Nitrobenzene	155
Palmatine Hydrochloride	25
Prednisolone	24
Timolol Maleate	92
Toluene	30
m-Tolunitrile	24, 155
o-Tolunitrile	24, 155
p-Tolunitrile	24, 155
Vanillylnonanamide	116

**SL-II**

Acetaminophen	110
Acetanilide	110
Benzoic Acid	36
β-Carotene	123, 127
Cetirizine Dihydrochloride	89
Cinnamyl Alcohol	36
Dantrolene Sodium	91
4-Dimethylaminoantipyrine	89
Fluocinolone Acetonide	101
Fluoxymesterone	100
Hydrocortisone	98
Methyl Prednisolone	100
1-Naphthol	37
2-Naphthol	37
p-Nitrobenzyl Alcohol	36

Phenacetin	110
L-α-Phosphatidyl Choline	134
Prednisone	98
all-trans-Retinol Acetate	123
Salicylamide	36
Salicylic Acid	36
Theophylline	91
α-Tocopherol	127
β-Tocopherol	127
γ-Tocopherol	127
δ-Tocopherol	127
Tolazamide	94
Tolbutamide	94
Triamcinolone Acetonide	101
Vitamin A	127
p-Xylene	36

**HILIC**

Allantoic Acid	112
Allantoin	112
Ammelide	39
Ammeline	39
Angiotensin II, [Asn1,Val5]	138
Angiotensin II, [Sar1,Ile8]	138
Angiotensin II, [Val5]	138
Angiotensin II, Des-Asp1-[Ile8]	138
Ascorbic Acid	125, 126
L(+)-Ascorbic Acid	125
L-Citrulline	120
Cyanuric Acid	39
Diethylene Glycol	145
meso-Erythritol	146
Ethylene Glycol	145
D-Fructose-6-phosphate	132
Glucose	132
α-D-Glucose-1-phosphate	132
D-Glucose-6-phosphate	132
Glyceric Acid	146
Glycerol	145
Glycine	38
Glycylglycine	38
Isoascorbic Acid	125
Isoleucine	39
Leucine	39
Malic Acid	120, 126
Melamine	39
Nicotinamide	125
Nicotinic Acid	125
Oxalic Acid	146
Oxamic Acid	146
D-Pantothenic Acid	125
Pyridoxine (Vitamin B6)	125

Riboflavin (Vitamin B2)	125
Sorbic Acid	125
Taurine	137
Trimethylene Glycol	145
Tris (hydroxymethyl) aminomethane	146
Uracil	145
Urea	146
Uridine	145
Valine	39

**Sugar-D**

Arabinose	41
α-Cyclodextrin	130
β-Cyclodextrin	130
γ-Cyclodextrin	130
meso-Erythritol	41, 130, 131
Fructooligosaccharides	42
Fructose	41, 42, 130, 131
Galactose	41
Glucitol (Sorbitol)	41, 131
Glucose	40, 41, 42, 130, 131
Glucuronic Acid	42
Glycerol	41
Guanosine	111
Inositol	41, 131
Isomaltoligosaccharides	130
Lactose	43, 131
Maltitol	41, 130, 131
Maltoheptaose	42
Maltohexaose	42
Maltopentaose	42
Maltose	40, 41, 42, 43, 130, 131
Maltotetraitol	131
Maltotetraose	42
Maltotriitol	131
Maltotriose	42, 131
Mannitol	131
Mannosamine	42
Mannose	42
Palatinit	131
Palatinose	130
cis-Platin (CDDP)	111
Raffinose	41, 130
Rhamnose	41, 42
Sucrose	41, 43, 130, 131
Trehalose	43, 130
Xylitol	41, 130, 131
Xylooligosaccharides	130
Xylose	42

<b>NH<sub>2</sub>-MS</b>	
Arabinose	41
Fructose	42
Galactose	41
Glucose	40, 42
Lactose	43
Maltitol	41
Maltose	40, 43
Mannose	42
Rhamnose	42
Sucrose	43
Trehalose	43
Xylose	42
<b>Protein-R</b>	
Actinase E	139
Albumin (Bovine)	44, 139
Albumin (Goat)	140
Albumin (Human)	139
Albumin (Ovalbumin)	140
Albumin (Rabbit)	140
Alcohol Dehydrogenase	139
α-Amylase	139
Angiotensin I (Human)	137
Angiotensin II (Human)	137
Carbonic Anhydrase	139
Catalase	140
Choline Oxidase	139
α-Chymotrypsinogen A	140
Conalbumin	140
Concanavalin A	140, 141
Cytochrome C	44, 140
Deoxyribonuclease I	140
Diaphorase	141
Fibrinogen	140
Glucose Oxidase	139
L-Glutamic Dehydrogenase	139
Lysozyme	44, 140
Milk	139
Myoglobin	44, 140
Myosin	138
Ovalbumin	44
Oxytocin	137
pd(T)12-18	135
Peroxidase	140
Poly-L-Glutamic Acid Sodium Salt	45
Pyruvate Kinase	140
Ribonuclease A	44, 140
Soybean Milk	139
Substance P	137
Thermolysin	139

Thyroglobulin	140
Transferrin	139
<b>C<sub>18</sub>-AR-300</b>	
[Ala2]-Leu-Enkephalin	138
[Ala2]-Met-Enkephalin	138
Angiotensin II	138
Bombesin	138
Bradykinin	138
Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Leu-Enkephalin	138
Met-Enkephalin	138
Neurotensin	138
Substance P	138
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141
<b>C<sub>8</sub>-AR-300</b>	
Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141
<b>C<sub>4</sub>-AR-300</b>	
Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141
<b>Ph-AR-300</b>	
Carbonic Anhydrase	141
α-Chymotrypsinogen A	141
Hemoglobin, Bovine	46
Insulin	141
Transferrin	141
Trypsin Inhibitor, Soybean	141
Trypsinogen	141

<b>Diol-120-II</b>	
Aprotinin	48
Bovine Serum	142
Conalbumin	48
Egg White	142
Glucose Oxidase	48
Human Serum	142
Myoglobin	48
Peroxidase	48
Polyethylene Glycol	49
Thyroglobulin	48
<b>Diol-300-II</b>	
Aprotinin	48
Bovine serum	142
Conalbumin	48
E.coli proteins	49
Egg White	142
Glucose Oxidase	48
Human Serum	142
Myoglobin	48
Peroxidase	48
Pullulan	49
Sodium Polystyrene Sulfonate	49
Thyroglobulin	48
<b>DEAE</b>	
Lysozyme	51
Myoglobin	51
Transferrin	51
Trypsin Inhibitor	51
<b>QA</b>	
Lysozyme	51
<b>CM</b>	
Lysozyme	51
<b>SP</b>	
Lysozyme	51
<b>HIC</b>	
Albumin (Bovine)	53
β-Glucosidase	53
Hemoglobin	53
β-Lactoglobulin	53

Myoglobin	53
<b>Buckyprep</b>	
C60	54, 56, 59
C60O	56
C70	54, 56
C76	56
C84	56
C86	57
Crude fullerenes	58
Fullerene toluene extract	55
Gd@C82(I)	57
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59
Sc2@C76(I)	57
Sc2@C78	57
Sc2@C80(I)	57
<b>Buckyprep-M</b>	
C60	54
C70	54
C86	57
Gd@C82(I)	57
Sc2@C76(I)	57
Sc2@C78	57
Sc2@C80(I)	57
<b>PBB</b>	
C60	54, 58, 59
C70	54, 58
C76	58
C84	58
Crude fullerenes	58
[6,6]-Phenyl-C61 Butyric Acid Butyl Ester	59
[6,6]-Phenyl-C61 Butyric Acid Methyl Ester	59
<b>PBB-R</b>	
Allylbenzene	149
1,1'-Binaphthyl	149
Dibenzosuberenone	149
Dibenzosuberone	149
Diphenylmethane	149
Fluorene	149
Perylene	149
Polyoxyethylene Hexadecyl Ether	148
Polyoxyethylene Lauryl Ether	148
Propylbenzene	149



I. COSMOSIL HPLC column	<p><b>CNT-300</b></p> <p>DNA-wrapped CNT ..... 60</p>
II. Liquid chromatography related product	<p><b>CNT-1000</b></p> <p>DNA ..... 135</p> <p>DNA-wrapped CNT ..... 60</p> <p>Hyaluronic Acid ..... 132</p> <p><b>CNT-2000</b></p> <p>DNA ..... 135</p> <p>DNA-wrapped CNT ..... 60</p>
III. Application data	
IV. Technical note	
V. Index	

# Genore chromatografia

ul. Dynarska 1/23, 01-493 Warszawa  
tel. 22 40 107 34 do 35, fax: -36  
e-mail: [info@genore.pl](mailto:info@genore.pl)  
[www.genore.pl](http://www.genore.pl)

## **Warranties and Disclaimers:**

Nacalai Tesque warrants that its products shall conform to the description of such products as provided by Nacalai Tesque through its catalog, analytical data or other literature. Nacalai Tesque makes no other warranty, express or implied, as to the fitness of these products for any particular purpose. Nacalai Tesque shall not in any event be liable for incidental or consequential damages that may result from any use or failure of the products.

For more information on products and pricing, please contact your local distributor.

---

## **NACALAI TESQUE, INC.**

Nijo Karasuma, Nakagyo-ku  
Kyoto 604-0855 JAPAN  
TEL :+81-(0)75-251-1730  
FAX :+81-(0)75-251-1763  
Web site :<http://www.nacalai.com>  
E-mail :[info.intl@nacalai.com](mailto:info.intl@nacalai.com)